



Nov 26,
2019

Aichivirus 3C3D RT-PCR

Judy Northill¹


¹Public Health Virology, Forensic and Scientific Services

1 Works for me dx.doi.org/10.17504/protocols.io.2i7gchn

Judy A. Northill
Public Health Virology, Forensic and Scientific Services

ABSTRACT

This RT-PCR will detected Aichivirus A from human samples. It spans the junction region of 3C and 3D and is used for genotyping.



Schematic of Aichivirus A with the 3C3D RT-PCR target region.

GUIDELINES

Method assumes the user is familiar with the thermocycler and software used to run the protocol.
Mix preparation should occur in a different laboratory or room to the amplification or post-PCR area.

STEPS MATERIALS

NAME	CATALOG #	VENDOR
SuperScript™ III One-Step RT-PCR System with Platinum™ Taq High Fidelity DNA Polymerase	12574035	Thermo Fisher Scientific
5ml Ethidium Bromide Solution [0.625mg/ml]	R041	G-Biosciences
Agarose low EEO (Agarose Standard)	A21140100	AppliChem
100bp DNA Ladder, 250ul (50 lanes)	G2101	Promega

Oligonucleotide sequences

1	Name	5'-3'
	AiV-6213F	ACTGGGCCACCCTCCAGACG
	AiV-7044R	GGTTGATTTCAGCTTGAGTTC

Reaction set-up

2

- Prepare sufficient for number of reaction plus a 'dead volume', usually 2 extra. Adjust as necessary if using a robotic dispenser.



SuperScript™ III One-Step RT-PCR
System with Platinum™ Taq High
Fidelity DNA Polymerase
by Thermo Fisher Scientific
Catalog #: 12574035

Reagent	Vol (μL) x1	Final reaction concentration
Nuclease free water	3.6	
Primer AiV-6213F (20pmol/μl)	0.5	500nM
Primer AiV-7044R (20pmol/μl)	0.5	500nM
2 X Reaction mix	10	1X
Superscript III RT/Platinum Taq HiFi enzyme mix	0.4	
TOTAL VOLUME	15	

Dispense 15μL to each reaction well.

Add 5μL of template, extracted RNA, controls or NTC (nuclease-free water).
Total reaction volume is 20μL

Amplification

- 3 The assay has been used with Eppendorf thermocyclers.

PCR cycling times

1 cycle	40 cycles	1 cycle	Hold
50°C 30 minutes	94°C 15 seconds	68°C 5 minutes	15°C
94°C 2 minutes	50°C 30 seconds		
	68°C 60 seconds		

4

Amplified products are analysed by gel electrophoresis or equipment such as a QIAxcel.

For gel electrophoresis, a 1.5% agarose gel with ethidium bromide was made using the following recipe.

Reagent	Volume
Agarose	1.5g
0.5 X TBE buffer	100ml

Boil in microwave for 1.5 -2 minutes until agarose powder is dissolved.

Add 2 drops of ethidium bromide (0.625mg/ml) and mix before pouring warm into the gel form with a comb.

Leave for approx 30minutes at room temperature to set.



For the novice user to gel electrophoresis there a detailed protocol at <https://dx.doi.org/10.17504/protocols.io.s38egrw> that will give you further guidance.

Experienced users may run their favorite gel recipe with appropriate dye. The aim is to visualise the amplified products. Gel running times and voltage can be adjusted depending on the equipment used.



5ml Ethidium Bromide Solution [0.625mg/ml]

by G-Biosciences

Catalog #: R041



Agarose low EEO (Agarose Standard)

by AppliChem

Catalog #: A21140100

CAS Number: 9012-36-6



100bp DNA Ladder, 250ul (50 lanes)

by Promega


Catalog #: G2101

Place gel (on a tray) into the tank and cover with 0.5X TBE with ethidium bromide.

Mix 2µl of loading dye with 10µl of amplified product and add to a well in the gel.

Mix 2µl of loading dye with 5µl of 100bp marker and add to a well in the gel. Ideally, the first and last lane but this is not essential.

Run the gel for 60-90 minutes at 80 volts or until the bands on the marker have separated adequately without run off the end of the gel.



Example of 3C-3D products analysed using a QIAxcel.



Amplified products will produce a band of 832bp in size. No template controls (NTC) should be not detected.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited