

ChIP protocol

Kseniya Obraztsova

Abstract

This is a ChIP protocol for mammalian cells using Millipore ChIP assay kit.

(Catalog # 17-295
Lot # 2135690)

Citation: Kseniya Obraztsova ChIP protocol. **protocols.io**
dx.doi.org/10.17504/protocols.io.njqdcmw

Published: 01 Jun 2018

Protocol

Step 1.

After passaging 2-3 times, plate cells on 10 cm dishes with appropriate density for them to reach 10^6 cells in two days. Have one control plate just for counting.

Step 2.

Cross link histones to DNA by adding formaldehyde directly to culture medium to a final concentration of 1% and incubate for 10 minutes at 37°C. (For example, add 270 μ L 37% formaldehyde into 10 mL of growthmedium on plate).

Step 3.

Aspirate medium, removing as much medium as possible. Wash cells twice using ice cold PBS containing protease inhibitors (1 mM phenylmethylsulfonyl fluoride (PMSF), 1 μ g/mL aprotinin and 1 μ g/mL pepstatin A). Note: Add protease inhibitors to PBS just prior to use. PMSF has a half-life of approximately 30 minutes in aqueous solutions.

Step 4.

Scrape cells into conical tube.

Step 5.

Pellet cells for 4 minutes at 2000 rpm at 4°C. Warm SDS Lysis Buffer (Catalog # 20-163) to room temperature to dissolve precipitated SDS and add protease inhibitors (inhibitors: 1mM PMSF, 1 μ g/mL aprotinin and 1 μ g/mL pepstatin A).

Step 6.

Resuspend cell pellet in 200 µL of SDS Lysis Buffer (Catalog # 20-163) and incubate for 10 minutes on ice.

Step 7.

Sonicate lysate to shear DNA to lengths between 200 and 1000 basepairs being sure to keep samples ice cold (Note: Once sonication conditions have been optimized following steps 1 to 9, proceed to Part B, step 1 below).

Step 8.

Add 8 µL 5 M NaCl (Catalog # 20-159) and reverse crosslinks at 65°C for 4 hours.

Step 9.

Recover DNA by phenol/chloroform extraction and run sample (example 5 µL, 10 µL and 20 µL samples) in an

agarose gel to visualize shearing efficiency.

Part B. Experimental protocol.

If sonication conditions have been optimized (Part A), complete steps 1 through 7 and continue with the protocol below. For a negative/background control, prepare a sample to use as a no-antibody immunoprecipitation control in step 5 below. Additionally, transcriptionally unactivated DNA samples should be prepared as controls for PCR in section II.

10)

Centrifuge samples (from part A, step 7) for 10 minutes at 13,000 rpm at 4°C, and transfer the supernatant to a new 2 mL-microcentrifuge tube. Discard pellet.

Step 10.

Dilute the sonicated cell supernatant 10 fold in ChIP Dilution Buffer (Catalog # 20-153), adding protease inhibitors as above. This is done by adding 1800 µL ChIP Dilution Buffer to the 200 µL sonicated cell supernatant for a final volume of 2 mL in each immunoprecipitation condition. Note: If proceeding to PCR a portion of the diluted cell supernatant 1% (20 µL) can be kept to quantitate the amount of DNA present in different samples at the PCR protocol, Part B, section II, step 6. This sample is considered to be your input/starting material and needs to have the Histone-DNA crosslinks reversed by adding 1 µL of 5 M NaCl and heating at 65°C for 4 hours (see section II, step 3).

28820 Single Oak Drive • Temecula, CA 92590

Technical Support: T: 800 437-7500 • F: 800 437-7502 www.millipore.com

Page Four of Four Catalog # 17-295 Lot # 2135690

12)

To reduce nonspecific background, pre-clear the 2 mL diluted cell supernatant with 75 µL of Protein A Agarose/Salmon Sperm DNA (50% Slurry), (Catalog # 16-157C) for 30 minutes at 4°C with agitation.

Step 11.

Pellet agarose by brief centrifugation and collect the supernatant fraction.

Step 12.

Add the immunoprecipitating antibody (the amount will vary per antibody) to the 2 mL supernatant fraction and incubate overnight at 4°C with rotation. For a negative control, perform a no-antibody immunoprecipitation by incubating the supernatant fraction with 60 µL of Protein A Agarose/Salmon Sperm DNA (50% Slurry), (Catalog # 16-157C) for one hour at 4°C with rotation and proceed to step 7.

Step 13.

Add 60 µL of Protein A Agarose/Salmon Sperm DNA (50% Slurry), (Catalog # 16-157C) for one hour at 4°C with rotation to collect the antibody/histone complex.

Step 14.

Pellet agarose by gentle centrifugation (700 to 1000 rpm at 4°C, 1min). Carefully remove the supernatant that contains unbound, non-specific DNA. Wash the protein A agarose/antibody/histone complex for 3-5 minutes on a rotating platform with 1 mL of each of the buffers listed in the order as given below:

17)

a) Low Salt Immune Complex Wash Buffer (Catalog # 20-154), one wash

Step 15.

b) High Salt Immune Complex Wash Buffer (Catalog # 20-155), one wash

Step 16.

c) LiCl Immune Complex Wash Buffer (Catalog # 20-156), one wash

Step 17.

d) TE Buffer (Catalog # 20-157), two washes

After step 7 above, the sample is now a protein A/antibody/histone/DNA complex ready for either an Immunoprecipitation/Immunoblot assay (Section I) or Polymerase Chain Reaction (PCR) assay (Section II):

Section I. Immunoprecipitation/Immunoblot protocol to detect histone.

1. Following washing of the beads in part B, step 7, immunoprecipitated histones can be analyzed by immunoblot analysis. Add 25 µL of 1X Laemmli buffer per sample and boil for 10 minutes. Load 20 µL

per lane and perform immunoblot procedure as described per appropriate antibody.

Section II. PCR protocol to amplify DNA that is bound to the immunoprecipitated histone.

21)

Freshly prepare elution buffer (1% SDS, 0.1 M NaHCO₃).

Step 18.

Elute the histone complex from the antibody by adding 250 µL elution buffer to the pelleted protein A agarose/antibody/histone complex from step 7d above. Vortex briefly to mix and incubate at room temperature for 15 minutes with rotation. Spin down agarose, and carefully transfer the supernatant fraction (eluate) to another tube and repeat elution. Combine eluates (total volume = 500 µL).

Step 19.

Add 20 µL 5 M NaCl (Catalog # 20-159) to the combined eluates (500 µL) and reverse histone-DNA crosslinks by heating at 65°C for 4 hours. At this step the sample can be stored at -20°C and the protocol continued the next day.

Note: Include the input/starting material (the sample saved from Part B, step 2, which has had the histone-DNA crosslinks reversed by adding 1 µL of 5 M NaCl per 20 µL sample and heating to 65°C for 4 hours) as well as a transcriptionally-unactivated DNA sample as negative and background controls for the PCR reaction. Previously, a 5 µL sample has been used in a nested PCR reaction. However, the amount of sample used per reaction must be determined empirically (e.g., titrate the sample at this step by using 1, 2, 5, or 10 µL per PCR reaction). If PCR results are poor, complete steps 4, 5 and 6 below to purify the DNA sample. NOTE: Handle the samples carefully; some DNA may be lost during the purification steps.

Step 20.

Add 10 µL of 0.5 M EDTA (Catalog # 20-158), 20 µL 1 M Tris-HCl, pH 6.5 (Catalog # 20-160) and 2 µL of 10 mg/mL Proteinase K to the combined eluates and incubate for one hour at 45°C.

Step 21.

Recover DNA using ChIP DNA concentrator kit. (Zymo Research).

Step 22.

Resuspend pellets in an appropriate buffer for PCR or slot-blot reactions. PCR or slot-blot conditions must be determined empirically.