

# Quick-DNA™ FFPE Kit Protocol

# Zymo Research

# **Abstract**

The **Quick-DNA™ FFPE Kit** provides a simple and reliable method for high yield/quality DNA isolation from formalin-fixed, paraffin embedded (FFPE) tissue samples and sections. The unique chemistries of the product have been optimized for maximum recovery of non-crosslinked, ultrapure DNA without RNA contamination. Simply digest deparaffinized tissues using the provided **Proteinase K**, heat, and then purify the DNA with the Fast-Spin columns in the kit. DNA >50 bp or >500 bp can be selectively isolated by altering the lysis buffer conditions as given in the protocol. PCR inhibitors are effectively removed during the isolation procedure, and eluted DNA is ideal for PCR, Next-Gen library prep, enzymatic manipulation, etc.

Citation: Zymo Research Quick-DNA™ FFPE Kit Protocol. protocols.io

dx.doi.org/10.17504/protocols.io.giybufw

Published: 21 Dec 2016

### **Guidelines**

#### **Product Contents**

Quick-DNA <sup>™</sup> FFPE Kit (Kit Size)	<b>D3067</b> (50 Preps.)	Storage Temperature
Deparaffinization Solution	20 ml	Room Temp.
Proteinase K & Storage Buffer <sup>1</sup>	2 x 5 mg	-20°C (after mixing)
2X Digestion Buffer	5 ml	Room Temp.
Genomic Lysis Buffer <sup>2</sup>	50 ml	Room Temp.
Genomic DNA Wash 1	25 ml	Room Temp.
Genomic DNA Wash 2 <sup>3</sup> (concentrate)	12 ml	Room Temp.
DNA Elution Buffer	10 ml	Room Temp.
RNase A <sup>4</sup>	2 mg	4°C
Zymo-Spin™ IIC Columns	50	Room Temp.
Collection Tubes	100	Room Temp.
Instruction Manual	1	-

<sup>&</sup>lt;sup>1</sup> The Proteinase K is stable as shipped. Add 260 μl **Proteinase K Storage Buffer** to each **Proteinase K** tube prior to use. The final concentration of **Proteinase K** after the addition of **Proteinase K** storage Buffer is ~20 mg/ml. Store at -20° C.

# **Specifications**

- **Sample Size** Up to 25 mg tissue from paraffin block or up to four (4) tissue sections (≤20 μm thick) with a total surface area ~20 cm². It is recommended to use 1-2 sections if performing the protocol for the first time. Compatible with fresh/frozen tissue specimens.
- **DNA Recovery**  $\square$  High quality total DNA (A<sub>260/</sub>A<sub>280</sub> >1.8) can be eluted into small volumes (i.e.,  $\ge 25 \ \mu$ l) allowing for highly concentrated samples. The maximum DNA binding capacity of the provided spin column is  $\sim 25 \ \mu$ g.
- **Processing Time** As little as 4 hours when processing large amounts of tissue. For maximum yields of the highest quality DNA, it is recommended to process samples overnight.
- **Equipment/Reagents** Microcentrifuge, thermomixer or heat block/bath capable of 55°C and 90°C, isopropanol, beta-mercaptoethanol (optional).


# **Product Description**

The **Quick-DNA™ FFPE Kit** provides a simple and reliable method for high yield/quality DNA isolation from formalin-fixed, paraffin embedded (FFPE) tissue samples and sections. The unique chemistries of the product have been optimized for maximum recovery of non-crosslinked, ultra-pure DNA without RNA contamination. Simply digest deparaffinized tissues using the provided **Proteinase K**, heat, and then purify the DNA with the Fast-Spin columns in the kit. DNA >50 bp or >500 bp can be selectively isolated by altering the lysis buffer conditions as given in the protocol. PCR inhibitors are effectively removed during the isolation procedure, and eluted DNA is ideal for PCR, Next-Gen library prep, enzymatic manipulation, etc. Shown below is a schematic and performance overview of the procedure.

<sup>&</sup>lt;sup>2</sup> Recommended: Add beta-mercaptoethanol to 0.5(v/v) i.e., 250 μl per 50 ml or 500 μl per 100 ml.

<sup>&</sup>lt;sup>3</sup> Before starting, add 48 ml 100% ethanol (52 ml 95% ethanol) to the 12 ml **Genomic DNA Wash 2** concentrate.

<sup>&</sup>lt;sup>4</sup> Re-suspend lyophilized RNase A in 300 μl of ddH<sub>2</sub>O. Store at 4° C.






Figure 1. Equivalent amounts of DNA isolated using Zymo and Supplier Q procedures were used for real time PCR analysis. DNA isolated using the *Quick-DNA™ FFPE* Kit consistently yielded lower *Ct* values values as depicted by the amplification curves above.




Figure 2. Equivalent amounts of DNA resolved in a 1% agarose/TAE/EtBr gel show binding conditions may be adjusted with the *Quick*-DNA™ FFPE Kit to selectively isolate DNA >50 bp or >500 bp. 100 bp DNA ladder and 1 kb DNA ladder from Zymo Research

# **Buffer Preparation**

Add 260  $\mu$ l **Proteinase K Storage Buffer** to reconstitute lyophilized **Proteinase K** at 20 mg/ml. Vortex to dissolve. Store at -20° C.

Before starting, add 48 ml 100% ethanol (52 ml 95% ethanol) to the 12 ml **Genomic DNA**Wash 2 concentrate.

Resuspend lyophilized **RNase A** in 300  $\mu$ l of ddH<sub>2</sub>0. Store at 4° C.

Recommended: Add beta-mercaptoethanol (user supplied) to the Genomic Lysis Buffer to a final

dilution of 0.5%(v/v) i.e., 250  $\mu$ l per 50 ml.

# **Before start**

# **Buffer Preparation**

Add 260  $\mu$ l **Proteinase K Storage Buffer** to reconstitute lyophilized **Proteinase K** at 20 mg/ml. Vortex to dissolve. Store at -20° C.

Before starting, add 48 ml 100% ethanol (52 ml 95% ethanol) to the 12 ml **Genomic DNA**Wash 2 concentrate.

Resuspend lyophilized **RNase A** in 300  $\mu$ l of ddH<sub>2</sub>0. Store at 4° C.

Recommended: Add beta-mercaptoethanol (user supplied) to the **Genomic Lysis Buffer** to a final dilution of 0.5%(v/v) i.e.,  $250~\mu l$  per 50~m l.

# **Materials**

Quick-DNA™ FFPE Kit <u>D3067</u> by <u>Zymo Research</u>

# **Protocol**

# Deparaffinization

# Step 1.

Remove (trim) excess paraffin wax from sample and transfer the sample to a 1.5 ml microcentrifuge tube.

#### NOTES

Shea Biondi 25 Nov 2016

**Note:** Up to 25 mg tissue from a paraffin block or up to four (4) tissue sections ( $\leq$ 20  $\mu$ m thick) with a 2 total surface area  $\sim$ 20 cm . It is recommended to use start with 1-2 sections.

#### Shea Biondi 25 Nov 2016

**Note:** If using fresh/frozen tissue specimens proceed directly with Proteinase K Digestion & DNA Isolation.

#### Deparaffinization

#### Step 2.

Add 400  $\mu$ l of **Deparaffinization Solution** to the sample. Incubate at 55°C for 1 minute. Vortex briefly.

© DURATION

00:01:00

NOTES

Shea Biondi 25 Nov 2016

**Note:** Xylene may also be used for deparaffinization. See the <u>Appendix</u> for more instructions.

### Deparaffinization

# Step 3.

Remove **Deparaffinization Solution** from the sample and proceed to next section.

# **Tissue Digestion**

# Step 4.

To the deparaffinized tissue sample ( $\leq$  25 mg) in a microcentrifuge tube, add the following mixture:

H <sub>2</sub> O	45µl
2X Digestion Buffer	45µl
Proteinase K	10μΙ

#### **P** NOTES

Shea Biondi 25 Nov 2016

**Note:** If the tissue sample is too large for the digestion volume, scale up the digestion to 200  $\mu$ l while keeping the amount of Proteinase K the same. Double the reagent volumes indicated in Step 8 & 9 of the DNA Purification section.

# Tissue Digestion

# Step 5.

Digest 1-4 hours or overnight at 55 °C.

# **Rapid Digestion**

# Standard Digestion

Incubate at 55°C for 1-4 hours Incubate at 55°C overnight (12-16 hrs)

### NOTES

Shea Biondi 25 Nov 2016

**Note:** The Rapid Digestion is recommended for processing slide tissue sections. The Standard Digestion ensures maximum yields of DNA from tough-to-lyse (collagen-rich, fibrous, etc.) or large tissue samples.

# Tissue Digestion

#### Step 6.

Transfer the digestion to 94°C and incubate for 20 minutes.

#### **O** DURATION

00:20:00

#### Tissue Digestion

#### Step 7.

Once done, add 5 µl of **RNase A**, mix, and incubate an additional 5 minutes at room temperature.

# **O DURATION**

00:05:00

#### **DNA Purification**

#### Step 8.

Add 350 µl of **Genomic Lysis Buffer** to the tube and mix thoroughly by vortexing.

#### **DNA Purification**

# Step 9.

Add 135 µl of isopropanol<sup>1</sup> (user supplied) to the sample and mix thoroughly.

#### NOTES

#### Shea Biondi 25 Nov 2016

This procedure will isolate total DNA > 50 bp. To isolate only DNA > 500 bp, skip this step.

FFPE DNA may be highly degraded and DNA >500 bp may not be present in sample.

# **DNA Purification**

#### Step 10.

Centrifuge at  $\geq$  12,000 x g for 1 minute to remove insoluble debris.

**O DURATION** 

00:01:00

<sup>&</sup>lt;sup>1</sup>ssDNA will also be purified if present in the sample upon the addition of isopropanol.

NOTES

Shea Biondi 25 Nov 2016

This procedure will isolate total DNA > 50 bp. To isolate only DNA > 500 bp, skip this step.

#### **DNA Purification**

**Step 11.** 

Transfer the supernatant to a **Zymo-Spin™ IIC Column**<sup>2</sup> in a **Collection Tube**.

#### NOTES

Shea Biondi 25 Nov 2016

The maximum loading volume for the Zymo-Spin™ Column is ~700 μl.

# **DNA Purification**

**Step 12.** 

Centrifuge at 10,000 x g for 1 minute.

**O DURATION** 

00:01:00

#### **DNA Purification**

**Step 13.** 

Add 400 µl of **Genomic DNA Wash 1** to the spin column in a <u>new</u> **Collection Tube**.

# **DNA Purification**

#### **Step 14.**

Centrifuge at 10,000 x g for 1 minute. Discard the flow-through.

© DURATION

00:01:00

#### **DNA Purification**

**Step 15.** 

Add 700 µl of **Genomic DNA Wash 2** to the spin column.

# **DNA Purification**

**Step 16.** 

Centrifuge at  $\geq$  12,000 x g for 1 minute. Discard the flow-through.

**O DURATION** 

00:01:00

#### **DNA Purification**

#### **Step 17.**

Add 200 µl of **Genomic DNA Wash 2** to the spin column.

# **DNA Purification**

Step 18.

Centrifuge at  $\geq$  12,000 x g for 1 minute.

**O DURATION** 

00:01:00

# **DNA Purification**

# Step 19.

Transfer the **Zymo-Spin**<sup>™</sup> **IIC Column** to a clean microcentrifuge tube. Add  $\geq$  50  $\mu$ l **DNA Elution Buffer**<sup>3</sup> or water (add  $\geq$ 100  $\mu$ l if sampling 25 mg tissue) to the spin column.

#### NOTES

Shea Biondi 25 Nov 2016

<sup>3</sup> Elution of DNA from the column is dependent on pH and temperature. If water is used, ensure the pH is >6.0. Also, the total yield may be improved by eluting the DNA with Elution Buffer or water pre-equilibrated to 60-70° C or by performing and pooling sequential elutions.

#### **DNA Purification**

#### Step 20.

Incubate 2-5 minutes at room temperature.

**O** DURATION

00:05:00

#### **DNA Purification**

#### Step 21.

Centrifuge at top speed for 30 seconds to elute the DNA.

**O DURATION** 

00:00:30

#### **DNA Purification**

# Step 22.

The eluted DNA can be used immediately for molecular based applications or stored ≤-20°C for future use.