

Re-amplification of CRISPRa and CRISPRi libraries

Pei-Chun Lin

Abstract

Re-amplification of CRISPRa and CRISPRi Libraries

CRISPRi/a sublibraries

		Dilute to 50 ng/ul			
		sgRNA #	[ul]	MegaX [ul]	Large plate #
TssLib_Sub1	Apoptosis+Cancer+Other_Cancer	37354	3	36	1
TssLib_Sub2	Drug_Targets+Kinase_Phosphatase	28340	2	24	1
TssLib_Sub3	Gene_Expression	28430	2	24	1
TssLib_Sub4	Membrane_Proteins	13417	1	12	1
TssLib_Sub5	Stress_Proteostasis	34556	3	36	1
TssLib_Sub6	Trafficking+Mitochondria+Motility	26621	2	24	1
TssLib_Sub7	Unassigned	24000	2	24	1

Citation: Pei-Chun Lin Re-amplification of CRISPRa and CRISPRi libraries. **protocols.io**

dx.doi.org/10.17504/protocols.io.dkv4w5

Published: 13 Aug 2015

Guidelines

Reference: Genome-Scale CRISPR-Mediated Control of Gene Repression and Activation Giltbert LA et al. Cell. 2014 Oct 23; 159(3):647-61.doi:10.1016/j.cell.2014.09.029.PMID:25307932

[http://www.cell.com/cell/abstract/S0092-8674\(14\)01178-7](http://www.cell.com/cell/abstract/S0092-8674(14)01178-7)

Protocol

Step 1.

Dilute each sub-library to 50 ng/ul in water or EB

Step 2.

Electroporate the library

Step 3.

Pre-chill 0.1 cm cuvettes, megaX cells, 10% glycerol on ice

Step 4.

CRISPRi/a sublibraries

Dilute to 50 ng/ul					
		sgRNA #	[ul]	MegaX [ul]	Large plate #
TssLib_Sub1	Apoptosis+Cancer+Other_Cancer	37354	3	36	1
TssLib_Sub2	Drug_Targets+Kinase_Phosphatase	28340	2	24	1
TssLib_Sub3	Gene_Expression	28430	2	24	1
TssLib_Sub4	Membrane_Proteins	13417	1	12	1
TssLib_Sub5	Stress_Proteostasis	34556	3	36	1
TssLib_Sub6	Trafficking+Mitochondria+Motility	26621	2	24	1
TssLib_Sub7	Unassigned	24000	2	24	1

Step 5.

Follow the table for the amounts of sub-library plasmid DNA and MegaX competent cells, mix gently and incubate on ice for 30 min



REAGENTS

MegaX DH10B [C6400-03](#) by [Thermo Scientific](#)



DURATION

00:30:00

Step 6.

Add pre-chilled 10% glycerol to the MegaX-library mix for a final 75 ul, transfer the mix to a prechilled 0.1 cm cuvette

Step 7.

Electroporate at 2.0 kV, 200 ohms, 25 uF (Gene Pulser Xcell, Bio-rad)

Step 8.

Transfer cells to a culture tube

Step 9.

Use 1 ml pipettes and gel loading tips

Step 10.

Wash cells out gently with 300 ul S.O.C. twice (total 600 ul)



AMOUNT

300 µl Additional info:



PROTOCOL

. [SOC Media](#)

CONTACT: [New England Biolabs](#)

Step 10.1.

SOB Media



PROTOCOL

. [SOB Media](#)

CONTACT: [New England Biolabs](#)

Step 1.1.

2% tryptone

Step 1.2.

0.5% yeast extract

Step 1.3.

10 mM NaCl

Step 1.4.

2.5 mM KCl

Step 1.5.

10 mM MgCl₂

Step 1.6.

10 mM MgSO₄

Step 10.2.

20 mM glucose

Step 11.

Incubate at 37°C, 250 rpm, 1.5 hour

 **DURATION**

01:30:00

Step 12.

Plate the transformations

Step 13.

Plate all in one large square plate per sub-library, use autoclaved beads

Step 14.

Incubate at 37°C for 18 hours

 **DURATION**

18:00:00

Step 15.

Collect all colonies with LB and do one maxiprep per plate, elute in 500 ul EB; an ideal concentration is about 23 ug/ul

Sequencing

Step 16.

To sequence the library, you can PCR the sgRNA region with the following primers (a common 3' primer, with 2 different 5' primers for CRISPRa or CRISPRi)

- CRISPRi TSS common 3' caagcagaagacggcatatcgaCGACTCGGTGCCACTTTTTC
 - oCRISPRi TSS_1 (TruSeq Index 12 CTTGTA),
aatgatacggcgaccaccgaGATCGGAAGAGCACACGTCTGAACTCCAGTCACCTTGTAgcacaaaaggaaactcacct
 - oCRISPRi TSS_2 (TruSeq Index 6 GCCAAT)
aatgatacggcgaccaccgaGATCGGAAGAGCACACGTCTGAACTCCAGTCACGCCAATgcacaaaaggaaactcacct
- a. Pool sub-libraries proportionally (based on the number of sgRNAs) to have the CRISPRa or CRISPRi library, measure the pooled concentration and dilute it to 400 ng/ul for PCR
- b. Run 3 tubes of 100 ul PCR reactions

Sequencing

Step 17.

index primer

CRISPRi oCRISPRi TSS_1

CRISPRa oCRISPRi TSS_2

	x3
	uL
Library (400ng/ul)	1
5x Q5 Buffer	20
5x Q5 GC Buffer	20
dNTPs (10mM)	2
Index Primer (10uM)	2.5
Common Primer (10uM)	2.5
Q5 Hoststart polymerase	1
H2O	51

PCR

- 1, 98C 30s
- 2, 98C 15s
- 3, 58C 15s
- 4, 72C 15s
- 5, 2-4x23 cycles
- 6, 72C 10 m
- 7, 12C hold

Sequencing

Step 18.

After PCR, combine the 3 tubes and proceed to purify the PCR product (size 270 bp)

Sequencing

Step 19.

Sequencing primer, oCRISPRi TSS_seq V2, gtgtgttttgagactataaGtatcccttgagaaCCAcctGTTGG