

# Staining of survival motor neuron (SMN) protein in peripheral blood mononuclear cells (PBMC) protocol

Noriko Otsuki

## Abstract

Functional SMN protein of peripheral blood-derived mononuclear cells was detected with an anti-SMN antibody labeled with a fluorescent dye and analyzed semiquantitatively using intracellular expression intensity and SMN spot formed in cell nucleus as an index, Consider the relationship with the clinical condition and motor function.

Briefly,

- 1) In order to identify the cell fraction, the fluorescently labeled cell surface antigen-specific antibody is added to the whole blood sample and stained.
- 2) Hemolyze peripheral blood and fix mononuclear cells.
- 3) Intracellular staining is performed with SMN protein and nuclear-specific fluorescently labeled antibody.
- 4) SMN protein expression analysis was performed on a strongly CD33 cell population detected as a cell surface marker, SMN spot detection algorithm was used to analyze the proportion of cells in which intracellular SMN expression and SMN protein aggregated in the nucleus.

**Citation:** Noriko Otsuki Staining of survival motor neuron (SMN) protein in peripheral blood mononuclear cells (PBMC) protocol. **protocols.io**

<https://www.protocols.io/view/staining-of-survival-motor-neuron-smn-protein-in-p-qkdvuw>

**Published:** 06 Jun 2018

## Materials

- ✓ Clear Back MTG-001 by Contributed by users
- BD Phosflow™ Lyse/Fix Buffer (5 x) 5558049 by [BD Biosciences](#)
- BD Phosflow™ Perm/Wash Buffer I (10 x) 557885 by [BD Biosciences](#)
- BD Pharmingen™ Stain Buffer (FBS) 554656 by [BD Biosciences](#)
- BD™ CompBeads 51-90-9001229 by [BD Biosciences](#)

- ✓ AF488-anti human SMN mAb by Contributed by users
- 🧪 AF488-MOPC21 by [BioLegend](#)
- 🧪 BV610-anti human CD3 mAb by [BioLegend](#)
- 🧪 R-PE-anti human CD19 mAb by [BioLegend](#)
- ✓ PE-Cy5-anti human CD33 mAb by Contributed by users
- ✓ Hoechst 33342 H3570 by Contributed by users

## Protocol

### Staining cell surface molecules

#### Step 1.

- (1) Dispense 1.5 mL of heparin-treated whole blood into 15 mL\_conical tube
- (2) Add 30  $\mu$ L of normal human Ig solution (Clear Back, MTG-001, MBL, Nagoya, Japan)
- (3) Gently shake, leave in the dark at room temperature for 15 minutes
- (4) Add 5  $\mu$ L each of PE-HIB19, PE-Cy5-WM53, BV610-UCHT1 (Biolegend)
- (5) Gently shake, leave in the dark at room temperature, for 30 minutes

### Hemolysis of erythrocytes and immobilization of cell membranes

#### Step 2.

- (6) Warm lysis/fixation buffer and PBS (-) to 37 °C in a water bath
- (7) Add 10 mL of pre-warmed lysis/fixation buffer (37 °C) and mix immediately upside down
- (8) Incubate in a water bath at 37 °C for 10 minutes
- (9) Mix by inverting, centrifuge (900  $\times g$ , 5 min)
- (10) Discard the supernatant using an aspirator
- (11) Resuspended in PBS (-) (37 °C)
- (12) Centrifuge (900  $\times g$ , 5 min)
- (13) Discard supernatant

### Permeabilization of transmembrane of cells

#### Step 3.

- (14) Resuspend with 1.5 mL of permeabilization buffer
- (15) Incubate at room temperature in the dark for 30 min
- (17) Discard supernatant

(18) Resuspended with 1.5 mL of permeabilization buffer

(19) Centrifugation ( $900 \times g$ , 5 min)

(20) Discard supernatant

#### Cell count and Intracellular staining

##### **Step 4.**

(21) Cell count

(22) Dispensing the  $1 \times 10^6$  cells/50  $\mu$ L in Permeabilization buffer into two microcentrifuge tubes (for stained with specific Ab and isotype control)

(23) Tube for SMN staining: Add 1  $\mu$ g of AF488-2B1

(24) Tube for isotype control staining: Add 1 test of AF488-MOPC21 (5  $\mu$ L)

(25) Incubate stained cells (both 23 and 24) at room temperature in the dark for 45 min

(26) Add 500  $\mu$ L of permeabilization buffer and centrifuge ( $500 \times g$ , 5 min)

(27) Repeat (26) twice, and discard supernatant

#### Nuclear staining

##### **Step 5.**

(28) Suspend in 50  $\mu$ L of Hoechst 33342 (final concentration: 0.25  $\mu$ g/mL)

(29) Incubate in the dark at room temperature for 15 min

(30) Add 500  $\mu$ L of PBS (-) and centrifuged ( $500 \times g$ , 5 min)

(31) Discard supernatant

(32) Suspend with 50  $\mu$ L of PBS (-)

(33) Analysis