



Feb 27, 2019

Working

Splitting 96 Well Plates for gDNA Extraction and Freezing Down

In 1 collection

Celeste Karch¹, Rita Martinez¹, Jacob Marsh¹¹Washington University in St Louis

dx.doi.org/10.17504/protocols.io.x79frr6

Neurodegeneration Method Development Community

Tech. support email: ndcn-help@chanzuckerberg.comCeleste Karch
Washington University in St LouisComprehensive Genomic
Editing and Screening
Protocol Updated
02142019.docx

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

This protocol is part of the [Screening Edited iPSC Clones collection](#).





SAFETY WARNINGS

Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.

BEFORE STARTING

Split cells upon reaching maximum 80% confluence and minimum 40%.

- 1 Aspirate media from original plate.
- 2 Wash with 200 μ l PBS and aspirate.
- 3 Add 25 μ l of 0.05% Trypsin.
- 4 Incubate at 37 °C for 00:05:00
- 5 Tap to lift cells from plate.

- 6 Check under microscope to ensure that cells have detached from plate.
- 7 Add  **45 µl** mTesR1 to plate and tap to mix.
- 8 Transfer  **15 µl** to a 96 well PCR plate, while maintaining the location of each sample (this plate will be used for gDNA extraction).
- 9 Add  **50 µl** of 2X Freezing Media (20% DMSO in FBS) to cell plate and tap to mix.
- 10 Wrap plate in parafilm. Add tape over parafilm. Label the plate on the outside of the tape.
- 11 Place plate in Styrofoam box. Fill any open areas with diapers, paper towels, or Kimwipes. Cover box and make sure that lid is closed completely.
- 12 Store in  **-80 °C** for up to 4 months .



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited