CviRI Purification From XZ-6E Virus Infected NC64A Chlorella

David Dunigan and Irina Agarkova

Abstract

Citation: David Dunigan and Irina Agarkova CviRI Purification From XZ-6E Virus Infected NC64A Chlorella. protocols.io

dx.doi.org/10.17504/protocols.io.er4bd8w

Published: 13 Jun 2016

Guidelines

MATERIALS AND ASSAY CONDITIONS:

- 1) 3-4 hour XZ-6E virus infected NC64A, pellets frozen at -80°C
- 2) Buffer A:
 - 10 mM Tris-HCl, pH 8.0
 - 10 mM 2-ME
 - 50 μg/mL PMSF
- 3) Buffer B:
 - 20 mM Tris-HCl, pH 7.8
 - 0.5 mM EDTA
 - 7 mM 2-ME, 10% Glycerol
- 4) Buffer C:
 - 10 mM KHPO₄, pH 7.4
 - 0.5 mM EDTA
 - 7 mM 2-ME, 10% Glycerol
- 5) Storage buffer:
 - 20 mM Tris-HCl, pH 8.0
 - 0.5 mM EDTA
 - 0.1 mM DTT
 - 50 mM KCl
 - 5 mM MgCl₂, 50% Glycerol
- 6) 4 M NaCl
- 7) 28% Polyethyleneglycol (PEG) 8000
- 8) Assay buffer and conditions:

1X CviRI assay buffer:

- 10 mM Tris-HCl, pH 7.8
- 150 mM KCl
- 10 mM MgCl2

• 100 μg/mL BSA

All assays are carried out in 20.0 μ L volumes with 1 μ g of pUC19 DNA as substrate for 60 to120 min at 25°C. The assays are electrophoresed on 1.2% agarose gels (100 mL gels with double 20 lane combs) for 1 hour at 100 volts in 1X TPE buffer. Gels are stained with 0.5 μ g/mL ethidium bromide for 30 min and visualized and photographed on a UV light box.

Before start

3-4 hour XZ-6E virus infected NC64A, pellets frozen at -80°C

Protocol

Step 1.

Prepare Buffer A:



. CviRI Buffer A

CONTACT: Irina Agarkova

Step 1.1.

10 mM Tris-HCl, pH 8.0

Step 1.2. 10 mM 2-ME

Step 1.3.

50 μg/mL PMSF

Step 2.

Prepare Buffer B:



. CviRI Buffer B

CONTACT: Irina Agarkova

Step 2.1.

20 mM Tris-HCl, pH 7.8

Step 2.2.

0.5 mM EDTA

Step 2.3.

7 mM 2-ME, 10% Glycerol

Step 3.

Prepare Buffer C:

₽ PROTOCOL

. CviRI Buffer C

CONTACT: Irina Agarkova

Step 3.1.

10 mM KHPO₄, pH 7.4

Step 3.2.

0.5 mM EDTA

Step 3.3.

7 mM 2-ME, 10% Glycerol

Step 4.

Prepare Storage Buffer:



. CviRI Storage Buffer

CONTACT: Irina Agarkova

Step 4.1.

20 mM Tris-HCl, pH 8.0

Step 4.2.

0.5 mM EDTA

Step 4.3.

0.1 mM DTT

Step 4.4.

50 mM KCI

Step 4.5.

5 mM MgCl₂, 50% Glycerol

Step 5.

Prepare 1X CviRI Assay Buffer:



. 1X CviRI Assay Buffer

CONTACT: Irina Agarkova

Step 5.1.

10 mM Tris-HCl, pH 7.8

Step 5.2.

150 mM KCl

Step 5.3.

10 mM MgCl₂

Step 5.4.

 $100 \, \mu g/mL \, BSA$

Step 6.

Thaw the 3-4 hour XZ-6E virus-infected NC64A *chlorella* and suspend in MSK flasks with Buffer A. Suspend with 20 mL per flask per $1.0-1.5 \times 10^{11}$ infected cells.

Step 7.

Homogenize the cells in the MSK mechanical homogenizer with 15 gm of 0.3 mm glass beads at 4,000 rpm for 90 sec (2 X 45 sec) with CO₂ cooling.

O DURATION

00:01:30

Step 8.

Recover the homogenate to clean tubes.

Step 9.

Wash the glass beads 3X with 5 mL of Buffer A and combine with the homogenate.

Step 10.

Centrifuge the homogenate in the Sorvall SS34 rotor at 10,000 rpm,20 min, 4°C.

O DURATION

00:20:00

Step 11.

Save the supernatant.

Step 12.

Adjust the homogenate supernatant to 70% saturation with $(NH_4)_2SO_4$ at 4°C with gentle stirring. Add the $(NH_4)_2SO_4$ gradually.

Step 13.

Incubate at 4°C for 60-90 min without stirring.

O DURATION

01:30:00

Step 14.

Centrifuge the material in the Sorvall SS34 rotor at 10,000 rpm, 10 min, 4°C.

O DURATION

00:10:00

Step 15.

Save the pellet.

Step 16.

Suspend the pellets with Buffer A.

Step 17.

Per mL of suspension add: 0.45 mL of 4 M NaCl and 0.45 mL of 28% PEG 8000 (heated to 65°C).

Step 18.

Mix gently by inversion for 5-10 min.

O DURATION

00:10:00

Step 19.

Centrifuge the material in the Sorvall SS34 rotor at 10,000 rpm, 10 min, 4°C.

© DURATION

00:10:00

Step 20.

Save the supernatant.

P NOTES

Irina Agarkova 29 Mar 2016

If there is small particulate matter in the supernatant, centrifuge a second time as before and save the supernatant.

Step 21.

Dilute the supernatant with 3-5 volumes of Buffer B to reduce the NaCl concentration.

Step 22.

Load the material overnight onto a Heparin-Sepharose column equilibrated with Buffer B in the cold room.

O DURATION

18:00:00

Step 23.

Elute the Heparin-Sepharose column with Buffer B using a 0.2-1.2 M NaCl gradient.

Step 24.

Assay the column fractions and pool the active fractions.

₽ PROTOCOL

. CviRI Assay Conditions

CONTACT: Irina Agarkova

Step 24.1.

All assays are carried out in 20.0 μ L volumes with 1 μ g of pUC19 DNA as substrate for 60 to 120 min at 25°C.

O DURATION

02:00:00

Step 24.2.

The assays are electrophoresed on 1.2% agarose gels (100 mL gels with double 20 lane combs) for 1 hour at 100 volts in 1X TPE buffer.

O DURATION

01:00:00

Step 24.3.

Gels are stained with 0.5 μ g/mL ethidium bromide for 30 min and visualized and photographed on a UV light box.

O DURATION

00:30:00

Step 25.

Dilute the pooled fractions with 3 volumes of Buffer B to reduce the salt concentration.

Step 26.

Load the material overnight onto a Blue-Sepharose column (or an Affi-Gel-Blue column) equilibrated with Buffer B in the cold room.

O DURATION

18:00:00

Step 27.

Elute the Blue-Sepharose column (or Affi-Gel-Blue column) with Buffer B using a 0.2-2.0 M KOAc gradient.

Step 28.

Assay the column fractions and pool the active fractions.

₽ PROTOCOL

. CviRI Assay Conditions

CONTACT: Irina Agarkova

Step 28.1.

All assays are carried out in 20.0 μL volumes with 1 μg of pUC19 DNA as substrate for 60 to 120 min at 25°C.

O DURATION

02:00:00

Step 28.2.

The assays are electrophoresed on 1.2% agarose gels (100 mL gels with double 20 lane combs) for 1 hour at 100 volts in 1X TPE buffer.

O DURATION

01:00:00

Step 28.3.

Gels are stained with 0.5 μ g/mL ethidium bromide for 30 min and visualized and photographed on a UV light box.

O DURATION

00:30:00

Step 29.

Dilute the pooled fractions with 2 volumes of Buffer C to reduce the salt concentration.

Step 30.

Load the material overnight onto an Hydroxylapatite column equilibrated with Buffer C in the cold room.

O DURATION

18:00:00

Step 31.

Elute the Hydroxylapatite column with Buffer C using a 0-1.0 M KHPO₄ gradient.

Step 32.

Assay the column fractions and pool the active fractions.

PROTOCOL

. CviRI Assay Conditions

CONTACT: Irina Agarkova

Step 32.1.

All assays are carried out in 20.0 μL volumes with 1 μg of pUC19 DNA as substrate for 60 to 120 min at 25°C.

O DURATION

02:00:00

Step 32.2.

The assays are electrophoresed on 1.2% agarose gels (100 mL gels with double 20 lane combs) for 1 hour at 100 volts in 1X TPE buffer.

O DURATION

01:00:00

Step 32.3.

Gels are stained with 0.5 μ g/mL ethidium bromide for 30 min and visualized and photographed on a UV light box.

O DURATION

00:30:00

Step 33.

Dilute the pooled fractions with 3-4 volumes of Buffer B.

Step 34.

Load the material overnight onto a small Heparin-Sepharose column equilibrated with Buffer B.

O DURATION

18:00:00

Step 35.

Elute the column with Buffer B containing 2.0 M NaCl.

Step 36.

Collect small fractions (approximately 35 drops per fraction).

Step 37.

Assay the column fractions (only about the first 10-15 fractions) and pool the active fractions.

. CviRI Assay Conditions

CONTACT: Irina Agarkova

Step 37.1.

All assays are carried out in 20.0 μL volumes with 1 μg of pUC19 DNA as substrate for 60 to 120 min at 25°C.

O DURATION

02:00:00

Step 37.2.

The assays are electrophoresed on 1.2% agarose gels (100 mL gels with double 20 lane combs) for 1 hour at 100 volts in 1X TPE buffer.

O DURATION

01:00:00

Step 37.3.

Gels are stained with 0.5 μ g/mL ethidium bromide for 30 min and visualized and photographed on a UV light box.

© DURATION

00:30:00

Step 38.

Concentrate the pooled enzyme by dialysis overnight into storage buffer at 4°C.

O DURATION

18:00:00

Step 39.

Store the enzyme at -20°C.