

Mar 28,
2019

Working

El-cheep-o bacterial DNA isolation / miniprep

Magdalena Julkowska¹¹King Abdullah University of Science and Technology[dx.doi.org/10.17504/protocols.io.pendjde](https://doi.org/10.17504/protocols.io.pendjde)Magdalena Julkowska
King Abdullah University of Science and Technology ⚡

ABSTRACT

This is a protocol for cheep miniprep isolation. Be aware that the DNA quality might not be sufficient to use this protocol for cloning purposes. It works very well for the digestion analysis though.

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

MATERIALS

| NAME ▾ | CATALOG # ▾ | VENDOR ▾ |
|-------------------|--------------|-----------------|
| RNase A | 19101 | Qiagen |
| Tris | | |
| EDTA | | |
| Isopropanol | | |
| Potassium acetate | 1.04820.1000 | Merck Millipore |
| Sodium hydroxide | 1064981000 | Merck Millipore |
| SDS | | Ambion |
| Acetic acid | 695092 | Sigma Aldrich |

SAFETY WARNINGS

BEFORE STARTING

Prepare following buffers and keep P3 at 4C:

P1: into 250 ml H₂O

EDTA 0.093g 1mM

TRIS 1.515g 50mM(pH 8.0)

AUTOCLAVE

P2: into 250 ml H₂O

NaOH 2g 0.2N

SDS 2.5g 1%

AUTOCLAVE

P3: into 250 ml H₂O

K-acetate 73.63g

Acetic Acid 27.5ml

pH 5.5

AUTOCLAVE

- 1 Spin 3ml of the overnight culture to pellet cells (1.5 ml per time) @13000 rpm for 5 min
- 2 Remove the supernatant
- 3 Add 250 ul of P1 solution + Rnase 5ul and vortex to resuspend the cells
- 4 Incubate 15 min @ room temperature
- 5 Add 250 ul of P2 solution and mix gently
- 6 Incubate 5 min @ room temperature
- 7 Add 350 ul of P3 solution and mix gently
- 8 Leave for 10 min on ice
- 9 Spin @ the max speed for 10 min
- 10 Take 750 ul of the supernatant and put in a new eppendorf tube. Add 750 ul of isopropanol
- 11 Mix well
- 12 Spin for 30 min @ 14000 rpm
- 13 Remove the supernatant. Wash with 80% EtOH
- 14 Spin again for other 10 min
- 15 Dry the pellet @ room temperature for 10 min. Remove the remaining EtOH

16 Add 50 ul of sterile / MiliQ water

17 Measure the DNA concentration using NanoDrop and store the DNA @ -20°C



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited