



Oct 16, 2019

Harvesting OP50 for *C. elegans* liquid culture

Vidur Sabharwal¹¹Tata Institute of Fundamental Research

In Development

[dx.doi.org/10.17504/protocols.io.xcmfiu6](https://doi.org/10.17504/protocols.io.xcmfiu6)

Vidur Sabharwal ⚡ 🌈

Harvesting OP50

- 1 After 8 hr incubation, use 2 ml of *E. coli* OP50 LB culture to inoculate each of 3 x 250 ml LB. Place flasks in shaking incubator (180 rpm) O/N. ⌚ 17:00:00

🌡 37 °C

- 2 Weigh 3 x 50 ml centrifuge tubes and write weight of the empty tube on the tube.

- 3
Day 2. Using a serological pipette, aliquot 40 ml of O/N *E. coli* OP50 culture into the pre weighed centrifuge tubes. Centrifuge in swing bucket centrifuge at 3,220 x g, 10°C, 10 min. ⌚ 00:10:00

🌡 10 °C

- 4

Carefully decant the supernatant, repeat step 4 till all the culture is used up. Keep the tube inverted with lid on and leave to stand for several minutes. Using a pipette remove any excess supernatant that may have collected in the lid. Weigh the tube and calculate the weight of the pellet. ➡ [go to step #4](#)

- 5 Calculate the volume needed to provide a suspension of 30 g/L and mark this volume on the tube. Date, label and place tubes at -20 °C for use within 1-3 months or at - 80 °C if storing for longer than 3 months.
- 6 Prepare bacterial suspension for growing nematodes; allow bacterial pellet to thaw out and add required volume of S complete medium to each tube, vortex gently to resuspend pellet, pool contents of different tubes to obtain the volume required. Work under sterile conditions.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited