



Working

## Immunofluorescence Assay to detect respiratory syncytial virus [↗](#)

PLOS One

Maria Lucia Silva<sup>1</sup>, Talita Gagliardi<sup>2</sup>, Eurico Arruda<sup>1</sup><sup>1</sup>School of Medicine University of São Paulo, <sup>2</sup>University of Maryland, College Park[dx.doi.org/10.17504/protocols.io.w8ufhww](https://doi.org/10.17504/protocols.io.w8ufhww) Talita B. Gagliardi  
University of Maryland, College Park 

### ABSTRACT

IF Protocol used to detect RSV in nasopharyngeal samples.

### EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0217744>

### THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Matsuno AK, Gagliardi TB, Paula FE, Luna LKS, Jesus BLS, Stein RT, Aragon DC, Carlotti APCP, Arruda E (2019) Human coronavirus alone or in co-infection with rhinovirus C is a risk factor for severe respiratory disease and admission to the pediatric intensive care unit: A one-year study in Southeast Brazil. PLoS ONE 14(6): e0217744. doi: [10.1371/journal.pone.0217744](https://doi.org/10.1371/journal.pone.0217744)

- 1 Fix the cells immerse them into a pure cold acetone for 5min.
- 2 Wash 3X the sample in PBS for 3min each and at room temperature (r.t.).
- 3 Permeabilize the cells incubate with 0.1% Triton-X100 for 15min at r.t.
- 4 Wash 3X the sample in PBS for 3min each and at r.t.
- 5 Incubate the cells with 1:100 dilution of "anti-RSV blend monoclonal antibody" (MAB858-4; Millipore, MA, USA) in PBS with 1% Bovine Serum Albumin (BSA) and 0.2% Tween 20. Keep it inside a humidified chamber for 30min at r.t.
- 6 Wash 3X the sample in PBS for 5min each and at r.t.
- 7 Incubate the sample for 15min at r.t in n a Block solution (10% of specific serum diluted in PBS 1X).

Wash 3X the sample in PBS for 3min each and at r.t.

8

9 Incubate the cells with the secondary antibody solution for 30min at r.t. and inside a humidified chamber. It was used "Donkey anti-mouse IgG (H+L) secondary antibody conjugated with Alexa Fluor® 488" (Life Technologies, Thermo Scientific, Carlsbad, CA, USA) diluted 1:200 in PBS with 1% BSA.

10 Wash 3X the sample in PBS for 5min each and at r.t.

11 Incubate the cells with 5µg/ml of DAPI solution (diluted in PBS) for 5min at r.t. and inside a humidified chamber.

12 Wash 2X the sample in PBS for 5min each and at r.t.

13 Mount coverslip with a drop of mounting medium (90% glycerol diluted in PBS).

14 Seal coverslip with nail polish to prevent drying and movement under microscope.

15 Store in dark at -20°C or +4°C.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited