

# Flag M2 immunoprecipitation and Western Blot in *Drosophila*

Caroline Grabbe, Ingrid Dacklin, Behzad Khoshnood

## Abstract

**Citation:** Caroline Grabbe, Ingrid Dacklin, Behzad Khoshnood Flag M2 immunoprecipitation and Western Blot in *Drosophila*. **protocols.io**

dx.doi.org/10.17504/protocols.io.jivcke6

**Published:** 24 Aug 2017

## Protocol

### Step 1.

#### Immunoprecipitation with Flag M2 magnetic beads in *Drosophila*

1. Grow flies expressing transgenic constructs under control of the UAS/GAL4 system at 25°
2. Procedure to prepare protein lysate from different developmental stages in *Drosophila*;

##### 1. Embryo collection

- To retrieve a collection of all embryonic stages, harvest 0-22 hour old embryos in a mesh-basket and wash away the yeast paste with deionized H<sub>2</sub>O
- Dechorinate the embryos in a 1:1 dilution of Sodium hypochlorite 14% (f.c. 7%) and deionized H<sub>2</sub>O for 2-4 minutes, rinse embryos for 1 min in deionized water.
- Transfer the embryos in 100 µl portions to 1.5 ml microtubes and mince with pestle in 200 µl Lysis buffer A (50 mM HEPES pH7.4, 150 mM NaCl, 1 mM EDTA, 1 mM EGTA, 1% Triton X-100, 10% Glycerol, 25 mM NAF, 10 µM ZnCl<sub>2</sub>), supplemented with 1 tablet Complete Inhibitors (Roche) per 50 ml, 1mM PMSF and 10 mM N-ethyl maleimide (NEM).
- Add 800 µl Lysis buffer A to each embryo-portion and continue grinding.
- Incubate all samples on ice for 30 minutes.
- Spin 14 000 x g for 15 min at 4°C and transfer and pool supernatants to fresh tubes.
- Measure protein concentration.

##### 1. Larvae collection

- Collect wandering third instar *Drosophila* larva and rinse them twice in 1xPBS.
  - Add 25 larva per aliquot to 1.5 ml microtubes with 125 µl Lysis buffer A (see above) and grind the larvae with a pestle.
  - Add 375 µl fresh Lysis buffer A and continue grinding.
  - Incubate all samples on ice for 30 minutes.
  - Spin 14 000 x g for 15 min at 4°C and transfer and pool supernatants to fresh tubes.
  - Measure protein concentration.
- 
- Adult fly collection
  - Collect 0-4 day old adults flies under CO<sub>2</sub> anesthesia, 40 flies per aliquot, in 1.5 ml microtubes.
  - Grind the flies in 200 µl Lysis buffer A (see above), add another 200 µl and continue grinding.
  - Incubate all samples on ice for 30 minutes.
  - Spin 14 000 x g for 15 min at 4°C and transfer and pool supernatants to fresh tubes.
  - Measure protein concentration.
- 
3. Take out 60 µl from the lysate, mix with 60 µl 2x Sample buffer and incubate 7 min at 98°C (=INPUT).
  4. Wash the Flag M2 magnetic beads (Sigma M8823) 3 times in Lysis Buffer A before usage. Make sure that the beads are completely defrosted before aliquoting the slurry with a cut pipet tip. Collect magnetic beads with a magnetic separator.
  5. Set up immunoprecipitations with Flag M2 magnetic beads at 4°C, by mixing 1-2 mg total protein with 30 µl Flag M2 magnetic beads (slurry, equals 15 µl packed beads) in 1.5 ml microtubes.
  6. Incubate for 2 hours at 4°C on rotator. Collect beads with magnetic separator.
  7. Wash 6 times with 1 ml Lysis buffer A (collect beads with a magnetic separator between each wash).
  8. Add 25 µl 2X Sample buffer with max 5% β-Mercaptoethanol to the magnetic beads. Incubate the samples at 98°C for 4 min (alternatively the elution can be performed with 0.1 M Glycine-HCl, pH 3.0).

9. Separate the samples for loading using a magnetic separator and transfer the eluted proteins to fresh tubes with a narrow ended tip.
10. Analyze the samples by separating the immunoprecipitated proteins on SDS-PAGE, followed by Western Blot.

### **Western Blot in *Drosophila***

1. Separate the eluted immunoprecipitates and protein lysates by SDS-PAGE (use a percentage of acrylamide suitable to detect your protein of interest or a gradient gel).
2. Transfer the separated proteins onto a 0.2 mm PVDF membrane (Merck Millipore) using a semi-dry Trans-Blot Turbo transfer system (Bio-Rad).
3. Block the membranes in 5 % BSA in 1X TBS (50 mM Tris-HCl, pH 7.5, 150 mM NaCl) for 1 hour at room temperature.
4. Incubate with primary antibodies, diluted in 1X TBS with 5 % BSA, overnight at 4°
5. Wash the membranes 5 times for 10 min in TBST (1X TBS with 0.075% Tween-20).
6. Incubate with HRP-linked secondary antibodies (GE Healthcare) for 1 hour at room temperature.

7. Wash the membranes 5 times 10 min in TBST at room temperature, and finally 10 minutes with 1X TBS.
8. Develop the membranes using standard techniques and according to the manufacturers recommendations; (i) SuperSignal™ West Pico Chemiluminescent Substrate (Thermo Scientific), ii) ECL™ Prime Western blotting Detection reagent (GE Healthcare).
9. Visualize by standard autoradiography on film or with a C-Digit Blot Scanner analysis (LI-COR Biosciences).