

Jul 17, 2019

Calibrating 10-AU fluorometer to quantify extracted chl-a

Shelby Whitehead¹, Samantha Coy¹, Steven Wilhelm¹¹The University of Tennessee, Knoxville

1

Works for me

dx.doi.org/10.17504/protocols.io.2ubgesn

The Aquatic Microbial Ecology Research Group - AMERG (The Buchan, Zinser and Wilhelm labs)



Shelby Whitehead ⚡

ABSTRACT

This method is for use for calibrating the 10-AU fluorometer using chl-a extracts. The method also provides guidelines for extracting chl-a and calculating absorbance using a spectrophotometer, making dilutions of extracts that can be used to create different chl-a concentrations

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

S.W. Jeffrey, G.F. Humphrey, New spectrophotometric equations for determining chlorophylls a, b, c1 and c2 in higher plants, algae and natural phytoplankton, *Biochemie und Physiologie der Pflanzen*, Volume 167, Issue 2, 1975, Pages 191-194, ISSN 0015-3796, [https://doi.org/10.1016/S0015-3796\(17\)30778-3](https://doi.org/10.1016/S0015-3796(17)30778-3).



Shelby's Notes for Fluorometer.pdf



Sam's Notes for Fluorometer.pdf

MATERIALS TEXT

97% acetone

spinach or green plant as source for chlorophyll extracts

SAFETY WARNINGS

97% acetone is flammable.

See MSDS safety and hazard warnings.



BEFORE STARTING

- *Fixed Application* can be found under "general tests" on the spectrophotometer.
- Make sure the cuvette is in the correct position for measurements. Misplacement may cause an error such as a negative absorbance measurement.
- Spinach only has chlorophyll a. If using an organism that has more pigments, will have to subtract those from total absorbance (i.e. chlorophyll a - chlorophyll c).

Chlorophyll Extraction

- 1 Extract chlorophyll from 1 g of ground spinach with 5 ml of 90% acetone at -20 °C overnight.

Spectrophotometry

- 2 Use the Biomate 5 spectrophotometer to measure the adsorption of extracted chlorophyll by adding  1 ml to a quartz cuvette.
 - Make sure the cuvette is in the correct position in the machine.
 - Zero the machine by blanking with  1 ml of 90% acetone.
 - Measure "multi λ " with $\lambda_1 = 664 \text{ nm}$ and $\lambda_2 = 647 \text{ nm}$. (Jeffrey 1975)
- 3 Based on the absorbance value from [go to step #2](#) : dilute the chlorophyll extract so that the absorbance is in optimum range for the fluorometer. Optimum range is Abs λ_1 and $\lambda_2 = 0.3-0.85$.

Example:

Undiluted

Abs $\lambda_1 = 1.88$

Abs $\lambda_2 = 1.324$

20-fold dilution:

Abs $\lambda_1 = 0.736$

Abs $\lambda_2 = 0.479$

Fluorometry

- 4 Calculate the concentration of diluted chlorophyll extract using the absorbance values from the spectrophotometer and the equation:
concentration of chlorophyll a = $11.93(\lambda_1) - 1.93(\lambda_2)$

(Jeffrey 1975)



If using an organism that has more pigments, will have to subtract those from total absorbance (i.e. chlorophyll a - chlorophyll c). The spinach we used did not require subtraction of chlorophyll c.

Example:

$$11.93(0.736) - 1.93(0.479) = 7.85601 \text{ ug/ml}$$

Calibrate Fluorometer

- 5 Convert answer from [go to step #4](#) to ug/L, because this is the readout for the fluorometer.
- 6 Turn on the 10-AU Fluorometer and allow the instrument to warm up for at least 30 minutes.

- 7 Create solutions of varying concentrations to test the fluorometer readings. Use $C1V1=C2V2$ to make solutions at different concentrations based on expected field concentrations, for example: 20, 10, 5, 1 and 0.5 ug/L. Also, create an acetone blank (with no chlorophyll).



For example:

7,856 ug/L (C1)

For a 5 ug/L concentration (C2) in 7 mls acetone (V2)

$$(7,856 \text{ ug/L})(V1) = (5 \text{ ug/L})(7 \text{ ml})$$

$$V1 = 4.45519 \text{ ug}$$

- 8 Use the concentration most expected in the field to set as the solid standard.

Chlorophyll (ug/l)	Rhodamine WT (ppb)	Range	% FS (+/-5)
240 - 180	100 - 75	High	80%
180 - 160	75 - 65	High	70%
160 - 130	65 - 55	High	60%
130 - 110	55 - 45	High	50%
110 - 80	45 - 35	High	40%
80 - 60	35 - 25	High	30%
60 - 40	25 - 15	High	20%
24 - 18	10 - 7.5	Medium	80%
18 - 16	7.5 - 6.5	Medium	70%
16 - 13	6.5 - 5.5	Medium	60%
13 - 11	5.5 - 4.5	Medium	50%
11 - 8.0	4.5 - 3.5	Medium	40%
8.0 - 6.0	3.5 - 2.5	Medium	30%
6.0 - 4.0	2.5 - 1.5	Medium	20%
2.4 - 1.8	1.0 - 0.75	Low	80%
1.8 - 1.6	0.75 - 0.65	Low	70%
1.6 - 1.3	0.65 - 0.55	Low	60%
1.3 - 1.1	0.55 - 0.45	Low	50%
1.1 - 0.8	0.45 - 0.35	Low	40%
0.8 - 0.6	0.35 - 0.25	Low	30%
0.6 - 0.4	0.25 - 0.15	Low	20%

Table 1 is an example of how to set the scale for a given linear range of the 10-AU.

For this example we are using a 20µg/L standard.

Table 1 is an example of how to set the sensitivity scale for the 10-AU. Use the table to determine the appropriate range and %FS for the calibration standard.

- 9 Clear past calibrations on screen #2.6.
- 10 Access screen #2.43 and set range control to "manual" and range to the appropriate level for the standard chosen. For example, we set range control to "Man" and "Med."
- 11 Access screen #3.2. Loosen the sensitivity locking screw using an allen wrench. (The sensitivity lock is the small screw to the left of the number pad.)
- 12 Insert and cover your selected standard concentration. Use a coin or tool to slowly turn the sensitivity knob (located to the right of the number pad, under the power button) to adjust the %FS to the value appropriate for the standard concentration according to Table 1 in [go to step #8](#). A deviation of 5% less than or greater than the table is acceptable. When %FS reaches the designated value, tighten the sensitivity lock.

At this point, sensitivity is set.



Example:
For 5 ug/L.
FS = 21%

13 Access screen #2.43 again, and set range control to "auto."

14 Access screen #2.1. Insert and cover the acetone blank created in [go to step #7](#). Press "1" to Run Blank. Wait for the reading to stabilize for about 8 seconds, and then press "0." The instrument will take 15 seconds to complete blanking.

At this point, blank is set.

15 Calibrate with your standard:

15.1 Access screen #2.2 and input the actual concentration of your standard. Return to screen #2.0

15.2 Insert and cover the standard.

15.3 Run standard by pressing "3." Wait about 8 seconds for reading to stabilize.

15.4 Press the "*" to set the calibration point. After about 15 seconds, a message will appear that says "Finished".

16 Test the 10-AU fluorometer using the concentrations created in [go to step #7](#) to assure appropriate readout.



Example:

Expected	Observed
5 ug/L	5.01 ug/L
1 ug/L	1.11 ug/L
Blank	0.0 ug/L



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited