



Jun 12,  
2019

## RNAi Plasmid Construction using pFGC5941 [↗](#)

Yaowu Yuan<sup>1</sup>

<sup>1</sup>University of Connecticut

Working

[dx.doi.org/10.17504/protocols.io.2w2gfge](https://doi.org/10.17504/protocols.io.2w2gfge)

Mimulus



Andrea Sweigart

### EXTERNAL LINK

[http://mimubase.org/FTP/Protocols/Plasmid\\_Construction/RNAi%20plasmid%20construction%20using%20pFGC5941.pdf](http://mimubase.org/FTP/Protocols/Plasmid_Construction/RNAi%20plasmid%20construction%20using%20pFGC5941.pdf)



RNAi plasmid construction  
using pFGC5941.pdf

### GUIDELINES

This protocol is based on the vector pFGC5941 (ABRC Stock CD3-447).

To avoid off-target effect, make sure no other regions in the interested genome perfectly match the RNAi fragment (150-500 bp) for a contiguous block longer than 16 bp. Also, make sure there are no restriction sites for the enzymes NcoI, Ascl, BamHI, or XbaI within the RNAi fragment.

When designing primers to amplify the RNAi fragment. Add "GTTCTAGACCATGG" at the 5' end of the Forward primer and add "GTGGATCCGGCGCGCC" at the 5' end of the Reverse primer.

Make sure you have digested the pFGC5941 vector using NcoI/Ascl before the first ligation.

### Primer sequences:

pFGC5941\_2372F: CTTTCATCGAAAGGACAGTAGAA

pFGC5941\_3082R: CCAAACAGGCTCATAGATACT

pFGC5941\_3930F: TGTACATCAGAATGTTTCTGAC


pFGC5941\_4430R: CGCTCTATCATAGATGTCGCTA

### SAFETY WARNINGS

For Safety Warnings and Hazard Information please refer to the SDS (Safety Data Sheet).

Amplifying insert from cDNA or gDNA using Phusion PCR

## 1 Amplify insert from cDNA or gDNA (if the fragment contains no intron) using Phusion PCR

Make **TWO**  **20 µl** reactions of the following in separate tubes:

Amount (µL)	Reagent
4 µL	5x Phusion Buffer
0.5 µL	10 mM dNTPs
0.6 µL	DMSO
1.0 µL	Template
0.2 µl	Phusion enzyme
11.0 µL	dH2O
1.5 µL	5 µM Forward Primer
1.5 µL	5 µM Reverse Primer
<b>20 µL</b>	<b>Total</b>

## 2 Run Phusion PCR program:

Cycle	Repeats	Temperature	Time
Cycle 1		98°C	0:30
Cycle 2	(32x)	98°C	0:10
		58°C (or the ideal annealing temperature)	0:20
		72°C	0:30
Cycle 3		72°C	5:00
Cycle 4		12°C	for ever

### Digestion

## 3 Digest one insert with **NcoI/Ascl** and the other one insert **BamHI/XbaI**.



See step 12 for BamHI/XbaI digestion.

Amount (µL)	Reagent
2.5 µL	10x CutSmart Buffer
4.5 µL	dH2O
1.5 µL	NcoI
1.5 µL	Ascl
15 µL	PCR Product
<b>25 µL</b>	<b>Total</b>

### 3.1 Incubate samples for **01:00:00** at **37 °C**.

### 3.2 Gel purify digests and save the BamHI/XbaI digested insert for the second ligation.

### First ligation

## 4 First Ligation (Would like insert to vector molar ratio to 2:1 to 6:1)

Amount (μL)	Component
2 μL	Linearized pFGC5941 digested with <i>AscI</i> / <i>NcoI</i> . (~175ng; adjust volume as needed)
4 μL	Insert (digested with <i>AscI</i> / <i>NcoI</i> ) (~15-30ng)
2 μL	T4 Ligase Buffer
1 μL	T4 Ligase
11 μL	dH <sub>2</sub> O
<b>20 μL</b>	<b>Total</b>

4.1 Incubate at  **Room temperature** for  **00:30:00**.

4.2 Transform  **10 μl** into *E. coli* competent cells (homemade) and plate on Kan plates.

Colony PCR to check for first insert

## 5 Colony PCR to check for first insert

Amount (μL)	Component
8.0 μL	dH <sub>2</sub> O
1.0 μL	10x buffer
0.125 μL	dNTPs
0.5 μL	pFGC5941 2372 F
0.5 μL	pFGC5941 3082 R
0.05 μL	Taq
<b>10 μL</b>	<b>Total</b>

5.1 Run Colony PCR

Cycle	Repeats	Temperature	Time
Cycle 1		95°C	3:00
Cycle 2	32x	95°C	0:15
		55°C	0:15
		72°C	1:00
Cycle 3		72°C	7:00
Cycle 4		12°C	forever

6 Circle the biggest colonies on your plate and label them 1-8.

7 Make a replica plate for your colonies.

8 PCR across the first insert using primers on the vector to check for an insert:  
**An empty vector will give a band of 700bp**

Picking Colonies and Plasmid Prep

9 Pick two correct colonies and inoculate into  **3 ml LB+Kan broth**.

10 Incubate in  **37 °C** shaker overnight.

11 The next day, do a plasmid prep (mini-prep kit) with 1 of the colonies that grew well.

#### Digest Plasmid with BamHI/XbaI

12

Amount	Component
5 µl	10x CutSmart Buffer
12 µl	dH2O
1.5 µl	XbaI
1.5 µl	BamHI
30 µl	Plasmid*
<b>50 µl</b>	<b>Total</b>

\* adjust volume based on concentration; you want 2000-5000 ng of plasmid

12.1  **37 °C** for  **01:00:00**.

12.2 Gel purify digest.

#### Ligation #2

13  **2 µl** vector that contains the first insert, digested with BamHI/XbaI (~**175 ng**; **adjust volume based on concentration**).

Amount	Component
4 µl	insert digested with BamHI/XbaI (done in step 3) (want ~15-30 ng)
2 µl	T4 ligase buffer
1 µl	T4 ligase
11 µl	dH2O
<b>20 µl</b>	<b>Total</b>


13.1 Incubate for  **00:30:00** at  **Room temperature**.

13.2 Transform  **10 µl** into *E. coli* competent cells (homemade) and plate on Kan plates.

#### Colony PCR to check for second insert

14 pFGC5941 **3930 F** & pFGC5941 **4430 R**

Vector without insert will give a band of 500bp

15 Pick two correct colonies and inoculate into  3 ml LB+Kan broth.

15.1 Incubate in  37 °C shaker overnight.

15.2 Plasmid prep (mini-prep kit)

Check plasmid for inserts

16 PCR to check for both inserts:  
2372F/3082R or RNAi\_R (insert specific)  
3930F/4430R or RNAi\_F (insert specific)

Sequence to verify

17 Use 4 primers:  
2372F, 3082R, 3930F, 4430R



**Note:** in the sequencing reaction, add DMSO to aid in the sequencing across the restriction enzyme digest sites (the chromatogram peaks usually drop off dramatically right after the digest sites; an alternative strategy is to PCR the final plasmid with 2372F&3082R for the left insert and 3930F&4430R for the right insert and then sequence the PCR product)

18 Transform into agrobacterium for infiltration.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited