

# Non-Centrifuged Plaque Assay Sample Steps

**Bonnie Poulos**

## Abstract

For [Transcriptomics During One-Step Growth Curves for Cellulophaga Phages protocol](#).

**Citation:** Bonnie Poulos Non-Centrifuged Plaque Assay Sample Steps. **protocols.io**

dx.doi.org/10.17504/protocols.io.ek2bcye

**Published:** 23 Feb 2016

## Protocol

### Step 1.

Pipet 10  $\mu$ l from the flask into 990  $\mu$ l of MSM in a 1.5 ml tube

#### NOTES

**Bonnie Poulos** 23 Feb 2016

You are diluting your sample 100x:  $10^{-2}$ .

### Step 2.

Vortex briefly

### Step 3.

Prepare  $10^{-4}$  and  $10^{-5}$  serial dilutions. Pipet 10  $\mu$ l of the  $10^{-2}$  dilution made from the growth flask into 990  $\mu$ l MSM for the  $10^{-4}$  dilution. Pipet 100  $\mu$ l from the  $10^{-4}$  dilution into 900  $\mu$ l MSM for the  $10^{-5}$  dilution.

### Step 4.

Plate 100  $\mu$ l each from the  $10^{-4}$  and  $10^{-5}$  dilutions in duplicate or triplicate.