

# Gateway BP recombination of attB tailed PCR products into pDONR/zeo (5 μl assay) Version 2

#### **Johannes Debler**

## **Abstract**

This is a slightly modified version of Thermo Fisher's gateway BP protocol which uses less enzyme and is therefore more economical.

**Citation:** Johannes Debler Gateway BP recombination of attB tailed PCR products into pDONR/zeo (5 μl assay).

protocols.io

dx.doi.org/10.17504/protocols.io.g5rby56

Published: 07 Feb 2017

## **Protocol**

## Calculate amount of attB PCR product needed

#### Step 1.

ng of attB PCR product = size of attB tailed PCR product in bp \* 0.0165

eg. 778 bp \* 0.0165 = 12.9 ng

If your purified PCR product has a concentration of 25 ng/ul you would use 0.5 μl in the next step.

#### NOTES

Johannes Debler 07 Feb 2017

Your attB PCR product should be purified either by gel purification or with a PCR clean-up kit

## Combine in 0.2 ml PCR tube

## Step 2.

attB PCR product	x µl (as calculated in step 1)
pDONR/zeo (150 ng/μl)	0.5 μΙ
TE buffer pH 8.0	to a total volume of 4 μl

#### BP Clonase II enzyme mix

## Step 3.

Remove from freezer and vortext for 2 seconds.

#### add BP Clonase II

Step 4.

Add **0.5** µI of BP Clonase II enzyme mix and mix well.

## Incubate at room temperature

Step 5.

Incubate for **1 hour** at **room temperature** (2-3 hours or overnight for more colonies).

**O DURATION** 

01:00:00

## Stop reaction with Proteinase K

Step 6.

Add **0.5** µl Proteinase K, mix well and incubate for **10 min at 37°C.** 

**O** DURATION

00:10:00

## Transform into Omnimax 2 E.coli

Step 7.

Combine **1**  $\mu$ I BP product with **50**  $\mu$ I Omnimax 2 competent *E.coli* cells in 1.5 ml tube.

### Incubate on ice

Step 8.

**O DURATION** 

00:10:00

Heat shock

Step 9.

30 seconds at 42°C

## Reco on ice

**Step 10.** 

**O DURATION** 

00:02:00

# Add SOC

Step 11.

Add **500 μl** SOC.

#### Incubate

**Step 12.** 

# 1 hour at 37° C at 250 rpm on shaking incubator.

**O DURATION** 

01:00:00

Plate out

**Step 13.** 

Plate 20 ul (or more) on LB plates containing 50 ug/ml zeocin and incubate at 37° C over night.

## **P** NOTES

Johannes Debler 07 Feb 2017

I have observed that the older your BP Clonase II enzyme mix is, the more of the BP reaction you need to plate out. You can therefore concentrate your cells by centrifuging at  $1000 \times g$  for a couple of minutes, decant most of the supernatant and resuspend the cell pellet in about  $50-100 \mu l$  of the leftover supernatant. Then plate all of that on the LB + zeocin plate.