

# AQUA cloning

Anna Behle, Nicolas Schmelling

## Abstract

Protocol for seamless, scarless, homology-based cloning of vectors using an enzyme-free technique.

The method was published by Zurbriggen et al., 2015, doi: [10.1371/journal.pone.0137652](https://doi.org/10.1371/journal.pone.0137652)

It relies on the endogenous homologous recombination machinery of *E. coli* lab strains.

**Citation:** Anna Behle, Nicolas Schmelling AQUA cloning. **protocols.io**

[dx.doi.org/10.17504/protocols.io.gs9bwh6](https://dx.doi.org/10.17504/protocols.io.gs9bwh6)

**Published:** 19 Dec 2016

## Before start

Make sure to recover high concentrations of fragments.

## Protocol

### Materials and equipment

#### Step 1.

Material	Volume/Mass	Company	Serial no.
Q5 Polymerase	Dependent on reaction	NEB	M0491L
Nucleo-Spin Gel & PCR clean up	250 reactions	Macherey-Nagel	740609.250
NEB5a competent cells	25 µL	NEB	C2987I
Sterile reaction tubes	1.5 mL		

### Primer design:

#### Step 2.




PROTOCOL

## **Primer design for assembly-based cloning**

CONTACT: [Anna Behle](#)

Experimental design/overview

### **Step 2.1.**



Black arrows: Annealing part of primer.

Dotted lines: overhang part of primer (color coded).

### **Step 2.2.**

Create in silico vector maps for template and target vectors using SnapGene or a similar tool.

Annealing part of primer:

### **Step 2.3.**

- 18-21 bp in length, 40-60 % GC-content,  $T_m = 55-65$  °C.  $T_m$ - values of two primer pairs (annealing part) should not differ by more than 3 °C.
- Single nucleotide repeats of four or more should be avoided when possible.
- A 3'-GC-clamp (1 or more G's or C's at the 3'-end of the primer) can be beneficial for annealing.

Overhang part of primer:

### **Step 2.4.**

- 16-32 bp in length - the longer, the better! 40-60 % GC-content (total primer)
- For  $T_m > 72$  °C, an annealing temperature of 72 °C can be used.

Quality control

### **Step 2.5.**

- Check primers for  $T_m$ , hairpin structures and primer dimers using online tool of choice (but stick with the same one after that, e.g. <http://eu.idtdna.com/calc/analyzer>)

- Hairpin structures: Check the  $T_m$  of the structure. It should be lower than the annealing temperature used in the PCR reaction.
- Primer dimers: Recommended  $\Delta G$  is -9 kcal/mol or more positive values. However, for primers containing palindromic sequences, such as restriction sites, a more negative  $\Delta G$  is not uncommon. These primers can still work.

#### Quality control

##### Step 2.6.

Before ordering, make sure your primers **bind in the template** you plan on using.

#### Quality control

##### Step 2.7.

After PCR, fragments should be separated using gel-electrophoresis and extracted using a kit of choice. This minimizes contamination with template backbones that could yield false-positive clones.

#### Preparation of PCR fragments

##### Step 3.

Create PCR products of all the fragments that are to be assembled, using the appropriate primers. Homologous overhangs should range between 16 and 32 bp.

#### Optional: Preparation of plasmid backbone

##### Step 4.

Optional when using pSHDY or other large backbone: Digest backbone using appropriate restriction enzymes. (*Note: AQUA assembly using digested backbone is less efficient.*)

#### Recovery of fragments

##### Step 5.

Gel-extract PCR and backbone fragments using the Nucleo-Spin kit. Make sure to recover high amounts. (*Note: PCR clean-up alone is insufficient due to backbone templates yielding false positives.*)

#### Assembly reaction

##### Step 6.

In a 10 to 15  $\mu\text{L}$  reaction, mix fragments with molar ratios of 3:1 (insert:vector) or more for 2-fragment assembly or equimolar ratios for 3 or more fragments.

#### ■ ANNOTATIONS

**Dennis Dienst** 10 Oct 2017

How about the absolute amount of DNA? According to Beyer et al. - '12 ng of linearized vector per 1 kb'. Further experiences?

In my hands, assembly of divergent fragments (1500bp: 2700bp: 4700 bp) works well with gradual/linear ratios (e.g. 1:2:4), based on that 12ng-of-backbone-suggestion

### Step 7.

Incubate reaction at RT for 1 hour.

 DURATION


01:00:00

### Step 8.

Transform appropriate reaction volume into chemically competent DH5α cells. (*Note: Best results are achieved using commercially available competent .*)



REAGENTS

 NEB 5-alpha Competent E.coli (High Efficiency) - 20x0.05 ml [C2987H](#) by [New England Biolabs](#)

### Step 9.

Centrifuge transformation mix at 2500 rpm for 2 minutes. Remove supernatant; resuspend pellet in 100 μL LB.

Plate entire amount of cells on LB-plate containing appropriate antibiotics. Incubate over night at 37 °C.

 DURATION

00:02:00

Control

### Step 10.

Test clones using restriction analysis or colony PCR. (*Note: For colony PCR, use primers that amplify a region spanning two assembled fragments to avoid false positives.*)