



Sep 08, 2019

## Generating chemically competent (*E. coli*) cells V.3

Su Yin Phua<sup>1</sup>, Nay Chi Khin<sup>1</sup>, Estee Tee<sup>1</sup>, [Diep R. Ganguly](#)<sup>1</sup>, Kai Xun Chan<sup>1</sup>, Veronica Albrecht-Borth<sup>1</sup>, Barry J Pogson<sup>1</sup>

<sup>1</sup>The Australian National University

**1** Works for me [dx.doi.org/10.17504/protocols.io.666hhhe](https://doi.org/10.17504/protocols.io.666hhhe)

[Pogson Genomics Group](#)



**Diep R. Ganguly**  
The Australian National University



### ABSTRACT

Protocol for generating chemically competent (ie. for heat shock transformation) *E. coli* cells for transformation (e.g. heterologous protein expression, cloning).

### GUIDELINES

Make sure to use aseptic technique to avoid contamination. Check for any growth conditions required for *E. coli* strain being grown. Prepare autoclaved 2x1 L Luria broth (25 g/L). Solutions A and B can be made prior to starting (ie. with overnight culture) or can be done in parallel with bacterial growth with Mg<sup>2+</sup>. Make sure labelled eppies and -80°C storage box is organised beforehand.

### MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
MOPS	<a href="#">View</a>	P212121
Luria-Bertani (LB) broth, makes 1L	K488	<a href="#">Amresco</a>
Magnesium Chloride	AC223210010	<a href="#">Fisher Scientific</a>
37°C Incubator		
100 Assay Cuvettes, 1.0ml	786-009A	<a href="#">G-Biosciences</a>
Refrigerated centrifuge for conical tubes 50 ml		
Glycerol	GB0232.SIZE.500ml	<a href="#">Bio Basic Inc.</a>
BL21(DE3) or BL21-Star(DE3) or Rosetta2(DE3) or etc for protein purification		
Falcon® Conical Tubes, 50 mL 500 Tubes	38010	<a href="#">Stemcell Technologies</a>
8-10 sterile 50 mL conical tubes		
Sodium hydroxide	1064981000	<a href="#">Merck Millipore</a>
EMD Millipore™ Stericup™ Sterile Vacuum Filter Units	SCGPU05RE	<a href="#">Fisher Scientific</a>
UV/Vis spectrophotometer	<a href="#">View</a>	

### BEFORE STARTING

Plate out desired *E. coli* strain (e.g. DH5α, BL21 star, Top10) to have colonies for an initial culture.

## Prepare solutions A and B

- 1 In a beaker, prepare 250 ml of Solution A: 1m
  - 100 mM  $\text{MgCl}_2$  (hexahydrate = 5.08 g; anhydrous = 2.38 g)
  - 10 mM MOPS (0.52g)
  - Up to ~250 mL  $\text{H}_2\text{O}$Adjust to pH 7 using 1 M NaOH (slowly; ~900  $\mu\text{L}$ )  
Transfer to measuring cylinder and ensure final volume = 250 mL
- 2 In a beaker, prepare 300 ml of Solution B 1m
  - 10 mM  $\text{MgCl}_2$  (hexahydrate = 0.61 g; anhydrous = 0.29 g)
  - 10 mM MOPS (6.78g)
  - 50 mM  $\text{CaCl}_2$  (dihydrate = 2.21g)
  - Fill to ~300 ml with  $\text{H}_2\text{O}$Adjust to pH 6.5 with 1 M NaOH (~6.8 ml)  
Transfer to measuring cylinder and ensure final volume = 300 mL
- 3 Filter sterilize (stericup) both solutions and store @ 4°C. MOPS degrades into unknown yellow byproduct from autoclaving! 1m


## Generate competent cells

- 4 Culture bacteria at 37°C overnight in 5 mL LB media 1h
- 5 Inoculate desired of LB media using an aliquot, or entire 5 mL, overnight culture. Volumes can range from 50 - 500 mL depending on a variety of factors (e.g. lab time required, equipment available, amount of cell desired). Ensure a proper vessel is used to allow efficient aeration of cultures. 1m

e.g. Add 5 mL starter culture to 500 mL LB media (use 1-2 L connical flask) for a large batch of competent cells. Or, add 500  $\mu\text{L}$  - 1 mL starter culture to 50 mL LB media (use 250 mL connical flask).

Add autoclaved  $\text{MgSO}_4$  to final concentration 10 mM.
- 6 Grow to  $\text{OD}_{600} = 0.4 - 0.5$  (do not exceed 0.6). This should take between 1 - 4 hrs depending on volumne of starter culture used (ie. 5 mL starter culture ~ 1 hr to reach  $\text{OD}_{600} \sim 0.5$ ). 1m
- 7 Cool on ice for 10 mins (make sure centrifuge is refridgerated and ready) 1m
- 8 Pellet at 4°C for 10 mins @ ~4200 rcf 1m
- 9 Remove supernatant and resuspend cells in cold solution A (volume can be adjusted to minimize the volume required for centrifugation but to adequately resuspend colonies e.g. 250 mL for 500 mL culture or 25 mL for 50 mL culture). 1m
- 10 Split resuspended cells equally across centrifuge bottles or falcon tubes. Incubate 10 - 20 minutes at 4°C. 2m

- |    |   |    |
|----|---|----|
| 11 | Pellet at 4°C for 10 mins @ ~4200 rcf   | 1m |
| 12 | Remove supernatant and resuspend cells in cold solution B (repeat as per solution A).   | 1m |
| 13 | Incubate 30 minutes at 4°C.   | 3m |
| 14 | Pellet at 4°C for 10 mins @ ~4200 rcf   | 1m |
| 15 | Remove supernatant and resuspend cells in 4.3 mL of cold Solution B + 700 µL glycerol (pipette slowly).   | 1m |
| 16 | Make aliquots of 20, 50, and/or 200 µL (variable volumes based on plasmid to be used). Freeze in liquid N <sub>2</sub> before storing at -80°C. | 4m |

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited