

# Sugar content in haemolymph

Christine Damrau, Brembs' lab members

## Abstract

This is a very simple protocol showing how to extract haemolymph from adult *Drosophila melanogaster*. (Based on protocols from Sigma Aldrich).

**Citation:** Christine Damrau, Brembs' lab members Sugar content in haemolymph. **protocols.io**

dx.doi.org/10.17504/protocols.io.c8bzsm

**Published:** 28 Jul 2015

## Guidelines

### Materials:

- 0.5 µl capillaries with adequate suction cup
- 0.5ml and 1.5ml Eppendorf cap
- Ice
- Peaked stylus
- 2 forceps (or scissors) to clip wings

## Protocol

### Fly preparation

#### Step 1.

Sting 3 holes in a 0.5ml Eppendorf cap and put into 15ml Eppendorf cap with removed lid.

#### 🔗 NOTES

**Björn Brembs** 22 Jul 2015

Size of the holes is important: not too small so that enough haemolymph can go through, not too big so that no fat or other dirt is going through.

### Fly preparation

#### Step 2.

Remove the flies' wings, spear the fly's thorax with the peaked stylus.

### Fly preparation

#### Step 3.

Collect 20 speared flies in the 0.5ml Eppendorf cap with holes, on ice.

### Fly preparation

#### Step 4.

Centrifuge the 0.5ml Eppendorf cap within the 1.5ml one (1 min, 5000 rpm, at 4°C).

#### 🕒 DURATION

00:01:00

### Fly preparation

#### Step 5.

Discard the 0.5ml Eppendorf cap, soak the pellet with a capillary.

#### Fly preparation

##### Step 6.

Record the amount of soaked haemolymph (to fill up the 0.5µl you need around 50 flies).

#### Fly preparation

##### Step 7.

The haemolymph from the capillary can be transferred with the suction cup anywhere.

#### Enzymatic procedure

##### Step 8.

Add 19.5µl cold PBS to 0.5µl haemolymph.

#### Enzymatic procedure

##### Step 9.

Add 10µl of this mixture to 30µl Citrate Acid Buffer and 10µl of a 3% Trehalase-Citrate acid buffer solution.

#### Enzymatic procedure

##### Step 10.

Incubate over night at 37°C.

 DURATION

18:00:00

#### Enzymatic procedure

##### Step 11.

Add 50µl Tris Buffer.

#### Enzymatic procedure

##### Step 12.

80µ of this mixture are added to 156.8µ Glucose oxidase (aliquot in the freezer) and 3.2µl o-Dianisidine (freshly added from the fridge).

#### Enzymatic procedure

##### Step 13.

Incubate for exactly 30 min at 37°C.

 DURATION

00:30:00

#### Enzymatic procedure

##### Step 14.

Stop reaction by adding 160µl Sulfuric Acid.

#### Enzymatic procedure

##### Step 15.

Measure at 540nm at the (nanoDrop) spectrometer.