



Sep 22, 2019

LB Agar (1 Liter) V.2

Alba Balletbó¹¹Wageningen University

1

Works for me

dx.doi.org/10.17504/protocols.io.7jvhkn6



Alba Balletbó



GUIDELINES

Mixture must be autoclaved before use.

MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

Distilled Water

Agar

NaCl

53014

Sigma Aldrich

Tryptone

BP1421-500

Fisher Scientific

Yeast Extract

211930

Thermo Fisher

- 1 Add LB agar components in a 1 L autoclavable. LB agar components (1 L):

NaCl	10 gram
Tryptone	10 gram
Yeast Extract	5 gram
Agar	15 gram

- 2 Add 1 L of distilled water into the bottle.
- 3 Vigorously shake the bottle.
- 4 Place the bottle into the autoclave and sterilize.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited