



Apr 16, 2019

Working

Western blot analysis [↗](#)

Version 2

PLOS One

Kiichi Hirota¹¹Kansai Medical University[dx.doi.org/10.17504/protocols.io.x9mfr46](https://doi.org/10.17504/protocols.io.x9mfr46)

Yoshiyuki Matsuo

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0215072>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Sumi C, Matsuo Y, Kusunoki M, Shoji T, Uba T, Iwai T, Bono H, Hirota K (2019) Cancerous phenotypes associated with hypoxia-inducible factors are not influenced by the volatile anesthetic isoflurane in renal cell carcinoma. PLoS ONE 14(4): e0215072. doi: [10.1371/journal.pone.0215072](https://doi.org/10.1371/journal.pone.0215072)

MATERIALS

| NAME | CATALOG # | VENDOR |
|--|------------|-----------------------------|
| Purified Mouse Anti-Human HIF-1α Clone 54/HIF-1α | 610959 | BD Biosciences |
| HIF-2 alpha/EPAS1 Antibody | NB100-122 | Novus Biologicals |
| HIF-1β/ARNT (D28F3) XP® Rabbit mAb | 5537 | Cell Signaling Technology |
| Anti α-Tubulin Monoclonal Antibody | 017-25031 | Fujifilm Wako Pure Chemical |
| Anti-Mouse IgG HRP-Linked Whole Ab Sheep | NA931 | Ge Healthcare |
| Anti-Rabbit IgG HRP-Linked Whole Ab Donkey | NA934 | Ge Healthcare |
| RIPA buffer | 16488-34 | Nacalai Tesque |
| cOmplete™ Protease Inhibitor Cocktail | 4693116001 | Roche |
| Blocking One | 03953-95 | Nacalai Tesque |
| Chemi-Lumi One Super | 02230-14 | Nacalai Tesque |
| DC™ Protein Assay Kit | 500-0112 | BIO-RAD |

MATERIALS TEXT

Antibody dilutions

Primary antibodies

- Anti-HIF-1α 1:1000
- Anti-HIF-2α 1:1000
- Anti-HIF-1β 1:1000
- Anti-α-tubulin 1:2000

Secondary antibodies

- Anti-mouse-IgG 1:10,000
- Anti-rabbit-IgG 1:10,000

Cell lysis

- 1 Wash the cells with ice-cold PBS and add ice-cold RIPA buffer with cOmplete™ Protease Inhibitor Cocktail.
- 2 Scrape cells and transfer the suspension to a 1.5 ml microcentrifuge tube.
- 3 Incubate the lysate on ice for 15 minutes.
- 4 Centrifuge at 10,000 x *g* for 5 minutes at 4°C and collect the supernatant in a new microcentrifuge tube.
- 5 Determine protein concentration by using DC™ Protein Assay Kit.

Separation of proteins by gel electrophoresis

- 6 Load 35 µg of total protein on a SDS-PAGE gel.
- 7 Run the gel for 1 h at 100 V.

Semi-dry membrane transfer

- 8 Transfer proteins onto a PVDF membrane using Trans-Blot Turbo™ Transfer System.

Western blotting

- 9 Place the blot in Blocking One and incubate with agitation for 20 minutes.
- 10 Incubate the blot with appropriate dilutions of primary antibody at 4°C overnight.
- 11 Place the blot in TBST and wash three times, 5 minutes each.
- 12 Incubate the blot with the species-matched horseradish peroxidase-conjugated secondary antibody for 1 hour.
- 13 Place the blot in TBST and wash three times, 5 minutes each.

Chemiluminescent Detection

- 14 Add Chemi-Lumi One Super to cover the blot and incubate for 1 minute.

15 Detect chemiluminescence using LAS-4000mini image analyzer.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited