

# Modified Filter-Transfer-Freeze ("FTF") Technique for Raman Microspectroscopic Analysis of Single Cells Version 2

Gordon T. Taylor, Zhuo-Qun "Emma" Li, Elizabeth Suter, Stephanie Chow

## Abstract

A method for concentrating aquatic microbes for single-cell analysis by Raman microspectroscopy. Modified from Hewes and Holm-Hansen (1983): A Method for Recovering Nanoplankton from Filters for Identification with the Microscope: The Filter-Transfer-Freeze (FTF) Technique. *Limnology and Oceanography*. 28(2): 389-394.

Many materials commonly used in microscopy to examine individual microorganisms contribute contaminating Raman scattered emissions, e.g., most types of membranes and filters, glass slides, cover slips, immersion oil, glycerol, etc. and are therefore unsuitable for this application. We developed a protocol that enables single-cell analysis while minimizing extraneous Raman scattered photons. The technique transfers cells from popular polycarbonate membranes to mirror-finished 304 stainless steel slides (1x3x0.0235") supplied by Stainless Supply® (Monroe, NC USA).

**Citation:** Gordon T. Taylor, Zhuo-Qun "Emma" Li, Elizabeth Suter, Stephanie Chow Modified Filter-Transfer-Freeze ("FTF") Technique for Raman Microspectroscopic Analysis of Single Cells. **protocols.io**

dx.doi.org/10.17504/protocols.io.ikqccvw

**Published:** 26 Jun 2017

## Before start

Additionally, a frozen aluminum block (previously stored at -80°C) is necessary as well as typical lab supplies such as pipettors and forceps.

**Note:** Polycarbonate membrane filters of various pore sizes can be used.

## Materials

● Methanol [PA-33900HPLCCS4L](#) by [P212121](#)

✓ MilliQ water by Contributed by users

Acetone solution [48358 SUPELCO](#) by [Sigma Aldrich](#)

- ✓ Isopropanol by Contributed by users
- 304 Stainless Steel Sheet SS-304-8-24G by [Stainless Supply](#)
- Polycarbonate Membrane 0.2µm pore size GTTP02500 by [Emd Millipore](#)

## Protocol

### Prepare Stainless Steel Slides

#### Step 1.

To begin, completely peel off protective film applied to the polished side of the steel slide by the manufacturer. This is the side where the sample will applied.

### Prepare Stainless Steel Slides

#### Step 2.

Acetone clean in ultrasonic bath for 15 minutes.

### Prepare Stainless Steel Slides

#### Step 3.

Rinse with DI water thoroughly.

### Prepare Stainless Steel Slides

#### Step 4.

Isopropanol clean in ultrasonic bath for 15 minutes.

### Prepare Stainless Steel Slides

#### Step 5.

Rinse with DI water thoroughly.

### Prepare Stainless Steel Slides

#### Step 6.

Methanol clean in ultrasonic bath for 15 minutes.

### Prepare Stainless Steel Slides

#### Step 7.

Rinse with DI water thoroughly.

### Prepare Stainless Steel Slides

#### Step 8.

MilliQ water clean in ultrasonic bath for 15 minutes.

### Prepare Stainless Steel Slides

## Step 9.

Rinse by MilliQ water thoroughly.

## Prepare Stainless Steel Slides

### Step 10.

Allow slide to air dry in lint-free environment.

## Filter

### Step 11.

Concentrate cells from fixed suspension onto white polycarbonate membranes by filtering appropriate volume.

#### 📌 NOTES

**Liz Suter** 23 Jun 2017

Some fixatives may be problematic for Raman analysis. Formaldehyde (2% final conc.) performs well, but 1-2% glutaraldehyde produces unacceptable fluorescence during Raman interrogation.

**Liz Suter** 23 Jun 2017

Other fixatives have not been rigorously evaluated.

## Rinse

### Step 12.

Rinse thoroughly with PBS (if seawater sample) or MilliQ H<sub>2</sub>O.

## Cut Wedges

### Step 13.

Cut filter into wedges of desired size with EtOH-sterilized scissors, indicating with pencil marks the top side of the filter.

## Hybridize

### Step 14.

If desired, subject filter wedge to standard fluorescent in situ hybridization (oligo-FISH) protocols with oligonucleotide probes of your choosing to selectively observe specific taxa.

#### 📌 NOTES

**Liz Suter** 23 Jun 2017

Samples prepared by catalyzed reporter deposition FISH (CARD-FISH) have **not** been evaluated for Raman microspectroscopic analysis and may present special challenges.

## Prepare for Transfer

### Step 15.

Place a droplet of sterile MilliQ H<sub>2</sub>O (2-5 µl depending on wedge size) onto polished side of stainless steel slide, near the edge.

## Transfer

### Step 16.

Place the filter wedge upside down on the droplet (see photo), making sure an edge of the wedge is hanging off the slide.

## Freeze

### Step 17.

Place the entire stainless steel slide on a frozen aluminum block (-80°C).

## Remove filter

### Step 18.

Once the droplet is frozen (after a few seconds), quickly peel the filter wedge away from surface with forceps. Cells are now transferred to stainless steel. Allow sample to dry at room temperature in the dark and store frozen until ready for Raman microspectroscopic surveys.

## 📌 NOTES

**Liz Suter** 23 Jun 2017

Transfer efficiency of cells from filter wedge will vary depending on a number of variables, e.g., cell types, debris on the filter, mode of preservation, etc. With *Synechococcus* sp. cells, for example, we have routinely achieved transfer efficiencies of 51 - 77% to mirror-finished stainless steel slides. Your results may vary.