

Cyanobacteria Trace Metal Mixture (CTMM)

Matthew Sullivan Lab

Abstract

Preparation of the trace metal mixture for addition to seawater for the cultivation of marine cyanobacteria, *Prochlorococcus* and *Synechococcus*

Citation: Matthew Sullivan Lab Cyanobacteria Trace Metal Mixture (CTMM). [protocols.io](https://doi.org/10.17504/protocols.io.c8nzvd)

[dx.doi.org/10.17504/protocols.io.c8nzvd](https://doi.org/10.17504/protocols.io.c8nzvd)

Published: 30 Nov 2015

Guidelines

All chemicals used for cultivating cyanobacteria should be of the highest quality.

All re-useable glassware and plastic should be acid-washed and autoclaved.

Use only plastic or teflon spatulas or dust-free weigh paper for weighing out chemicals.

Before start

Prepare the six trace metal stock solutions as described in the Primary Trace Metal Stocks protocol:

- 0.08M $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ ((2.30 g/100 mL)
- 0.05M $\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$ (1.19 g/100 mL)
- 0.90M $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$ (17.81 g/100 mL)
- 0.03M $\text{Na}_2\text{MoO}_4 \cdot 2\text{H}_2\text{O}$ (0.726 g/100m mL)
- 0.10M Na_2SeO_3 (1.73 g/100 mL)
- 0.10M $\text{NiCl}_2 \cdot 6\text{H}_2\text{O}$ (2.38 g/100 mL)

Materials

Ethylenediaminetetraacetic acid disodium salt dihydrate [E4884](#) by [Sigma Aldrich](#)

Protocol

Step 1.

Weigh out 0.435 g $\text{Na}_2\text{EDTA} \cdot 2\text{H}_2\text{O}$ using dust-free weigh paper



REAGENTS

Ethylenediaminetetraacetic acid disodium salt dihydrate [E4884](#) by [Sigma Aldrich](#)

Step 2.

Transfer into 100 mL volumetric flask filled with 60 mL Milli-Q water

Step 3.

Dissolve Na₂EDTA by inserting the stopper and inverting flask several times

📌 NOTES

Bonnie Poulos 11 Aug 2015

May have to heat 5 min at 80°C to dissolve.

Step 4.

Weigh out 0.32 g FeCl₃·6H₂O using dust-free weigh paper

Step 5.

Dissolve FeCl₃ into same volumetric flask and mix by inverting several times

Step 6.

Individually add and dissolve 100 µl of each of the six primary trace metal stocks described in the Primary Trace Metal Stocks protocol. The six primary trace metals to add are ZnSO₄, CoCl₂, MnCl₂, Na₂MoO₄, Na₂SeO₃ and NiCl₂.

🔗 LINK:

<https://www.protocols.io/view/Primary-trace-metal-stocks-c8hzt5>

Step 7.

Adjust volume to 100mL mark with Milli-Q water

Step 8.

Using a polycarbonate syringe, filter through 0.2µm syringe filter into a sterile, acid-washed teflon or polycarbonate container in a laminar flow hood

Step 9.

Store sterile CTMM (cyanobacteria trace metal mixture) at 4°C