

# How to make a 1 M PMSF Stock Solution

Sean Seaver

## Abstract

How to make a 1 M PMSF Stock Solution

**Citation:** Sean Seaver How to make a 1 M PMSF Stock Solution. **protocols.io**

dx.doi.org/10.17504/protocols.io.ci8uhv

**Published:** 03 Nov 2014

## Guidelines

PMSF is hygroscopic and is highly unstable in aqueous solutions. A typical working solution of PMSF is 1mM.

[PMSF](#) is an anti-protease! Do not allow to come into contact with skin.

PMSF is inactivated in aqueous solutions. The rate of inactivation increases with increasing pH and is faster at 25°C than at 4°C. The half-life of a 20mM aqueous solution of PMSF is ~35 minutes at pH 8.0. This short half-life means that aqueous solutions of PMSF can be safely discarded after they have been rendered alkaline (pH >8.6) and stored for several hours at room temperature.

## Materials

 PMSF [RP-P20270](#) by [P212121](#)

## Protocol

### Step 1.

Weigh 1.74 g PMSF

 **AMOUNT**

2 g Additional info:

 **REAGENTS**

 PMSF [RP-P20270](#) by [P212121](#)

 **DURATION**

00:03:00

### Step 2.

Add DMSO to a final volume of 100 ml.

Dissolve completely.

 **REAGENTS**

 DMSO [GB-D-360](#) by [P212121](#)

 **DURATION**

00:03:00

### Step 3.

Stocks may be kept at -20°C for up to 6 months.

PMSF is used at a final concentration of 0.1- 1.0 mM



**REAGENTS**

 PMSF [RP-P20270](#) by [P212121](#)