



May 13,
2019

Working

U Michigan - Optokinetic Measurements of Visual Acuity and Contrast Sensitivity [↗](#)

David A. Antonetti¹

¹University of Michigan - Ann Arbor

[dx.doi.org/10.17504/protocols.io.x99fr96](https://doi.org/10.17504/protocols.io.x99fr96)

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

A virtual optometry system is used to quantify the spatial vision of laboratory animal.

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=306&docType=Protocol>

- 1 Animal is placed inside a square box displaying a rotating cylinder comprised of a vertical sine wave grating is calculated and drawn in virtual three-dimensional coordinate space on four computer monitors facing the animal to form a square.
- 2 Animal stands unrestrained on a platform in the center of the square
- 3 The animal's head movement is tracked for reflexive head and neck movements in response to the grating rotating around the animal
- 4 The spatial frequency of the grating is clamped at the viewing position by repeatedly re-centering the cylinder on the head in real time
- 5 Visual acuity is quantified by increasing the spatial frequency of the grating until an optomotor response could not be elicited
- 6 Contrast sensitivity is measured by identifying the minimum contrast that generates tracking over a range of spatial frequencies



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited