

SQK-MAP005 protocol for library preparation for Nanopore sequencing

Benjamin Istace, Anne Friedrich, Léo dAgata, Sébastien Faye, Emilie Payen, Odette Beluche, Claudia Caradec, Sabrina Davidas, Corinne Cruaud, Gianni Liti, Arnaud Lemainque, Stefan Engelen, Patrick Wincker, Joseph Schacherer, Jean-Marc Aury

Abstract

This protocol describes the library preparation for Nanopore sequencing according to the SQK-MAP005 protocol.

It accompanies the *GigaScience* publication:

Benjamin Istace, et al. (2017) De novo assembly and population genomic survey of natural yeast isolates with the Oxford Nanopore MinION sequencer. *GigaScience*...

Citation: Benjamin Istace, Anne Friedrich, Léo dAgata, Sébastien Faye, Emilie Payen, Odette Beluche, Claudia Caradec, Sabrina Davidas, Corinne Cruaud, Gianni Liti, Arnaud Lemainque, Stefan Engelen, Patrick Wincker, Joseph Schacherer, Jean-Marc Aury SQK-MAP005 protocol for library preparation for Nanopore sequencing. **protocols.io**

dx.doi.org/10.17504/protocols.io.gvubw6w

Published: 06 Jan 2017

Protocol

DNA fragmentation

Step 1.

Fragment DNA (2µg for a 8Kb or 6 to 10µg for a 20Kb size) by using Covaris gTube

DNA repair

Step 2.

Perform PreCR treatment (NEB PreCR® Repair Mix) of fragmented DNA

Clean Up

Step 3.

The DNA repair reaction was cleaned up with Agencourt AMPure XP beads (1x)

End Repair

Step 4.

DNA fragments were End-repaired by using NEBNext® End Repair Module

Clean Up

Step 5.

The end-repair reaction was cleaned up with Agencourt AMPure XP beads (1x)

dA-tailing

Step 6.

DNA fragments were dA-tailed by using the NEBNext® dA-Tailing Module

Clean Up

Step 7.

The dA-tail reaction was cleaned up with AMPure beads (1x)

Adaptors ligation

Step 8.

Nanopore adaptors were ligated to the dA-tailed DNA fragment by using NEB Blunt/TA Ligase Master Mix

Clean Up

Step 9.

The ligation reaction was cleaned up with His-Tag Dynabeads (1x)