



Jun 12,  
2019

## ASW+NO3 medium

Roscoff Culture Collection<sup>1</sup>

<sup>1</sup>CNRS-Sorbonne Université, Station Biologique, Place G. Tessier 29680 Roscoff FRANCE

Working

dx.doi.org/10.17504/protocols.io.s5zeg76

Roscoff Culture Collection



Daniel Vaultot

Station Biologique, Roscoff, France



### ABSTRACT

Medium used for cyanobacteria based on artificial seawater (ASW)

- 1
  - Dissolve 25 g of NaCl in MilliQ water
  - To this solution, add :

| Quantity              | Compound      | Stock Solution   | Concentration in medium (in mM) |       |
|-----------------------|---------------|--|---------------------------------|-------|
| <i>mL stock /LASW</i> | <i>g/LASW</i> |  |                                 |       |
| 10                    | 0,75          | Sodium nitrate (NaNO <sub>3</sub> )                                  | 75 g/L                          | 8.8   |
| 10                    | 2             | Magnesium chloride hexahydrate (MgCl <sub>2</sub> 6H <sub>2</sub> O) | 200 g/L                         | 9.8   |
| 5                     | 0,5           | Potassium chloride (KCl)   | 100 g/L                         | 6,7   |
| 5                     | 0,5           | Calcium chloride (CaCl <sub>2</sub> )                                | 100 g/L                         | 4.5   |
| 10                    | 3,5           | Magnesium sulfate heptahydrate (MgSO <sub>4</sub> 7H <sub>2</sub> O) | 350 g/L                         | 14.2  |
| 5,5                   | 1,1           | TRIS-Base  | 200 g/L                         | 9.08  |
| 2,5                   | 0,03          | Dipotassium phosphate (K <sub>2</sub> HPO <sub>4</sub> )             | 12 g/L                          | 0.172 |

- Adjust the pH to 8 with concentrated HCl
- Adjust to 999 mL with milliQ water
- Add 1 mL of trace metals (see receipe below)
- Autoclave the medium

## Trace metal stock solution

- 2
- Dissolve all these components separately in milliQ water :

| Quantity | Compound  |
|----------|---|
| 2.86g    | Boric acid ( $\text{H}_3\text{BO}_3$ )  |
| 1.81g    | Manganese (II) chloride tetrahydrate ( $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$ )  |
| 0.222g   | Zinc sulfate monohydrate ( $\text{ZnSO}_4 \cdot \text{H}_2\text{O}$ )               |
| 0.390g   | Sodium molybdate dihydrate ( $\text{Na}_2\text{MoO}_4 \cdot 2\text{H}_2\text{O}$ )  |
| 0.008g   | Copper sulfate pentahydrate ( $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ )           |
| 0.0494g  | Cobalt nitrate hexahydrate ( $\text{Co}(\text{NO}_3)_2 \cdot 6\text{H}_2\text{O}$ ) |
| 3.0g     | Ferric chloride hexahydrate ( $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$ )           |
| 0.5g     | EDTA magnesium disodium ( $\text{EDTA}(\text{Na}_2\text{Mg})$ )                     |

- Combine the various solutions after full dissolution
- Make final volume up to 1L with milliQ
- Store in refrigerator



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited