

Preparing AQUIL Media

Alexis Wasson

Abstract

Procedure in order to prepare EDTA media (also known as Aquil).

Citation: Alexis Wasson Preparing AQUIL Media. **protocols.io**

dx.doi.org/10.17504/protocols.io.fp9bmr6

Published: 16 Sep 2016

Guidelines

Reagents

- Sodium Silicate
- Sodium Nitrate:
- Sodium dihydrogen Phosphate:
- Vitamin B12
- Biotin
- Thiamine
- $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$
- $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$
- $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$
- $\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$
- $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$
- $\text{Na}_2\text{MoO}_4 \cdot 2\text{H}_2\text{O}$
- Na_2SeO_3

Protocol

Step 1.

Nutrient Stock

Add 1.42g Sodium Silicate to 50ml of UPW in a 50ml vial.

Add 0.43g Sodium Nitrate to 50ml of UPW in a 50ml vial.

Add 0.069g Sodium dihydrogen Phosphate to 50ml of UPW in a 50ml vial.

Step 2.

Vitamin Stock

Add 0.037g of Vitamin B12 to 50ml of UPW in a 50ml vial.

Add 0.0061g of Biotin to 50ml of UPW in a 50ml vial.

Add 1.5g Thiamine to 50ml of UPW in a 50ml vial.

Step 3.

EDTA Stock

Add all of the following dry chemicals into one clean 50ml vial in order to make EDTA stock.

- 46g of EDTA
- 0.0135g of $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$
- 0.0011g of $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$
- 0.002g of $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$
- 0.001g of $\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$
- 0.002g of $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$
- 0.004g of $\text{Na}_2\text{MoO}_4 \cdot 2\text{H}_2\text{O}$

- 0001g of Na_2SeO_3

Step 4.

Nutrient and vitamin stocks should be chelexed and filter-sterilized before use. EDTA(metals) stock should be filter-sterilized. To make AQUIL, add the following reagents to 1L of chelexed, microwave-sterilized ASW, or 1L microwave sterilized surface seawater. Add 1ml EDTA stock, 1ml sodium silicate, 1ml sodium nitrate, 1ml sodium dihydrogen phosphate, 100ul vitamin B₁₂, 100ul biotin and 100ul thiamine stock. Invert bottle twice for mixing.