

Conjugation of *Thalassiosira pseudonana* Version 3

Jernej Turnsek

Abstract

This protocol was used to express a nourseothricin (*nou/nat*) resistance gene and a gene encoding a silaffin precursor TpSil3p-APEX2 fusion protein in *Thalassiosira pseudonana* (*Tp*) strain [CCMP1335](#) using conjugation by largely following Karas et al. (2015). In addition to the protocol, I am listing some of my observations working with this diatom species. Please check the 'Before start' and 'Guidelines' sections for more information including a list of *Tp* episomes I've constructed and which are available upon request. For additional background on my overarching experimental aim, please refer to [this document](#).

1. E. V. Armbrust, The Genome of the Diatom *Thalassiosira Pseudonana*: Ecology, Evolution, and Metabolism. *Science*. **306**, 79-86 (2004).
2. B. J. Karas et al., Designer diatom episomes delivered by bacterial conjugation. *Nat. Commun.* **6**, 6925 (2015).

Keeping in mind that diatom episomes contain a yeast centromere, it would be very interesting to see if conjugation could be adapted to organisms beyond diatoms.

Citation: Jernej Turnsek Conjugation of *Thalassiosira pseudonana*. [protocols.io](#)

[dx.doi.org/10.17504/protocols.io.f55bq86](https://doi.org/10.17504/protocols.io.f55bq86)

Published: 19 Oct 2016

Guidelines

[1] 'Crashing phenomenon' in *Tp* cultures

In my experience *Tp* CCMP1335 has to be subcultured frequently - every ~2 weeks - to avoid 'crashing'. I usually dilute my cultures 1:10 in fresh L1 medium. By 'crashing' I mean sudden take over of what I believe is a naturally associated bacterium in the culture. This usually happens in the stationary phase over the course of 1 or 2 days during which brown diatom culture becomes milky due to bacteria. It could be a standard contamination issue, but I've pretty consistently observed this phenomenon even after ordering completely fresh stocks from the NCMA. I analyzed the bacterium via 16S sequencing and some simple growth experiments and confirmed it was a *Pseudomonas* species. You can read more about how I've done this [here](#) and find the paper associated with my 16S and *rpoD* (*Pseudomonas*-specific) analysis [here](#).

Interestingly, such *Tp-Pseudomonas* interplay was [already observed in the late 1970's](#).

[2] Episomes for *Tp* conjugation available upon request

All episomes used in my conjugation experiments including their maps are available upon request. Send an email to turnsek@fas.harvard.edu, turnsek@ucsd.edu or jturnsek@jcv.org and we'll go from there/I can provide more information. Here is a list of available episomes:

pTpPuc3 - empty nourseothricin resistance cassette-containing episome for *T. pseudonana* conjugation

eGFP - 2 different enhanced GFP constructs

[DreidelTeal](#) - blue chromogenic protein

TpSil3p-eGFP - frustule-associated silaffin precursor protein TpSil3p tagged with eGFP

TpSil3p-APEX2-FLAG - frustule-associated silaffin precursor protein TpSil3p tagged with APEX2 and FLAG

TpSil3p-APEX2-6xHis - frustule-associated silaffin precursor protein TpSil3p tagged with APEX2 and 6xHis

NOTE:

All the genes above were cloned into pTpPuc3 Kan^R backbone via Gibson Assembly starting 36 nucleotides downstream of the HIS3 gene. They are flanked by a constitutive fcp promoter and an associated terminator. Here are their sequences in FASTA format:

>Tp_fcp_promoter

```
GCGCTTTTTCCGAGAACTCCCCATAAGTCAACGGCTCCAATCAAGAATGTATCCGACAACGGCGAGCATAGCAA  
CACGTCCGTCTTTGGAGTAGAATCATCATGTTGTGGATGAATACACAGA
```

```
TGAATGACATTAAGCATGAACATGTTAGAGAGTAGGAGGTAGAGATTGATATGGTAGCATTGCGATGTTTGT  
TTTTGGTCAGCATATGATGAGTGGATACCAATATGATGAAAGTTGAATCT
```

```
CGCGTTTGAGCTCAGCGGTACGTTATTGATCGAAAGTAGCCTGATCAAAATCCTTGGAGAGTACAAGAGGATCA
```

AGAATCCAGTGGGGGCGATAACTCCAAGCTCGTTCTCAAAGAGGCAA
 TGGAGGTAGAACTCATCCCAGTTGAGAAGAAGTGAAGGCAGTGGCGGTGGCGAAAGCAGAGGCAACGAGGA
 CAGACTTCCTGTGGGTTGATGCAACGAATATTTCCAGAAGGAGAAG
 TTTAGAGAGTTGAACCGCTACCTACAATGACAAAGTATCGTATCGATTTTGATGTTGGTTGGTTATGAATTCAAA
 CTGTAAGTTGGATTGTGAGAAGATCAGAAGTTGAACGAACACATCTTTC
 CGATCATTACCTCCACACTGCAACAACACGGTACTTCTTCCGCGGCAGGTCTCTGTCGCCATTCTCTTGTCCTG
 TTGTTGGCTGTGAGACGAGGAAAGCAACGACAAGTTTCACAAAAGGG
 AGTTCCTTTAACGAGATATGTTTTTTATAAAGAGTCCCAATAGAAAGACAAATTGATTCTCCGTGCAAACGCGC
 AAATAAACACCACGTCCATTATATCCATATCTTTCAGAGTATCCAACAAGT
 GTTGAAGGACAGGTAGTTGAAGTAACGTATCTTCCCCCTCGACTGGATCCATCAACAAGGCGAACAAATCCATT
 CAACCTCTCATAAATTATCTGATTTACCAAACC

>Tp_fcp_terminator

ATACTGGATTGGTGAATCAATGAGCCGTAGCACAATGGTTACATTCGGCTAGCTAAGATCCAATGGCAAGGACC
 AAGTGCTGGAACCTGTTTTGCTTTAGCAGATCTTAGCGTGAGAGGTATT
 TGTCTCTGTCAGGAGTAGATAGTAGATGTTCTTTTTAACTAAAATGCTAACTGTTCCGAATTCCTCATCGCAGC
 TAATCCGTACATCAAAGACAAAATGCTAGGTATGTGTACTACATCTCC
 TGTTGCTAGATAAGACATATGATAGGAAACACACCATCAATAGTCATTGTAGCTTTACTTATACTACGCATTTGCA
 CTTTCCCCTGAGTGGCAGAGGCGCATTGAGAAAATCGATCTCAACATAG
 TTTATGTAGCATCCCCTAGATCCATTACTTTAAGTCTCCTTCGTCTTTGGTGTAGGCATGTTGGACACAACGAGG
 TAAACACAACACAACAATGTGTCCAGCAAAGTAGTAGCTGCTCCAGT
 TCT

Before start

[1] Preparation of donor bacteria

1.1. Any *E. coli* strain will probably work just fine, but I stuck with the one from the original publication: [EPI300](#) (Epicentre, catalog # EC300110).

1.2. Electroporate pTA-Mob into EPI300 cells and select on LB agar plates w/ 10 or 20 µg/mL gentamycin.

NOTE:

Please refer to the following publication to learn more about the mobilization plasmid pTA-Mob: T. A. Strand, R. Lale, K. F. Degnes, M. Lando, S. Valla, A new and improved host-independent plasmid system for RK2-based conjugal transfer. *PLoS One*. **9**, e90372 (2014).

1.3. Make pTA-Mob-containing strain chemically competent.

1.4. Transform your engineered pTpPuc3 cargo plasmid into pTA-Mob-containing cells.

NOTE:

If you already have a working *Tp* plasmid that you successfully transformed via biolistics or some other method, you can easily convert it into a conjugation-compatible vector. Refer to Karas et al. (2015) for details on how to perform this conversion.

1.5. Select double transformants on LB agar plates with 10 or 20 µg/mL gentamycin and 50 or 100 µg/mL kanamycin.

1.6. Store plates at 4°C.

[2] Preparation of L1 medium and L1 1% agar plates

I was using [L1 Medium Kit](#) and [Gulf of Maine Seawater](#), both from NCMA, to make my liquid medium and plates.

NOTE:

I also had success with [ESAW](#), [f/2](#) and [Daigo's IMK](#) media, but never tested them for the full conjugation protocol. Daigo's IMK can be ordered on [this website](#) by using a search term "IMK" (select catalog # 392-01331).

2.1 L1 medium

I would usually make batches of 1L by combining seawater, N, P and Si sources from the kit followed by autoclaving, cooling down to ~50°C and adding vitamins and antibiotics. The resulting medium

had pH ~8.5. Before starting a diatom culture I would adjust the pH with HCl to ~8.0. Starting pH anywhere between 7.8 and 8.2 should be fine. All my liquid cultures were grown in a batch mode at constant light, ~18°C and moderate shaking (~100 rpm). I never buffered the medium and/or aerate with air or CO₂.

NOTE:

The resulting medium contain some precipitates due to Si and seawater-derived compounds. You can use a combination of microwaving and filter sterilization instead which should prevent precipitation.

2.2 1/2xL1 1% agar plates w/ 5% LB

I would always prepare 0.5L and make ~12 plates 40 mL each. For 0.5L combine and mix seawater, N, P and Si sources from the kit, add agar and 5% LB. Autoclave, cool down to ~50 °C and add vitamins.

NOTES:

- 5% LB is **v/v not w/v**. I had success preparing plates w/v, but **v/v is what you should use**. In other words: 1.25 g LB powder per 1 L.
- 1/2xL1: cut all N, P, Si and vitamin volumes from the kit in half.
- I had success with Bacto Agar and [Sigma's high purity agar](#).

2.3.1/2xL1 1% agar plates w/ 100 µg/mL nourseothricin

Same as above except the LB part; add nourseothricin when cooled down to ~50 °C. My nourseothricin stocks are 200 mg/mL.

NOTE:

I observed a lot of background - false positive - colonies on plates with only 50 µg/mL nourseothricin so I switched to 100 µg/mL. You can read more about nourseothricin on the [following website](#) which I found very useful.

Protocol

Growth and preparation of *E. coli* donor

Step 1.

Pick bacterial colonies from your Gent+Kan plates and inoculate 10 mL LB medium. Grow overnight.

📌 NOTES

Jernej Turnsek 19 Oct 2016

Start growing colonies early in the afternoon the day before conjugation to get enough biomass. ~16-20 hrs should suffice. Expect OD₆₀₀ ~2 next morning.

Growth and preparation of *E. coli* donor

Step 2.

Measure OD₆₀₀ and start a 150 mL LB subculture (recommended starting OD₆₀₀ either 0.05 or 0.1).

Growth and preparation of *E. coli* donor

Step 3.

Grow at 37°C until OD₆₀₀ reaches 0.3-0.4.

Growth and preparation of *E. coli* donor

Step 4.

Centrifuge at 4,000 rpm, 10°C, for 10 min.

🕒 DURATION

00:10:00

Growth and preparation of *E. coli* donor

Step 5.

Decant supernatant and resuspend in 800 µL SOC.

Growth and preparation of *Thalassiosira pseudonana* cells

Step 6.

Tp was cultured in L1 medium as described in the 'Before start' section.

📌 NOTES

Jernej Turnsek 19 Oct 2016

- I've never tested if cell density before spinning cells down matters with respect to final conjugation outcome same as the Alverson Lab did for *P. tricornutum*. They observe ~8 x 10⁶ cells/mL to be a sweet spot. In my hands spinning cells down at ~4-8 x 10⁶ cells/mL seemed to work fine. According to my information harvesting *T. pseudonana* at ~0.8 x 10⁵ cell/mL works best for biolistic experiments.

- I counted cells with [BioRad's TC20 automatic cell counter](#).

Growth and preparation of *T. pseudonana* cells

Step 7.

Spin down 2×10^8 cells at 4000 rpm, 10°C, for 10 min.

DURATION

00:10:00

Growth and preparation of *T. pseudonana* cells

Step 8.

Decant supernatant and resuspend pellet in 1 mL L1 medium.

NOTES

Jernej Turnsek 19 Oct 2016

I usually saved some supernatant and used it for resuspension and later on for scraping cells before selection (see Steps 14 and 15).

Conjugation

Step 9.

Mix 200 μ L *T. pseudonana* cells and 200 μ L *E. coli* cells in a 1.5 mL tube.

Conjugation

Step 10.

Pipette up and down a few times.

Conjugation

Step 11.

Plate on 1/2xL1 1% agar plates w/ 5% LB.

NOTES

Jernej Turnsek 19 Oct 2016

Make sure the plates are dry. Leave them open at 37 °C for an hour or so before plating diatom-bacteria co-cultures.

Conjugation

Step 12.

Incubate in dark at 30°C for 90 minutes.

DURATION

01:30:00

Conjugation

Step 13.

Move plates to your standard *Tp* growth conditions - in my case 18°C and constant light - for 4 hours.

 DURATION

04:00:00

Conjugation

Step 14.

Add 1 mL L1 medium and scrape with a [cell scraper](#) or [L spreader](#).

Selection

Step 15.

Expect to recover 500 µL co-culture suspension after scraping.

Selection

Step 16.

Plate 200 µL of the resulting suspension on pre-dried 1/2xL1 1% agar plates w/ 100 µg/mL nourseothricin.

Selection

Step 17.

Leave at 18°C and constant light until colonies appear - 2 weeks.

Screening

Step 18.

Here are a few ways to confirm the presence and expression of your heterologous gene in resulting colonies:

1. Growth under selection pressure.

Make sure colonies are able to grow in liquid L1 medium with 100 µg/mL nourseothricin (Nou100). Pick colonies with a small tip or better a toothpick and inoculate 0.5 mL L1 medium. Once you observe growth subculture in larger volume.

2. PCR

- Use 1 µL of diatom culture as a template to amplify your expression cassette.
- Confirm the absence of donor DNA by amplifying *E. coli*-specific genes.

NOTE:

I used primers to amplify the [corC](#) gene. Make sure you always run *E. coli* EPI300 positive control.

- Confirm the absence of live donor cells by plating some diatom culture on LB plates w/o antibiotics.

NOTE:

Any remaining donors cells and donor DNA are gone after a few liquid subcultures. You can also patch colonies on fresh 1/2xL1 Nou100 plates.

3. RT-PCR

- Purify total RNA from Nou100 diatom culture, convert it to cDNA and use it to run a PCR with heterologous gene-specific primers.
- An example of my result with controls can be found [here](#).

4. Episome recovery

- Perform a diatom miniprep as described by Karas et al. (2015).
- Transform *E. coli* with diatom-derived DNA.
- Select on LB agar plates with 50 µg/mL kanamycin.
- Miniprep, digest and analyze on a gel.

NOTE:

You will learn more about the state of your episome from this analysis.

5. Western blot

6. Protein-specific assays

- enzymatic assay, microscopy etc.

I've had success with 1, 2 and 3; tried 4, 5 and 6 without success.