



Mar 20,
2019

Working

Growth Incremental analysis of thins sections from marine molluscs

Niklas Hausmann¹, Harry K. Robson¹, Chris Hunt²

¹University of York, BioArCh, ²Liverpool John Moores University, Natural Sciences and Psychology,
dx.doi.org/10.17504/protocols.io.y7efzje

Niklas Hausmann
University of York

ABSTRACT

Annual growth patterns in marine mollusc shells are valuable indicators of the condition of marine ecologies through time. In archaeological contexts, the mollusc’s time of death (i.e. the last season of growth) is an indicator of human exploitation patterns throughout the year, enabling the reconstruction of when and how often gathering occurred as well as when sites were occupied. Both pieces of information, growth rate and season of death, are vital for understanding exploitation pressure in the past, and building baselines for modern environmental policies that secure sustainable marine resources. Previously, these parameters have been determined by incremental growth-line or isotopic analyses, which are time consuming and often expensive techniques, thus restricting sample size and the overall robustness of palaeoecological interpretations.

Here, we apply Laser Induced Breakdown Spectroscopy (LIBS) to produce elemental maps (Mg/Ca) with the potential to trace and display growth patterns quickly, and at a reduced cost. We further compare the elemental maps with the results of incremental growth-line analysis to provide a structural context for the geochemical data, and demonstrate the utility of an integrated methodological approach.

Our study was undertaken on 12 European oysters (*Ostrea edulis*) from the Late Mesolithic shell midden at Conors Island, Co. Sligo in the Republic of Ireland. Our LIBS analysis enabled us to accurately and quickly determine repeating growth patterns, which were in agreement with the annual growth increments visible through the microscopic analysis. Based on this comparative dataset, including structural and geochemical patterns, the Late Mesolithic site of Conors Island had been occupied throughout the year. Moreover, our analyses highlight the applicability of LIBS to determine prehistoric seasonality practices as well as biological age and growth at an improved rate and reduced cost than was previously achievable.

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

STEPS MATERIALS

| NAME ▾ | CATALOG # ▾ | VENDOR ▾ |
|----------------------------|----------------|----------|
| Ethanol | | |
| Release Agent | 20-8185-016 | Buehler |
| Epoxy Resin | No 20-8140-12B | Buehler |
| Epoxy Hardener | 20-8142-016 | Buehler |
| Polycrystalline suspension | 40-6016 | Buehler |
| Adhesive | n/a | |

Producing thin sections

1

Section oysters following Milner (2001) and using a precision saw with a diamond blade. First section is at the hinge and in direction of growth, the second (optional) section is done to remove unnecessarily long portions of the shell section.

Milner, N., 2001. At the Cutting Edge: Using Thin Sectioning to Determine Season of Death of the European Oyster, *Ostrea edulis*. Journal of archaeological science, 28(8), pp.861–873.



ISOMET 1000 Precision Saw
Saw

Buehler 11-2180
Using a Buehler Diamond Wafering Blade
(Series 15LC Diamond No. 11-4276)



2 Clean cylindric plastic moulds with



Ethanol

leave to dry for 00:05:00 . Then coat the inside with



Release Agent

by Buehler

Catalog #: 20-8185-016

and dry again for 00:30:00

3

Place sectioned shells face down in 1- and 2-inch plastic moulds depending on the size of the cut hinge.
Then add resin (2 parts by volume of Epoxy Resin mixed with 1 part by volume Epoxy Hardener.



Epoxy Resin

by Buehler

Catalog #: No 20-8140-12B

Leave to cure in an environment of constant humidity and temperature overnight (or 12:00:00)



Epoxy Hardener

by Buehler

Catalog #: 20-8142-016

4

Remove resin block from plastic moulds once hardened, and grind hinge-surface using successively finer metallographic grit papers (P600, P1200 and P2500 grades respectively) on a Buehler Motorpol 2000.



Motorpol 2000
Polariser

Buehler 49-7268



Clean the surfaces of the resin blocks with water and polish further using a polishing cloth and



Polycrystalline suspension

by Buehler

Catalog #: 40-6016

polycrystalline
suspension of 3µm.

Leave to dry overnight (or ⌚ 12:00:00)

5

Attach the polished facets of the resin blocks to a glass slide using Loctite 322 Adhesive, exposed to an ultraviolet light source and leave to harden for

⌚ 24:00:00 .



Adhesive

Catalog #: n/a

6

Once the sample had bonded to the slide, cut the the bodies of the resin blocks from the slide using a Buehler ISOMET 1000 Precision Saw and leave a slice of the embedded hinge of approximately 200 to 500 µm in thickness on the glass slide.☒☒



ISOMET 1000 Precision Saw

Saw

Buehler 11-2180

Using a Buehler Diamond Wafering Blade
(Series 15LC Diamond No. 11-4276)

7

Repeat Step 4, iln order to lightly ground the shell and the resin to produce the final thin section of ~50 µm.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited