

Human Parechovirus A conventional nested genotyping RT-PCR ["Harvala assay"; 2008-2015] Version 3

Ian Mackay

Abstract

I and my team used this assay between 2008-2015; we dubbed it the "Harvala assay". It produces an amplicon that includes the 3' end of VP3 and the 5' end of VP1, spanning the junction.

In silico sequence alignments indicated the highly degenerate oligonucleotides could theoretically detect at least HPeV 1-7, 17 and 18.

Citation: Ian Mackay Human Parechovirus A conventional nested genotyping RT-PCR ["Harvala assay"; 2008-2015]. **protocols.io**

dx.doi.org/10.17504/protocols.io.krxcv7n

Published: 12 Nov 2017

Guidelines

- This protocol assumes the user is familiar with working in a laboratory, with PCR, the thermocycler and software used to run it
- This protocol should be re-evaluated if being used with different reagents, if the oligonucleotide sequences are changed or if the cycling conditions are altered

Protocol

Oligonucleotides...

Step 1.

Round	Name	5'-3' oligonucleotide sequence
Round 1. RT-PCR	HPeV_VP3/1_OS	GAYAATGCIATMTAYCAWATYTGTA
Round 1. RT-PCR	HPeV_VP3/1_OAS	ACWGTRAARATRTCHACATTSATDG
Round 2. nPCR	HPeV_VP3/1_IS	TTYTCMACTGGATGMGGAARAC
Round 2. nPCR	HPeV_VP3/1_IAS	DGGYCCATCATCYTGWGCTGA

OS-outer sense; OAS-outer antisense; IS-inner sense; IAS-inner antisense

Reagents...

Step 2.



REAGENTS

SensiFAST Probe no ROX one-step kit BIO-76005 by [Bioline](#)

MyTaq HS DNA Polymerase BIO-21113 by [Bioline](#)

Reaction setup...

Step 3.

Below is the reaction setup for a single RT-PCR reaction.

Ideally, this work is conducted in a laboratory separate to any space used to *perform* PCR, molecular cloning or the analysis of high concentration DNA.

This volume has been used in 0.2ml tubes.

Multiply this according to the number of reactions you will need, remembering to include a positive control and at least two non-template controls (NTCs).

You may also need to allow some extra volume, depending on the method used to pipette mix into tubes for the run. For example, some robot-loaded tubes can require two reaction 'dead volumes'.

Round 1: RT-PCR

Reagent (stock concentration)	Vol (μL) / reaction	Final concentration
Nuclease free water	7.28	
AN345_panHPeV/LV (200pmol/ul [200uM])	0.06	600nM
AN344_panHPeV/LV (200pmol/ul [200uM])	0.06	600nM
SensiFast OneStep Mix(2x)	10	1X
RNase inhibitor	0.4	
RT/Taq (6U/mL)	0.2	1X
Template extract RNA	2	
Final volume	20 μ l	

Round 2: nPCR

Reagent (stock concentration)	Vol (μL) / reaction	Final concentration
Nuclease free water	12.424	
AN345_panHPeV/LV (200pmol/ul [200uM])	0.038	380nM
AN344_panHPeV/LV (200pmol/ul [200uM])	0.038	380nM
MyTaq Reaction Buffer (5X)	4	1X
MgCl ₂ (25mM)	1.4	
MyTaq HS DNA Polymerase (5U/uL)	0.1	1X
Round 1 amplicon	2	
Final volume	20μl	

Cycling conditions...

Step 4.

This assay has been optimized and for use with a conventional block thermal cycler.

The cycling conditions for the one-step RT-PCR and the nested PCR (nPCR) are as follow:

Round 1: RT-PCR			
45°C	20min		
94°C	2min		
95°C	30s		
60°C	30s		40X
72°C	105s		
72°C	7min		
4°C	∞		

Round 2: nPCR			
94°C	1min		
94°C	30s		
50°C	30s		40X
72°C	105s		
72°C	7min		
4°C	∞		

Step 5.

A positive result is determined by a suitably sized band on an agarose gel after electrophoresis.

Run both reactions on a 1.5% agarose gel after Round 2 is complete. Only the second round amplicon needs to be sequenced, but if the viral load is high enough, there may be a useful, larger band from Round 1 which can be sequenced instead.