



May 10,
2019

Working

U Mass - Triglyceride

Jason Kim¹

¹University of Massachusetts

dx.doi.org/10.17504/protocols.io.x5pfq5n

Mouse Metabolic Phenotyping Centers

Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

This experiment involves a spectrophotometric measurement using Roche Cobas Clinical Chemistry Analyzer. Triglyceride levels may be measured in serum, plasma, and tissue samples. Serum and tissue triglyceride levels are affected by alterations in systemic lipid metabolism, lipid digestion/absorption, and lipid clearance. Serum and tissue triglyceride levels are altered in obesity, insulin resistance, type 2 diabetes, alcoholic steatohepatitis, non-alcoholic fatty liver disease, and non-alcoholic steatohepatitis.

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=167&docType=Protocol>

MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾	CAS NUMBER ▾	RRID ▾
Triglycerides	04657594 190	Roche		
Calibrator f.a.s.	10759350 360	Roche		
Precinorm U plus	12149435 160	Roche		
Precipath U plus	12149443 160	Roche		
NaCl Diluent 9 %	04774230 190	Roche		
Chimneys	11930630 001	Roche		
Cleaner	04774248 190	Roche		
Micro Sample cups	11406680 001	Roche		
NERL High Quality Water	9805	Fisher Scientific		

MATERIALS TEXT

Note:

Roche, [RRID:SCR_001326](#)

Fisher Scientific, [RRID:SCR_008452](#)

BEFORE STARTING

Notes:

✓ Try to use freshly prepared serum and plasma samples for this assay.

✓ No dilution or treatment of the sample is required, but plasma samples should be centrifuged to remove any fibrin/fibrinogen

clumps.

✓ Samples should be stored at 2-8°C for 24 hours prior to analysis. For longer periods, store samples at -70°C, and avoid repeated freeze/thaw cycles.

✓ A 50 µl dead volume is required in addition to sample volume for multi-protein analysis (typically 1-5 µl).

- 1 Perform daily quality control assessment of instrumentation before analysis.
- 2 Load each sample into a specialized micro-sample cup for the clinical chemistry analyzer.
- 3 Select Triglyceride test on display and run the analysis.
- 4 Collect and analyze the data.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited