

Post-stain protocol for GreenGlo™ Safe DNA Dye, 20,000X in Water

Denville Scientific

Abstract

GreenGlo $^{\text{TM}}$, 20,000X in Water, is a non-carcinogenic and non-toxic alternative to Ethidium bromide used for the detection of nucleic acids in agarose gels. It is as sensitive as Ethidium bromide. There is no toxic DMSO as GreenGlo $^{\text{TM}}$ is supplied in water.

GreenGlo [™] has fluorescence excitation maxima at 295 nm and 490 nm. The fluorescence emission maxima is similar to EtBr when bound to DNA – at 530 nm.

Citation: Denville Scientific Post-stain protocol for GreenGlo™ Safe DNA Dye, 20,000X in Water. protocols.io

https://www.protocols.io/view/post-stain-protocol-for-greenglo-safe-dna-dye-20-0-gqcbvsw

Published: 12 Dec 2016

Materials

GreenGlo[™] Safe DNA Dye, 20,000X in Water, 500µl <u>CA3600</u> by <u>Denville Scientific Inc.</u>

Protocol

POST-STAIN

Step 1.

Use 10 to 15 µl per 100 ml of staining solution. Same as when using ethidium bromide.

NOTES

Nicole Clouse 10 Dec 2016

If pre-staining follow: Pre-stain protocol for GreenGlo™ Safe DNA Dye.

POST-STAIN

Step 2.

For an average gel thickness of about 7 mm, stain 30 minutes.

O DURATION

00:03:00

POST-STAIN

Step 3.

Destain of 30 minutes in water.

© DURATION 00:30:00

NOTES

Nicole Clouse 10 Dec 2016

Protect gel and staining solution from light with aluminum foil or place in dark.

DETECTION

Step 4.

Detect bands under UV illumination (yellow or green gelatin- or cellophane filters is recommended for clearer bands) or non-UV LED illuminators such as Blue Light LED illumination.