

# Pulse Field Gel Electrophoresis (PFGE) Protocol for separation of *Chlorella* chromosomal DNA and *Chlorella* virus DNA

Irina Agarkova

## Abstract

**Citation:** Irina Agarkova Pulse Field Gel Electrophoresis (PFGE) Protocol for separation of *Chlorella* chromosomal DNA and *Chlorella* virus DNA. **protocols.io**

[dx.doi.org/10.17504/protocols.io.essbeee](https://dx.doi.org/10.17504/protocols.io.essbeee)

**Published:** 27 Apr 2016

## Guidelines

### Buffers and Solutions

#### Suspension Buffer (SB):

25 mM Tris pH 7.5

1 M Sorbitol

25 mM EDTA pH 8.0

#### Digestion Buffer (DB):

250 mM EDTA pH 9.5

1% N-Lauroylsarcosine

#### 1 × TAE Buffer:

40 mM Tris

20 mM acetic acid

1 mM EDTA

## Formulation of MBBM (Modified Bold's Basal Medium)

### Stock Solutions:

1. 25.0 g  $\text{NaNO}_3$  /1L d- $\text{H}_2\text{O}$
2. 2.5 g  $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$  /1L d- $\text{H}_2\text{O}$
3. 7.5 g  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$  /1L d- $\text{H}_2\text{O}$
4. 7.5 g  $\text{K}_2\text{HPO}_4$  /1L d- $\text{H}_2\text{O}$
5. 17.5 g  $\text{KH}_2\text{PO}_4$  /1L d- $\text{H}_2\text{O}$
6. 2.5 g  $\text{NaCl}$  /1L d- $\text{H}_2\text{O}$
7. 50.0 g disodium EDTA; 31.0 g  $\text{KOH}$  /1L d- $\text{H}_2\text{O}$
8. 4.98 g  $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$  / 1L acidified  $\text{H}_2\text{O}$  (Acidified  $\text{H}_2\text{O}$  is 999.0 mL d- $\text{H}_2\text{O}$  + 1.0 mL concentrated  $\text{H}_2\text{SO}_4$ )
9. 11.42 g  $\text{H}_3\text{BO}_3$  /1L d- $\text{H}_2\text{O}$
10. 8.82 g  $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ ; 1.44 g  $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$ ; 0.71 g  $\text{MoO}_3$ ; 1.57 g  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ ; 0.49 g  $\text{CoNO}_3 \cdot 6\text{H}_2\text{O}$  /1L d- $\text{H}_2\text{O}$

### MBBM preparation:

to 950 mL of d- $\text{H}_2\text{O}$  add:

10.0 mL of stock solutions 1; 2; 3; 4; 5 and 6

1.0 mL of stock solutions 7; 8 and 9

1.0 mL of stock solution 10

1.0 g of bacto-peptone

5.0 g of sucrose

## Formulation of FES Medium (for Chlorella Pbi growth)

### Stock Solutions:

- 1). 10.0 gm  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$  /1L d- $\text{H}_2\text{O}$
- 2). 1.0 gm  $\text{KNO}_3$  /1L d- $\text{H}_2\text{O}$
- 3). 1.0 gm  $\text{K}_2\text{HPO}_4$  /1L d- $\text{H}_2\text{O}$
- 4). 50.0 g disodium EDTA; 31.0 g  $\text{KOH}$  /1L d- $\text{H}_2\text{O}$
- 5). 4.98 g  $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$  / 1L acidified  $\text{H}_2\text{O}$  (Acidified  $\text{H}_2\text{O}$  is 999.0 mL d- $\text{H}_2\text{O}$  + 1.0 mL concentrated  $\text{H}_2\text{SO}_4$ )

6). 11.42 g  $\text{H}_3\text{BO}_3$  /1L d- $\text{H}_2\text{O}$

7). 8.82 g  $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ ; 1.44 g  $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$ ; 0.71 g  $\text{MoO}_3$ ; 1.57 g  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ ;

0.49 g  $\text{CoNO}_3 \cdot 6\text{H}_2\text{O}$  /1L d- $\text{H}_2\text{O}$

FES preparation:

to 950 mL of d- $\text{H}_2\text{O}$  add:

20.0 mL of stock solutions 1, 2, and 3.

1.0 mL of stock solutions 4, 5 and 6.

2.0 mL of stock solution 7

1.0 g of bacto-peptone

2.0 gm of Oxoid Lab-Lemco Powder

5.0 g of sucrose

## Protocol

### Step 1.

Harvest *Chlorella* NC64A (or Pbi) cells from 4 day old cultures ( $1.2 - 2.0 \times 10^7$  cells/ml) by centrifugation at 4000 x g for 5 minutes.

 DURATION

00:05:00

### Step 2.

Re-suspend in MBBM (NC64A cells) or FES (Pbi cells) at concentration of  $8.6 \times 10^7$  cells/mL.

### Step 3.

Add *chlorella* virus to a multiplicity of infection (MOI) = 10 plaque-forming units (pfu)/cell.

### Step 4.

Add 3.1 ml of 37% formaldehyde into 40 ml centrifuge tubes and place them on ice.

### Step 5.

Sample infected *chlorella* cells (25 mL) into prepared centrifuge tubes with formaldehyde (the final formaldehyde concentration is 4%) and place on ice.

### Step 6.

Centrifuge at 4000 x g for 5 minutes.

 DURATION

00:05:00

### Step 7.

Wash samples by re-suspending them in 10 mL of MBBM amended with 50 mM EDTA and following centrifugation at 4000 x g for 5 minutes. Repeat wash step 3 times.

#### **Step 8.**

Re-suspend washed infected cells in 0.5 mL of SB.

#### **Step 9.**

Add to the cells 0.5 mL of 2% low melting point agarose (BioRad) in SB (kept at 45°C), mix well (work quickly, try not to generate any air bubbles), and pour the mix into BioRad plug molds.

#### **Step 10.**

Place plug molds in refrigerator for 15 minutes to solidify.

 **DURATION**

00:15:00

#### **Step 11.**

Carefully remove agarose blocks from mold and place them into 2 mL of DB amended with 1mg/mL Proteinase K.

#### **Step 12.**

Incubate agarose blocks for 24 hrs at 50°C

#### **Step 13.**

After incubation, wash agarose blocks for 30 minutes with DB 4 times. Cut blocks in small pieces to fit gel wells.

 **DURATION**

00:30:00

#### **Step 14.**

Prepare 1% agarose gel (PFGE grade) in 1× TAE buffer using BioRad casting stand.

#### **Step 15.**

Load agarose blocks into gel wells and seal them with melted (45°C) 1% low melting point agarose in running buffer.

#### **Step 16.**

The chromosomes of *Hansenula wingei* (1.05-3.13 Mbp), cat#170-3667; *Schizosaccharomyces* (3.5-5.7 Mbp), cat#170-3633 (Bio-Rad, Hercules, CA, USA), and Yeast Chromosome PFG Marker (225-1,900 Kbp), cat#N0345S (New England Biolabs, Beverly, MA, USA) are used as molecular weight markers.

#### **Step 17.**

Separate chromosomal DNA in CHEF-DR (BioRad, Hercules, CA) electrophoresis unit with 1X TAE running buffer. Run electrophoresis at 3 V/cm (100 V) with pulse time ramping from 250 to 900 seconds for 60 hrs. Change buffer every 24 hrs.

 **DURATION**

00:00:00

#### **Step 18.**

Stain gel with 0.5 mg/L ethidium bromide for 20.

 **NOTES**

**Irina Agarkova** 19 Apr 2016

Alternatively, Sybr-Gold (Molecular Probes, Eugene, OR).