

DNA extraction from gram negative bacteria for ONT MinION sequencing

Louise Judd

Abstract

This is the protocol I have developed to generate High'ish' Molecular Weight gDNA from gram negative bacterial species to use for ONT MinION sequencing.

This protocol was initially developed for the high throughput extraction of bacterial gDNA for Illumina short read sequencing. To perform our high throughput work we utilise Beckman Coulter Biomek liquid handling robots and the Beckman Coulter gDNA extraction kit Genfind v2. This kit utilises magnetic particles and is relatively gentle. We have found that gDNA extracted using this method has a high'ish' molecular weight. On a fragment analyser the gDNA is typically a clean single peak >60 kb (the limits of the fragment analyser). ONT MinION libraries generated with this gDNA typically have N50s in excess of 20 kb and the data is suitable for the assembly of almost bacterial genomes we have tested to date. This protocol will not give you the ultra long 'whale' reads that others have reported but it is a good compromise between HMW DNA, speed and ease of use and capacity to scale up to high throughput.

NOTE

We have found that gDNA extracted using this method often have underrepresentation of smaller plasmids (<20kbp). We are not sure exactly why this happens but possible reasons are; the circular replicons do not bind as readily to the magnetic particles and so are excluded from the extraction or the plasmids remain circular and so are not available to adapter ligation in the library prep. I will upload an additional protocol for gDNA extraction that focuses on enrichment of plasmid sequences.

Citation: Louise Judd DNA extraction from gram negative bacteria for ONT MinION sequencing. **protocols.io** dx.doi.org/10.17504/protocols.io.p5mdq46

Published: 15 May 2018

Protocol

Bacterial culturing

Step 1.

- 1) Inoculate single bacterial colony into 2 ml of appropriate growth media
- 2) Incubate overnight 37°C with shaking
- 3) Transfer 1.5 ml of overnight culture into microfuge tube
- 4) Pellet bacteria
- 5) Remove supernatant

NOTE

At this stage bacterial pellet can be stored at -80°C for gDNA extraction at a later date with no compromise to extraction

Bacterial Lysis

Step 2.

- 1) Dissolve lyophilised RNase A (American Bioanalytical #AB14043) in 50% glycerol and 10 mM Tris (pH 8.0) to a final concentration of 100 mg/ml
(stock can be stored at -20°C for future use)
- 2) Dissolve lyophilised Genfind proteinase K in Genfind proteinase K buffer to a final concentration of 96 mg/ml
(stock can be stored at -20°C for future use)
- 3) For each bacterial isolate to be extracted make a lysis master mix containing 400 µl Genfind lysis buffer, 9 µl Genfind proteinase K and 1 µl RNase A (see step 1)
- 4) Add 400 µl of the lysis master mix to each bacterial pellet
- 5) Gently tip mix 10 times to ensure the bacterial pellet is well resuspended
- 6) Incubate at 37°C for 30 minutes
- 7) Proceed to DNA extraction

DNA extraction

Step 3.

- 1) Invert the Genfind binding buffer 20 times to ensure complete resuspension of the magnetic particles
- 2) Add 300 µl of the Genfind binding buffer to each of the samples and gently tip mix 10 times, avoid

generating bubbles.

The magentic particles will most likely clump together at this stage due to the high concentration of DNA in the sample.

- 3) Incubate samples at room temperature for 5 minutes
- 4) Place the sample on a magentic rack for 5 minutes or until solution completely clears
- 5) With the tubes on the magnet remove the supernatant avoiding the magentic particles
- 6) Remove samples from the magnet
- 7) Add 800 ul of Genfind wash buffer 1 to each of the samples
- 8) Gently tip mix each sample 10 times.

The idea is to try and resuspend the magentic beads but this will most likely not be possible as they will clump together.

I try to make sure the clumps of beads move in and out of the pipette tip with each mix

- 9) Place the sample on a magentic rack for 5 minutes or until solution completely clears
- 10) Repeat steps 7-10 for a total of two washes with Genfind wash buffer 1
- 11) Add 500 ul of Genfind wash buffer 2 to each of the samples
- 12) Gently tip mix each sample 10 times.

The idea is to try and resuspend the magentic beads but this will most likely not be possible as they will clump together.

I try to make sure the clumps of beads move in and out of the pipette tip with each mix

- 13) Place the sample on a magentic rack for 5 minutes or until solution completely clears
- 14) Repeat steps 11-13 for a total of two washes with Genfind wash buffer 2
- 15) Briefly centrifuge samples to bring all liquid to the bottom of the tube
- 16) Place the sample on a magentic rack for 5 minutes or until solution completely clears
- 17) Remove all traces of wash buffer 2 with a pipette, but do not allow beads to fully dry out.
- 18) Add 50-200 ul of 10 mM Tris pH 8.0 to each sample
(I usually start with 100 ul but it will depend on the groth dynamics of your baccteria)
- 19) Remove samples from the magnet
- 20) Gently tip mix each sample 10 times and incubate at 37°C for 5 minutes then tip mix again.

To ensure elution of DNA off the particles it is crucial that the magnetic particles are completely resuspended.

If it is not possible to entirely resuspend particles then add more 10 mM Tris.

21) Place the sample on a magnetic rack for 5 minutes or until solution completely clears

22) Transfer supernatant to a clean microfuge tube. Try to avoid carryover of magnetic particles into the supernatant.

NOTE

After extraction DNA should be stored at 4°C to minimise any damage. Storage at -20°C will result in ice crystal formation that will shear your DNA.

I try to do my DNA extractions on the same day as my library prep but have stored DNA at 4°C for a couple of weeks and still made successful libraries.

DNA QC

Step 4.

1) Visualise sample to estimate the molecular weight. We have a Fragment Analyzer and run a gDNA 1-60kbp gel.

Once you have extracted gDNA from the same species a number of times it is possibly not necessary to run a "gel" with each subsequent extraction. However, I would run QC gels each time you tackle a new bacterial species

2) Quantify your sample with a fluorimeter e.g. QUANTUS or Qubit.

3) Analyse 1 µL in a UV spectrophotometer (e.g. Nanodrop).