

Mar 28, 2019

Working

## P Arabidopsis flower dip transformation

Forked from Arabidopsis flower dip transformation

## Magdalena Julkowska<sup>1</sup>

<sup>1</sup>King Abdullah University of Science and Technology

dx.doi.org/10.17504/protocols.io.pemdjc6



**ABSTRACT** 

A standard protocol how to transform Arabidopsis plants using a flower-dip method.

For higher rates of transformation, plants may be dipped two times at 5-7 day intervals. We suggest one dip two days after clipping, and a second dip one week later.

PROTOCOL STATUS

## Working

We use this protocol in our group and it is working

SAFETY WARNINGS

- Grow healthy Arabidopsis plants until they are flowering.
- 2 Prepare Agrobacterium tumefaciens strain carrying gene of interest on a binary vector. Grow a large liquid culture @ 28C in LB with antibiotics to select for the binary plasmid, or grow in other media.
- 3 Spin down Agrobacterium, resuspend to OD600 = 0.8-1.0 (can be higher or lower) in LB with 5% Sucrose. You will need 250 ml for each four small pots to be dipped.
- 4 Before dipping, add Silwet L-77 to a concentration of 0.05% (125 ul/250mL) and mix well.
- 5 Dip above-ground parts of plant in Agrobacterium solution for 10-30 seconds, with gentle agitation. You should then see a film of liquid coating plant.
- 6 Place dipped plants under a cover in a dark O/N to maintain high humidity (plants should be laid on their side).
- 7 Water and grow plants normally. Stop watering as seeds become mature.
- 8 Harvest dry seed. Transformants are usually all independent, but are guaranteed to be independent if they come off of separate plants.

- 9 Select for transformants using antibiotic or herbicide selectable marker.
- 10 Transplant transformants to soil.

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited