

# **RNA Imaging with MERFISH - Probe Construction**

# Jeffrey R. Moffitt, Xiaowei Zhuang

# **Abstract**

Staining hundreds to thousands of RNA species each of which requires tens of unique encoding probes requires a very large number of unique oligonucleotide sequences. The sheer number of unique sequences makes traditional solid phase oligonucleotide synthesis prohibitively expensive in most cases. For example, assuming a modest cost of 0.10/base, an encoding probe length of 0.10 nt, and the need for 0.10000 unique oligos for a single MERFISH measurement, the cost of the needed oligos would be 0.10001. To circumvent this astronomical cost, we have developed a high-throughput approach to generating these probes which utilizes array-derived synthesis of complex oligonucleotide pools. These arrayderived complex oligonucleotide pools that contain 0.100000 custom designed sequences can be purchased for only a few thousand dollars. The challenge to using these pools is that each individual sequence is provided in quantities far too small to be used directly for labeling. Thus, we developed an enzymatic amplification protocol to generate the encoding probes in high quantity sufficient for RNA FISH experiments using the array-derived complex oligonucleotide pools as templates.

Citation: Jeffrey R. Moffitt, Xiaowei Zhuang RNA Imaging with MERFISH - Probe Construction. protocols.io

dx.doi.org/10.17504/protocols.io.meqc3dw

Published: 29 Mar 2018

# **Guidelines**

The basic protocol involves four steps, which are described in detail in the following sections. First, the oligopool is amplified via PCR to create template molecules for in vitro transcription. Second, an in vitro transcription produces RNA from these template molecules as well as another ~100-fold amplification. Third, a reverse transcription reaction generates single-stranded DNA from this RNA. Fourth, the RNA template is removed via alkaline hydrolysis. We utilize an RNA intermediate for two reasons: i) in vitro transcription can produce very large quantities of nucleic acid in very small volumes, reducing the amount of material that must be used, and ii) the alkaline susceptibility of RNA allows it to be easily removed from the final DNA probes. Because this reaction involves the use of an RNA intermediate, we recommend using the higher standards of laboratory cleanliness often required for handling RNA. Specifically, all surfaces and pipettes should be cleaned daily using RNase removal solutions and separate stocks of all buffers should be kept solely for RNA work.

### **Equipment:**

1. Table top centrifuge

- 2. qPCR machine or thermocycler
- 3. 37 °C incubator or water bath
- 4. 50 °C water bath
- 5. 95 °C water bath
- 6. Vacuum manifold (optional)
- 7. Gel electrophoresis equipment for poly-acrylamide gels (optional)
- 8. Vacuum concentrator (optional)

### The protocol workflow is as follows:

### Amplification of in vitro Template (Steps 1-8)

The first step in this protocol is to use PCR to amplify template molecules for the in vitro transcription. We recommend running this reaction as a limited-cycle PCR, i.e. monitor the status of the reaction in real time with a qPCR machine and remove the samples immediately before the final amplification plateau. We recommend limiting the number of PCR cycles because we have found that over amplification of complex libraries can produce molecules that miss-prime on other molecules, forming long concatemers that both reduce the yield of proper encoding probes and which could produce spurious signals.

### In vitro Transcription (Steps 9-14)

The second step of this protocol is a high yield in vitro transcription reaction that further amplifies the template molecules created in Section **Amplification of** *in vitro* **Template** as well as converts them into RNA.

### Reverse Transcription and Purification of Encoding Probes (Steps 15-25)

In this step of the protocol, the large quantities of RNA produced by the high yield in vitro transcription are converted to single-stranded DNA using a reverse transcription reaction. This RNA template is then removed via alkaline hydrolysis, and the final encoding probes are purified and concentrated.

### **Materials**

20X EvaGreen 31000 by Biotium

2X Phusion hot start polymerase master mix M0536S by New England Biolabs

Tris-EDTA (TE) pH 8 buffer AM9849 by Ambion

DNA binding buffer D4004-1-L by Zymo Research

DNA wash buffer C1016-50 by Zymo Research

Oligo binding buffer D4060-1-40 by Zymo Research

100-μg capacity silicon columns (Spin-V) D4003-2-48 by Zymo Research

RNA binding buffer (Optional) R1013-2-100 by Zymo Research

RNA prep buffer (Optional) R1060-2-100 by Zymo Research

RNA wash buffer (Optional) R1003-3-24 by Zymo Research

Quick HiScribe T7 polymerase kit E2050S by New England Biolabs

RNasin plus N2611 by Promega

Maxima H- reverse transcriptase EP0751 by Thermo Scientific

10 mM mix of dNTPs N0447S by New England Biolabs

0.5 M EDTA AM9261 by Ambion

1 N NaOH JT5635-2 by Vwr

Nuclease-free water AM9932 by Ambion

100% Ethanol (KOPTEC) 89125-186 by Vwr

D/RNAaseFree 47751-044 by Vwr

1.5 mL LoBind tubes 022431021 by Eppendorf

PCR tubes by Contributed by users

### **Protocol**

# Amplification of in vitro Template

# Step 1.

Design the primers. This PCR will not only amplify the library, it will also add the T7 promoter to these molecules to allow in vitro transcription of these templates. The sequence of the forward primer is the same as that designed in the protocol for Design of Oligonucleotide Probes. However, the sequence of the reverse primer, must include a T7 promoter, TAATACGACTCACTATAGGG, at the 5' end. Example primers can be seen in Table 3. The forward primer will also be used as the primer for reverse transcription in Reverse Transcription and Purification of Encoding Probes below; thus, it is recommended to order this primer at a relatively large synthesis scale, such as 100 nmol or 250 nmol. The T7 promoter ends in a G triplet, and the presence of a G quadruplet, i.e. four Gs in a row, often significantly lowers synthesis yields; thus, we recommend removing any 5'-terminal G nucleotides in the sequence of the reverse primer. The presence of the terminal G nucleotides in the T7 promoter region will serve as replacements for these nucleotides in the priming region. Resuspend forward primer to 200  $\mu$ M and the reverse primer to 100  $\mu$ M both in TE.

#### NOTES

#### Max Brown 29 Dec 2017

#### Table 3

Example primers. These primers are compatible with the encoding probe template listed in Table 2. Note that the 5' G in the reverse primer has been removed so as not to create a G quadruplet with the terminal GGG of the T7 promoter.

Forward primer	CGCGGGCTATATGCGAACCG
Reverse primer with T7 promoter	$\mathit{TAATACGACTCACTATAGGGCCGTGGAGGGCATACAACGCC}$

# Amplification of in vitro Template

# Step 2.

Prepare the PCR. In a 1.7 mL Eppendorf tube, mix the following: 40  $\mu$ L 20X Eva Green; 2  $\mu$ L 200  $\mu$ M forward primer; 4  $\mu$ L 100  $\mu$ M reverse primer; 1  $\mu$ L of 80  $ng/\mu$ L complex oligopool; 353  $\mu$ L nuclease free water; 400  $\mu$ L 2X Phusion hot start polymerase master mix. Aliquot 50  $\mu$ L volumes into 16 PCR tubes.

**■** AMOUNT

1 μl Additional info: 80 ng/μL complex oligopool

AMOUNT

353 µl Additional info: nuclease free water

AMOUNT

4 μl Additional info: 20X Eva Green

AMOUNT

2 μl Additional info: 200 μM forward primer

AMOUNT

4 μl Additional info: 100 μM reverse primer

AMOUNT

1 μl Additional info: 80 ng/μL complex oligopool

AMOUNT

353 µl Additional info: nuclease free water

**■** AMOUNT

400 µl Additional info: 2X Phusion hot start polymerase master mix

# Amplification of in vitro Template

# Step 3.

Amplify the template. Run the following protocol on a qPCR machine: 1) 98 °C for 3 minutes; 2) 98 °C for 10 s; 3) 65 °C for 10 s; 4) 72 °C for 15 s; 5) Measure the fluorescence of each sample. Repeat cycle steps 2 through 5 until the rate at which the sample amplifies decreases, which is a sign that it is approaching the final amplification plateau. Due to the complexity of these oligonucleotide pools, it is very unlikely that, once denatured, each molecule will find and rehybridize to its complement as opposed to partially hybridize with the common priming regions of a different molecule; thus, it is recommended that samples be removed after the elongation step—while the instrument is at 72 °C—and before it reaches the 98-°C-melting step of the next cycle. If a qPCR machine is not available, we recommend determining the appropriate number of cycles to run by quantifying the yield of small-scale PCR reactions run for different number of cycles.

# Amplification of in vitro Template

### Step 4.

Purify the template. We utilize column purification to remove enzyme, nucleotides, and primers. In a 15 mL Falcon tube, mix the following:  $800~\mu\text{L}$  of the PCR reaction generated in Step 3; and 4 mL of DNA binding buffer. Pull this mixture across a  $100-\mu\text{g}$  capacity column using either a vacuum manifold or a centrifuge.

**■** AMOUNT

800 µl Additional info: PCR reaction

**■** AMOUNT

4 ml Additional info: DNA binding buffer

### Amplification of in vitro Template

# Step 5.

Wash the column with 300  $\mu$ L DNA wash buffer, spinning the column in a table top centrifuge at maximum speed for 30 s. (1/2)

**■** AMOUNT

300 µl Additional info: DNA wash buffer

# Amplification of in vitro Template

# Step 6.

Wash the column with 300  $\mu$ L DNA wash buffer, spinning the column in a table top centrifuge at maximum speed for 30 s. (2/2)

**■** AMOUNT

300 µl Additional info: DNA wash buffer

### Amplification of in vitro Template

# Step 7.

Elute the template by adding 170  $\mu$ L nuclease-free water to the column, transferring the column to a fresh 1.7 mL Eppendorf tube, and spinning at maximum speed for 30 s. Set aside 10  $\mu$ L of this reaction for quality control.

**AMOUNT** 

170 µl Additional info: Nuclease-free water

# Amplification of in vitro Template

# Step 8.

(Optional) Quality control for template reaction. Two important quality control steps can be performed at this point. First, it is useful to measure the concentration of the template with a spectrophotometer, such as the Nanodrop. The concentration should be between 10 ng/ $\mu$ L to 50 ng/ $\mu$ L. The second quality control step is gel electrophoresis and will be described in Section 'Reverse Transcription and Purification of Encoding Probes' below.

### In vitro Transcription

# Step 9.

in vitro transcription. In a fresh 1.7 mL Eppendorf tube, mix the following: 160  $\mu$ L of the in vitro template created above; 176  $\mu$ L of nuclease free water; 250  $\mu$ L of the NTP buffer mix provided with the Quick HiScribe T7 polymerase kit; 25  $\mu$ L of RNasin Plus; and 25  $\mu$ L T7 polymerase (from the same

HiScribe kit). Incubate the reaction in a 37 °C incubator for 12–16 hours. Often the reaction is complete after 6–8 hours, but it is typically convenient to leave this reaction overnight. Remove 20  $\mu$ L for quality control.

**■** AMOUNT

160 µl Additional info: in vitro template

**■** AMOUNT

176 µl Additional info: nuclease free water

AMOUNT

250 µl Additional info: NTP buffer mix

**■** AMOUNT

25 µl Additional info: RNasin Plus

**■** AMOUNT

25 μl Additional info: T7 polymerase

37 °C Additional info: Incubation

# In vitro Transcription

# Step 10.

(Steps 10-14 are optional) Quality control for the *in vitro* transcription. To confirm that the *in vitro* transcription was successful, purify the reaction and measure its concentration with a spectrophotometer. To purify, mix 20  $\mu$ L of the *in vitro* reaction with 30  $\mu$ L nuclease-free water, 100  $\mu$ L RNA binding buffer, and 150 $\mu$ L 100% ethanol. Pass across a 100-ug-capacity spin column in a table top centrifuge.

**■** AMOUNT

20 μl Additional info: in vitro reaction

AMOUNT

30 µl Additional info: nuclease-free water

**■** AMOUNT

100 μl Additional info: RNA binding buffer

AMOUNT

150 µl Additional info: 100% ethanol

# In vitro Transcription

### **Step 11.**

Wash this column once with 400  $\mu L$  RNA prep buffer with a 30 s spin at the maximum speed of the table top centrifuge.

**■** AMOUNT

400 μl Additional info: RNA prep buffer

# In vitro Transcription

# Step 12.

Wash column with 200  $\mu$ L RNA wash buffer with a 30 s spin at the maximum speed of the table top centrifuge. (1/2)

AMOUNT

200 µl Additional info: RNA wash buffer

### In vitro Transcription

# Step 13.

Wash column with 200  $\mu$ L RNA wash buffer with a 30 s spin at the maximum speed of the table top centrifuge. (2/2)

# **AMOUNT**

200 µl Additional info: RNA wash buffer

### In vitro Transcription

# Step 14.

Elute the RNA with 100  $\mu$ L nuclease-free water. If successful, the concentration of the in vitro transcription should be between 0.5  $\mu$ g/ $\mu$ L to 2  $\mu$ g/ $\mu$ L. Purified RNA can also be run on a gel as described below (Step 24).

# **■** AMOUNT

100 μl Additional info: nuclease-free water

Reverse Transcription and Purification of Encoding Probes

# **Step 15.**

Reverse transcription. To the unpurified in vitro transcription created in the Steps above, add the following and mix well: 200  $\mu$ L 10 mM dNTP mix; 120  $\mu$ L 200  $\mu$ M forward primer; 240  $\mu$ L 5X Maxima buffer; 24  $\mu$ L RNasin Plus; 24  $\mu$ L Maxima H- reverse transcriptase.

# **■** AMOUNT

200 ul Additional info: 10mM dNTP mix

AMOUNT

120 μl Additional info: 200 μM forward primer

**■** AMOUNT

240 µl Additional info: 5X Maxima buffer

AMOUNT

24 µl Additional info: RNasin Plus

AMOUNT

24 μl Additional info: Maxima H- reverse transcriptase Reverse Transcription and Purification of Encoding Probes

# **Step 16.**

Incubate in a 50°C water bath for 1 hour.

### **▮** TEMPERATURE

50 °C Additional info: Water bath

#### NOTES

Max Brown 27 Dec 2017

It is important to use a water bath, not an air incubator, to insure that the temperature of the sample rises to 50°C quickly.

### Reverse Transcription and Purification of Encoding Probes

### **Step 17.**

Alkaline hydrolysis. Split the above reaction into two 1.7-mL Eppendorf tubes and add the following to each: 300  $\mu$ L 0.5 M EDTA and 300  $\mu$ L 1 N NaOH. Incubate in a 95°C water bath for 15 minutes.

AMOUNT

300 µl Additional info: 0.5 M EDTA

**■** AMOUNT

300 µl Additional info: 1N NaOH

**AMOUNT** 

300 µl Additional info: 0.5 M EDTA

**■** AMOUNT

300 µl Additional info: 1 N NaOH

**■ TEMPERATURE** 

95 °C Additional info: Water bath

Reverse Transcription and Purification of Encoding Probes

### **Step 18.**

Purification of ssDNA probe. Combine the two aliquots above into a single 50 mL Falcon tube and add the following: 4.8 mL Oligo binding buffer and 19.2 mL 100% ethanol. Mix well and split equally between eight 100-µg capacity spin columns. Pull the sample across the columns with a vacuum manifold or via centrifugation.

**■** AMOUNT

4.8 ml Additional info: Oligo binding buffer

**■** AMOUNT

19.2 ml Additional info: 100% ethanol

Reverse Transcription and Purification of Encoding Probes

Step 19.

Wash the columns once with 750 µL DNA wash buffer.

**■** AMOUNT

750 µl Additional info: DNA wash buffer

Reverse Transcription and Purification of Encoding Probes

Step 20.

Elute the columns using 100  $\mu$ L of nuclease-free water. Combine eluates and set aside 10  $\mu$ L for quality control.

**■** AMOUNT

100 μl Additional info: Nuclease-free water

Reverse Transcription and Purification of Encoding Probes

Step 21.

Concentration of probe. Use a vacuum concentrator to dry the samples.

NOTES

Max Brown 27 Dec 2017

This process could take several hours.

### Reverse Transcription and Purification of Encoding Probes

### Step 22.

Resuspend the dried pellet in 24  $\mu$ L of nuclease-free water, or if desired, the hybridization buffer described in the attached full manuscript (Section 5.1).

**■** AMOUNT

24 µl Additional info: Nuclease-free water

Reverse Transcription and Purification of Encoding Probes

# Step 23.

Store probe at -20 °C and avoid unnecessary freeze-thaw cycles. If a vacuum concentrator is not available, it is also possible to concentrate the probe using ethanol precipitation.

### **▮** TEMPERATURE

-20 °C Additional info: Storage

Reverse Transcription and Purification of Encoding Probes

# Step 24.

(Optional) Quality control of in vitro template, RNA, and probe. We recommend running the in vitro template, the RNA, and the final probe on a 15% TBE-urea polyacrylamide gel to identify both RNase contamination and low conversion of the reverse transcription primer to full length probe. Large smearing both in the RNA band and the probe band can indicate RNase contamination. Failure to efficiently convert the reverse transcription primer into probe is revealed by a bright band running at the 20-nt length corresponding to the primer. Increasing the amount of RNA template in the reverse transcription often improves the fraction of primer converted into probe. We routinely obtain 75% or greater incorporation of the reverse transcription primer into probe with the above protocol.