

# Protocol for STO Cell Transfection by FuGENE HD

Misha Gurevich, V. Katerov

## Abstract

Protocol for Transfection Mouse Embryonic Stem Cells.

**Citation:** Misha Gurevich, V. Katerov Protocol for STO Cell Transfection by FuGENE HD. **protocols.io**

dx.doi.org/10.17504/protocols.io.e8jbhun

**Published:** 23 Jun 2016

## Protocol

### Step 1.

#### Cell plating

STO cells were seeded the day before transfection with the density 15,000 cells per well in 100  $\mu$ l complete growth medium DMEM+10% Fetal Bovine Serum.

### Step 2.

Prepare 0.02 $\mu$ g/ $\mu$ l pCMV $\beta$  plasmid DNA solution in OptiMEM®.

### Step 3.

Add 6 $\mu$ l of reagent to 100  $\mu$ l of OptiMEM® /DNA solution.

### Step 4.

Mix carefully by pipetting (10-15 times).

### Step 5.

Incubate 5 min at room temperature.

### Step 6.

Add 5 $\mu$ l complex per well to the cells, and mix thoroughly.

### Step 7.

Place the cells into CO2 incubator for 26-28 hours.

### Step 8.

Remove the medium from the well and wash the cells once with 100 $\mu$ l per well PBS.

### Step 9.

Fix the cells in the well with 50 $\mu$ l solution of 4% formaldehyde in PBS for 5min at room temperature.

### Step 10.

Wash each well twice with 100 $\mu$ l PBS.

### Step 11.

Add 50µl per well of substrate/stain solution and incubate the plate overnight at 37°C.

**Step 12.**

Observe the cells under microscope and evaluate the proportion of blue (β-gal-positive) cells.