

Mate-pair large libraries preparation for assembly of the *Lateolabrax maculatus* genome

Chang Li

Abstract

This protocol is used to clarify the process of the mate-pair large libraries preparation for the *L. maculatus*.

Citation: Chang Li Mate-pair large libraries preparation for assembly of the *Lateolabrax maculatus* genome. **protocols.io**

dx.doi.org/10.17504/protocols.io.ss2eege

Published: 23 Aug 2018

Protocol

Genomic DNA interruption

Step 1.

The genomic DNA was fragmented using a Covaris E220 ultrasonicator (Covaris, Brighton, UK) to obtain 2kb (for 2 kb library) and a Hydroshear (GeneMachines, CA, USA) to obtain 5 kb, 10 kb and 20kb fragments (for 5 kb, 10 kb and 20 kb libraries respectively).

EQUIPMENT

Equipment brand:

Covaris

SKU:

1

Specifications:

E220 ultrasonicator

EQUIPMENT

Equipment brand:

GeneMachines

SKU:

1

Specifications:

Hydroshear

End-repair


Step 2.

Repaired using T4 DNA polymerase, (ENZYMATICS, Beverly, the U.S.) 30 min at 20 °C.

TEMPERATURE

20 °C :

REAGENTS

 T4 DNA polymerase by [Enzymatics](#)

 DURATION

00:30:00 :

Biotin Label

Step 3.

Add Biotin Label by Biotin dNTP Mix (5mM) 30 min at 20 °C.

 TEMPERATURE

20 °C :

 REAGENTS

 Biotin dNTP Mix by [Invitrogen - Thermo Fisher](#)

 DURATION

00:30:00 :

Fragment selection

Step 4.


These fragments were further selected into size ranges of 2-2.4 kb, 5-5.5 kb, 10-11 kb or 20-23 kb by agarose gel electrophoresis.

Fragment cyclizing

Step 5.

The T4 and T3 DNA ligase were used to connect the ring. And then, Covaris LE220 was used to cyclizing DNA fragments.

 REAGENTS

 T3 DNA ligase by [Enzymatics](#)

End-repair

Step 6.

Fragmented DNA labeled with biotin was captured on M280 streptavidin beads (Invitrogen, CA, USA), followed by end repair (30 min. at 20°C, 1000 rotation per minute, rpm, vibrate for 15 sec. per 2 min.), A-tailing (30 min. at 37°C, 1000 rpm vibrate for 15 sec. per 2 min.).

 TEMPERATURE

20 °C :

 TEMPERATURE

37 °C :

 REAGENTS

 M280 streptavidin beads by [Invitrogen - Thermo Fisher](#)

 DURATION

00:30:00 :

 DURATION

00:00:15 :

 DURATION

00:30:00 :

 DURATION

00:00:15 :

Add adapter

Step 7.

Adaptor ligation (1h at 20°C, 1000rpm vibrate for 15 sec per 2 min.).

 TEMPERATURE

20 °C :

 DURATION

01:00:00 :

 DURATION

00:00:15 :

PCR amplification

Step 8.

PCR amplifications on beads 95°C 3 min., (98 °C 20 sec., 60 °C 15 sec., 72 °C 45 sec.) for N cycles, 72 °C 10 min., 4°C hold (For 2 kb library, N=16; For 5 kb library, 10 kb library and 20 kb library, N=18)] using Enzymatics (MA, USA) and NEB (MA, USA) reagent.

 TEMPERATURE

95 °C :

 DURATION

00:03:00 :