



# Isolation, culturing, and cryopreservation of *Endozoicomonas* (Gammaproteobacteria: Oceanospirillales: Endozoicomonadaceae) from reef-building corals

Claudia Pogoreutz<sup>1</sup>, Christian Voolstra<sup>1</sup>

Working

Sep 30, 2018

<sup>1</sup>King Abdullah University of Science and Technology

[dx.doi.org/10.17504/protocols.io.t2aeqae](https://doi.org/10.17504/protocols.io.t2aeqae)

reefgenomics , Aiptasia-Symbiodinium Model System



Christian Voolstra

King Abdullah University of Sc...



## ABSTRACT

*Endozoicomonas* are gram-negative bacteria widely and often abundantly associated with marine invertebrates and fish (Yang *et al.*, 2010; Bayer *et al.*, 2013; Nishijima *et al.*, 2013; Hyun *et al.*, 2014; Katharios *et al.*, 2015; Ding *et al.*, 2016; Neave *et al.*, 2017a; Schreiber *et al.*, 2016; Pogoreutz *et al.*, 2018). Despite their ubiquitous distribution, only few cultured strains are available, as *Endozoicomonas* are supposedly difficult to isolate and to maintain in pure cultures (Neave *et al.*, 2017b). Here we detail a protocol that allowed us to reproducibly isolate *Endozoicomonas* from stony corals (*Acropora humilis*) from the Red Sea. This protocol should be useful in isolating abundant *Endozoicomonas* from other corals or marine invertebrates.

## References

- Bayer, T., Neave, M.J., Alsheikh-Hussain, A., Aranda, M., Yum, L.K., Mincer, T., et al. (2013) The microbiome of the Red Sea coral *Stylophora pistillata* is dominated by tissue-associated *Endozoicomonas* bacteria. *Appl. Environ. Microbiol.* **79**: 4759–4762.
- Ding, J.Y., Shiu, J.H., Chen, W.M., Chiang, Y.R., and Tang, S.L. (2016) Genomic insight into the host-endosymbiont relationship of *Endozoicomonas montiporae* CL-33T with its coral host. *Front. Microbiol.* **7**: 251.
- Hyun, D.W., Shin, N.R., Kim, M.S., Oh, S.J., Kim, P.S., Whon, T.W., and Bae, J.W. (2014) *Endozoicomonas atrina* sp. nov., isolated from the intestine of a comb pen shell *Atrina pectinata*. *Int. J. Syst. Evol. Microbiol.* **64**: 2312–2318.
- Katharios, P., Seth-Smith, H.M.B., Fehr, A., Mateos, J.M., Qi, W., Richter, D., et al. (2015) Environmental marine pathogen isolation using mesocosm culture of sharpnose seabream: Striking genomic and morphological features of novel *Endozoicomonas* sp. *Sci. Rep.* **5**: 1–13.
- Neave, M.J., Michell, C.T., Apprill, A., and Voolstra, C.R. (2017a) *Endozoicomonas* genomes reveal functional adaptation and plasticity in bacterial strains symbiotically associated with diverse marine hosts. *Sci. Rep.* **7**: 40579.
- Neave, M. J., Apprill, A., Ferrier-Pagès, C., & Voolstra, C. R. (2016) Diversity and function of prevalent symbiotic marine bacteria in the genus *Endozoicomonas*. *App. Microbiol. Biotech.*, **100**: 8315–8324.
- Neave, M.J., Rachmawati, R., Xun, L., Michell, C.T., Bourne, D.G., Apprill, A., and Voolstra, C.R. (2017b) Differential specificity between closely related corals and abundant *Endozoicomonas* endosymbionts across global scales. *ISMEJ* **11**: 186–200.
- Nishijima, M., Adachi, K., Katsuta, A., Shizuri, Y., and Yamasato, K. (2013) *Endozoicomonas numazuensis* sp. nov., a gammaproteobacterium isolated from marine sponges, and emended description of the genus *Endozoicomonas* Kurahashi and Yokota 2007. *Int. J. Syst. Evol. Microbiol.* **63**: 709–714.
- Pogoreutz, C., Rädcker, N., Cárdenas, A., Gärdes, A., Wild, C., and Voolstra, C.R. (2018) Dominance of *Endozoicomonas* bacteria throughout coral bleaching and mortality suggests structural inflexibility of the *Pocillopora verrucosa* microbiome. *Ecol. Evol.* **8**: 2240–2252.
- Schreiber, L., Kjeldsen, K.U., Funch, P., Jensen, J., Obst, M., Lopez-Legentil, S., and Schramm, A. (2016) *Endozoicomonas* are specific, facultative symbionts of sea squirts. *Front. Microbiol.* **7**: 1–15.
- Yang, C.-S., Chen, M.-H., Arun, A.B., Chen, C.A., Wang, J.-T., and Chen, W.-M. (2010) *Endozoicomonas montipora* sp. nov., isolated from the encrusting pore coral *Montipora aequituberculata*. *Int. J. Syst. Evol. Microbiol.* **60**: 1158–1162.



Endozoicomonas\_isolation\_protocol\_PogoreutzVoolstra.pdf

## PROTOCOL STATUS

### Working

We use this protocol in our group and it is working.

## GUIDELINES

For the Isolation of *Endozoicomonas*, several factors are important to consider: **homogenization, dilution, incubation temperature, replication, and colony morphology**:

(1) **Homogenization**: Thoroughly homogenize the coral tissue to make sure the *Endozoicomonas* aggregates are 'set free' from the host tissue (e.g., for 30 seconds at medium speed with an UltraTurrax or comparable tool). We didn't do that in a first attempt and all we isolated was *Vibrios*. Once we implemented the tissue homogenization, overall taxonomic diversity of isolates increased.

(2) **Dilution**: if the *Endozoicomonas* make up a high proportion of the microbiome (for instance, 60 - 90 % in Red Sea *Acropora humilis*), the dilution will help to select for abundant organisms. We have plated tissue slurry "undiluted" and in 1:10, 1:100, 1:1000 dilutions. The *Endozoicomonas* grew on the Marine Agar plates inoculated with 'undiluted' and 1:10 diluted slurries.

(3) **Temperature**: we settled for an incubation temperature low compared to ambient seawater temperatures (23°C in our case; the central Red Sea has an annual temperature average of 29°C). At this temperature, *Vibrios*, *Alteromonas*, and other fast-growing taxa would form colonies within the first 24 h, while the *Endozoicomonas* started to appear on day 4 post-inoculation, but would form many more colonies than other taxa.

(4) **Replication**: we had an agar plate replication of  $n = 3$  for each dilution of slurry and type of media (other media than Marine Agar were used, too, but no *Endozoicomonas* grew on them). This was not entirely necessary in our case because we had hundreds of *Endozoicomonas* colonies growing, but it might increase the chance of isolating *Endozoicomonas* strains (or other bacteria) that are not quite as abundant.

(5) **Colony morphology**: *Endozoicomonas* will form tiny round convex creamy colonies with entire margins. When plating the picked colonies, you will see that the average *Endozoicomonas* colony is rather sticky and will strongly adhere to loops or needles. These bacteria are certainly easier to handle in suspension culture than on plates. Once you have successfully isolated *Endozoicomonas* on plates and confirmed their identity with Sanger sequencing, you may transfer them in liquid culture (Difco 2216 Marine Broth).

- 1 Collect coral fragment(s) (approx. finger-sized) in a sterile zip-lock bag. Take notes on sampling conditions (site, sampling depth, habitat).
- 2 It is recommended to process coral fragment(s) right away. If not feasible, maintain at ambient reef water temperatures in flow-through aquaria (Temperature 28°C, salinity 40 PSU) or in closed aquarium systems with daily water exchange.
- 3 To obtain coral tissue slurry for inoculation, blast tissue off the coral skeleton using an air gun and autoclaved filtered seawater (AFSW; Whatman, 0.22 µm).

- 4 Homogenize a total volume of tissue slurry of no more than 15 ml per finger-sized fragment for 30 sec (e.g., using an IKA UltraTurrax).
- 5 Plate 50 ul of slurry on Marine Agar 2216 (MA; BD Difco) undiluted and in 1:10 dilution in triplicates. After incubation at 23 °C for 4 days, *Endozoicomonas* colonies will form and should be purified starting from a single colony.
- 6 After purification (minimum of 2 clean passages), confirm identity with Sanger sequencing.
- 7 Inoculate 8 ml Marine Broth 2216 (MB; BD Difco) with a 10 ul loop of cells from purified colonies on MA. Close the tube, vortex, and incubate under constant motion (e.g., 60 rpm) at 25 °C for 48 h.
- 8 Snap-freeze aliquots of *Endozoicomonas* strain in suspended culture as a 20% (v/v) glycerol suspension in MB and store at -80 °C or -140 °C.
- 9 Note: for the isolation of low abundance *Endozoicomonas* or other bacteria, (1) the amount of AFSW used should be kept to a minimum or no AFSW should be used at all to not unnecessarily dilute the slurry, (2) a selective medium might be required.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited