

Purification of nucleocapsid of tospoviruses from *Chenopodium quinoa* leaves

Tsung-Chi Chen

Abstract

Purification of nucleocapsid of tospoviruses from *Chenopodium quinoa* leaves

Citation: Tsung-Chi Chen Purification of nucleocapsid of tospoviruses from *Chenopodium quinoa* leaves. **protocols.io**
dx.doi.org/10.17504/protocols.io.iscceaaw

Published: 05 Jul 2017

Protocol

Step 1.

Leaves of tospovirus-infected *Chenopodium quinoa* (3-5 days post-inoculation with tospoviruses rely on daily observation) were homogenized with a blender (full speed) for 1 min with the TB buffer [0.01 M Tris-HCl (pH 8.0) containing 0.01 M sodium sulfite and 0.1% cysteine] (300 ml/100 g tissue).

Step 2.

After filtering through four layers of cheesecloth, the extract was centrifuged at 15300 *g* (10000 *rpm*) in a Beckman JA-14 rotor for 10 min.

Step 3.

The supernatant was stirred with 1-2% Triton X-100 (final concentration) at 4°C for 30 min.

Step 4.

Centrifugation at 79700 *g* (25000 *rpm*) in a Beckman 35Ti rotor for 2.5 hr through a 20% (w/w) sucrose cushion (sucrose within TB buffer, 15 ml/tube).

Step 5.

The pellets were then resuspended in the TBG buffer (TB buffer containing 0.01 M glycine) (8 ml/100 g tissue) further purified by isopycnic centrifugation through 35% (w/w) cesium sulfate (cesium sulfate within TBG buffer, 8-9 ml/each tube, and load 2 ml/each tube of samples) at 83500 *g* (30000 *rpm*) in a Beckman SW41 rotor for 18 hr.

Step 6.

The opalescent zones containing virus nucleocapsids were collected in 60 Ti centrifuge tube, fill full with TBG buffer, further centrifugation at 163600 *g* (45000 *rpm*) in a Beckman 70Ti rotor for 1 hr.

Step 7.

The pellets were resuspended in TBG buffer (3 ml/100 g tissue) and regarded as purified nucleocapsids.

Step 8.

Store at -20°C for further use.