



Protein preparation for LC-MS/MS analysis [↗](#)

Journal of Pain Research

Lihui Li¹

¹Department of Sport Medicine and Rehabilitation Center, Shanghai University of Sport, Shanghai, China; Rehabilitation Research Laboratory 2rLab, Department of Business Economics Health and Social Care, University of Applied Sciences and Arts of Southern Switzerland, Manno, Switzerland

[dx.doi.org/10.17504/protocols.io.upjevkn](https://doi.org/10.17504/protocols.io.upjevkn)



ABSTRACT

This protocol is with regard to quantitative proteomics analysis to identify biomarkers of chronic myofascial pain and therapeutic targets of dry needling in a rat model of myofascial trigger points. For further nanoscale liquid chromatography-tandem mass spectrometry (LC-MS/MS) analysis, the first step is to extract and prepare the protein of the sample.

EXTERNAL LINK

<https://doi.org/10.2147/JPR.S185916>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Li L, Huang Q, Barbero M, Liu L, Nguyen T, Beretta-Piccoli M, Xu A, Ji L, Quantitative proteomics analysis to identify biomarkers of chronic myofascial pain and therapeutic targets of dry needling in a rat model of myofascial trigger points. Journal of Pain Research doi: [10.2147/JPR.S185916](https://doi.org/10.2147/JPR.S185916)

PROTOCOL STATUS

Working

GUIDELINES

The animal experiments complied with the ARRIVE (Animal Research: Reporting of In Vivo Experiments) guidelines and were carried out in accordance with the U.K. Animals (Scientific Procedures) Act, 1986, and associated guidelines.

SAFETY WARNINGS

BEFORE STARTING

After the treatment was terminated on the 5th day, all rats were deeply anesthetized using intraperitoneal sodium pentobarbital at 60mg/kg, and the left GM tissue was quickly removed and rapidly frozen in liquid nitrogen for 5min and then immediately stored at -80°C refrigerator. The whole procedure for removing the GM in vivo was completed within 1 min from each rat, and air embolism was then administered.

- 1 The SDT lysis buffer was added to the sample and transferred to 2ml tubes with quartz sand consisting of 1/4-inch ceramic beads (MP 6540-424 for tissue samples).
- 2 The lysate was homogenized twice for 60 s using an MP homogenizer (24×2, 6.0M/S).
- 3 The homogenate was sonicated and then boiled for 15min.

- 4 After centrifuging at 14000 x g for 40min, the supernatant was filtered with 0.22µm filters. The filtrate was quantified with the BCA Protein Assay Kit (Bio-Rad Inc., Hercules, USA).
- 5 Protein samples of 20µg (each sample) were mixed with 5X loading buffer and the proteins were separated on 12.5% SDS-PAGE (sodium dodecyl sulfate-polyacrylamide gel electrophoresis) gel for 90min with a constant current of 14mA after boiling for 5min.
- 6 200µg of proteins (each sample) were added to a sample of 30µl of SDT buffer (150mM Tris-HCl, 4% SDS, 100mM DTT, pH8.0) and then removed from the dithiothreitol (DTT), the detergent and other low-molecular-weight components, using uric acid buffer (150mM Tris-HCl, 8M urea, pH8.0) and 100µl iodoacetamide to block the reduced cysteine residues for Filter-aided sample preparation (FASP).
- 7 The samples were then incubated for 30min in the dark. Then, the filters were washed in triplicate with 100µl uric acid buffer and washed twice in 100µl of 100mM tetraethyl-ammonium bromide (TEAB) buffer.
- 8 Finally, 4µg of trypsin in 40µl of TEAB buffer was used to digest the protein suspensions with overnight at 37°C.
- 9 Following sample digestion and immediately before use, the Tandem mass tag (TMT) labeled reagents were equilibrated to room temperature, and 4-plex TMT reagent was applied to label 100µg of the peptide mixture of each sample, according to the instructions of the manufacturer (Thermo Fisher Scientific Inc., Waltham, USA).
- 10 The protein peptides were labeled with TMT-126, TMT-127, TMT-128, and TMT-129 for the CG, MG, DG and NDG groups for three biological replicates.
- 11 Incubation for 60min was undertaken to allow the reaction to occur, and 8µl of 5% hydroxylamine was added to the sample, which was incubated for a further 15 minutes.
- 12 TMT-labeled samples were fractionated into 10 fractions using a Pierce high pH reversed-phase fractionation kit by an increasing acetonitrile step-gradient elution.
- 13 Before LC-MS/MS analysis, the fractions were dissolved in 0.1% formic acid, according to the instructions of the manufacturer (Thermo Fisher Scientific Inc.).



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited