



Nov 09,
2019

Transmission electron microscopy protocol for anaerobic ciliates

Johana Rotterova¹, Ivan Čepička¹

¹Charles University

1 Works for me [dx.doi.org/10.17504/protocols.io.85uhy6w](https://doi.org/10.17504/protocols.io.85uhy6w)

Johana Rotterova
Charles University

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Rotterová, J., Bourland, W. and Čepička, I., 2018. Tropidoatractidae fam. nov., a deep branching lineage of metopida (Armophorea, Ciliophora) found in diverse habitats and possessing prokaryotic symbionts. *Protist*, 169(3), pp.362-405. Beinart, R.A., Rotterová, J., Čepička, I., Gast, R.J. and Edgcomb, V.P., 2018. The genome of an endosymbiotic methanogen is very similar to those of its free-living relatives. *Environmental microbiology*, 20(7), pp.2538-2551.

MATERIALS

NAME	CATALOG #	VENDOR
Osmium tetroxide solution 2%	75633	Sigma Aldrich
25% Glutaraldehyde		
Acetone		
Sodium cacodylate trihydrate	View	ProSciTech
Ethanol		
Poly/Bed 812	08791-500	Polysciences

SAFETY WARNINGS

Please see Osmium Tetroxide safety data sheet here - [PrintMSDSAction.pdf](#)

Please see Sodium Cacodylate buffer safety data sheet here - <https://archive.proscitech.com/?navaction=pages&page=sds&catno=c0205>

Please see Glutaraldehyde safety data sheet here - [354400_SDS_ZA_EN.PDF](#)

BEFORE STARTING

Get a well growing culture of anaerobic ciliates following protocol listed here - [dx.doi.org/10.17504/protocols.io.85why7e](https://doi.org/10.17504/protocols.io.85why7e) or in Rotterová et al. 2018.

Rotterová J, Bourland W, Čepička I (2018). Tropidoatractidae fam. nov., a Deep Branching Lineage of Metopida (Armophorea, Ciliophora) Found in Diverse Habitats and Possessing Prokaryotic Symbionts.. *Protist*.
<https://doi.org/10.1016/j.protis.2018.04.003>

Fixation

- 1 Isolate 1 ml of thriving culture of ciliate cells, fix with 2.5% (v/v) glutaraldehyde (Polysciences) and subsequently centrifuge at 800×g for 5 min at 4 °C. 10m
- 2 Carefully replace supernatant with 1 ml of 2.5% (v/v) glutaraldehyde in 0.2M SCB (Sodium Cacodylate buffer, pH 7.2) and incubate 1 hour on ice. 1h 10m
- 3 Then, rinse the pellet carefully with SCB three times per 15 minutes. 50m
- 4 Postfix the cells with 1% OsO₄ in distilled water for 1 hour on ice. 3h 10m

Rinsing and embedding

- 5 Rinse the cells with distilled water (5 minutes), dehydrate them in a graded ethanol series from 30% to absolute ethanol: through 30%, 50%, 70%, 80%, 90%, 95%, and finally three times in absolute ethanol, each per eight minutes. 1h
- 6 Incubate cells two times in 1:1 ethanol-acetone mixture per 15 minutes.
- 7 Incubate cells once in 2 ml of 100% acetone per 15 minutes, and finally substitute acetone with acetone-resin mixture (EPON resin Poly/Bed 812, Polysciences) and incubate per 6-12 hours.
- 8 Substitute absolute EPON resin (Poly/Bed 812, Polysciences) three times and finally embed cells in the resin into capsules.
- 9 Let the embedded cells get polymerized at 70 °C for 48 hours. 48h 30m

Cutting ultrathin sections

- 10 Cut serial ultrathin sections on an Ultracut E ultramicrotome (Reichert) using a diamond knife and stain with lead citrate and uranyl acetate (2–3%). 3h

Protocol is adapted from Rotterova et al. (2018) and Beinart et al. (2018).



Rotterová J, Bourland W, Čepička I (2018). Tropidoatractidae fam. nov., a Deep Branching Lineage of Metopida (Armophorea, Ciliophora) Found in Diverse Habitats and Possessing Prokaryotic Symbionts.. Protist.

<https://doi.org/10.1016/j.protis.2018.04.003>



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author

and source are credited



Beinart RA, Rotterová J, Čepička I, Gast RJ, Edgcomb VP (2018). The genome of an endosymbiotic methanogen is very similar to those of its free-living relatives.. Environmental microbiology.
<https://doi.org/10.1111/1462-2920.14279>