



Mar 28,
2019

Working

Electrocompetent Agrobacterium transformation

Forked from [Electrocompetent Agrobacterium transformation](#)

[Magdalena Julkowska](#)¹

¹King Abdullah University of Science and Technology

[dx.doi.org/10.17504/protocols.io.pewdjfe](https://doi.org/10.17504/protocols.io.pewdjfe)



Magdalena Julkowska
King Abdullah University of Science and Technology 

ABSTRACT

Standard protocol for Agrobacterium transformation

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

SAFETY WARNINGS

- 1 Add a small amount of plasmid DNA (<100ng) to 40-80µl of competent cells in a tube on ice. Stir gently.
- 2 Pipette Agro cells with DNA into a chilled 0.2 cm electroporation cuvette. Shake cells to bottom of cell
- 3 Pulse cells at 2.5kV.
- 4 Add 1 ml YEP (of LB) media to cuvette, mix, and immediately transfer cells into sterile test tubes.
- 5 Allow cells to shake at 28°C for 2-4 hours.
- 6 Plate 10 µl and 100µl on selective media. Put the plates at 28°C.
- 7 Select positive cultures after 48h. If you use multiple antibiotics for selection, you might have to wait even longer.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited