

MBBM Stock Solutions Version 2

Dr. Steven Wilhelm, Samantha Coy

Abstract

Van Etten, J. (n.d.). Formulation of Modified Bold's Basal Medium (MBBM). Retrieved from <http://ncv.unl.edu/vanettenlab/>
Contact Dr. Steven Wilhelm (wilhelm@utk.edu) or Samantha Coy (srose16@vols.utk.edu) for additional information regarding this protocol.

Citation: Dr. Steven Wilhelm, Samantha Coy MBBM Stock Solutions. **protocols.io**

[dx.doi.org/10.17504/protocols.io.hgtb3wn](https://doi.org/10.17504/protocols.io.hgtb3wn)

Published: 29 Mar 2017

Protocol

Stock Solution #1

Step 1.

Add 25.0 g NaNO_3 per liter d- H_2O

Stock Solution #2

Step 2.

Add 2.5 g $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ per liter d- H_2O

Stock Solution #3

Step 3.

Add 7.5 g $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ per liter d- H_2O

Stock Solution #4

Step 4.

Add 7.5 g K_2HPO_4 per liter d- H_2O

Stock Solution #5

Step 5.

Add 17.5 g KH_2PO_4 per liter d- H_2O

Stock Solution #6

Step 6.

Add 2.5 g NaCl per liter d-H₂O

Stock Solution #7

Step 7.

Add 50.0 g disodium EDTA and 31.0 g KOH per liter d-H₂O

Stock Solution #8

Step 8.

Add 4.98 g FeSO₄*7H₂O per liter acidified H₂O

(Acidified H₂O is 999.0 mL d-H₂O + 1.0 mL concentrated H₂SO₄)

Stock Solution #9

Step 9.

Add 11.42 g H₃BO₃ per liter d-H₂O

Stock Solution #10

Step 10.

Add per liter d-H₂O

- 8.82 g ZnSO₄*7H₂O
- 1.44 g MnCl₂*4H₂O, 0.71 g MoO₃
- 1.57 g CuSO₄*5H₂O
- 0.49 g CoNO₃*6H₂O

Autoclave

Step 11.

Autoclave all stock solutions at 121°C for 20 min except solution #8 and #10 (stir to dissolve)