

Isolation of cell fractionation

Mengchen Shi

Abstract

Citation: Mengchen Shi Isolation of cell fractionation. protocols.io

dx.doi.org/10.17504/protocols.io.kcfcstn

Published: 18 Oct 2017

Protocol

Step 1.

Step 2.

A total of 5×106 cells were collected after brief trypsinization, followed by two more washes with PBS, and the cell pellet was resuspended in 200 μ l of mitochondria extraction buffer containing 0.02 mM phenylmethanesulfonyl fluoride (PMSF) and proteinase inhibitors (Keygentec, Nanjing, China).

Step 3.

After incubating on ice for 20 min, the cells were homogenized using a glass Dounce and pestle.

Step 4.

The homogenates were centrifuged at 600 g for 15 min at 4°C, and the resulting supernatant was collected and centrifuged at 11,000 g for 15 min at 4°C to separate the mitochondria (pellet) and cytoplasmic proteins (supernatant). The mitochondria pellet was lysed in mitochondria extraction buffer.