

Dec 20, 2019

## Chloral Hydrate Seed Clearing V.2

Gabrielle Sandstedt<sup>1</sup>, Andrea Sweigart<sup>1</sup><sup>1</sup>University of Georgia**1** Works for me dx.doi.org/10.17504/protocols.io.bakaicse

Mimulus

 Gabrielle Sandstedt 

### ABSTRACT

To characterize early seed development in *Mimulus* (1-5 days after pollination), we clear seeds with chloral hydrate and quickly obtain images of embryo and endosperm development using a Differential Interference Contrast Microscope.

### MATERIALS TEXT

Hoyer's Solution: 19g gum Arabic, 12g glycerol, 250g Chloral Hydrate, 75ml Water.

Diluted Hoyer's solution: 3 parts Hoyer's solution: 1 part 10% gum Arabic.

Forceps, glass microscope slides, coverslip, access to Differential Interference Contrast Microscope.

### SAFETY WARNINGS

Chloral Hydrate is DEA regulated and you either need a permit or access to a lab that has a permit with this substance.

This substance is also hazardous and acutely toxic. SDS: <https://www.caymanchem.com/msdss/21843m.pdf>

### BEFORE STARTING

Prepare Hoyer's solution: 19g gum Arabic, 12g glycerol, 250g Chloral Hydrate, 75ml Water.


Prepare a diluted Hoyer's solution: 3 parts Hoyer's solution: 1 part 10% gum Arabic.


### Prepare Hoyer's Solution

1 Hoyer's Solution:  **19 g** gum Arabic,  **12 g** glycerol,  **250 g** Chloral Hydrate,  **75 ml** Water.

 **10 ml** Diluted Hoyer's solution:  **7.5 ml** Hoyer's solution:  **2.5 ml** 10% gum Arabic

### Dissecting fruit and clearing

- 2 Emasculate a bud from some maternal plant. 2 to 3 days later pollinate by selfing/outcrossing or use an unfertilized fruit.
- 3 Remove the developing fruit 1 to 5 days after pollination or 2 to 3 days after emasculatation. In *Mimulus*, this protocol is useful for capturing early seed development (0- to 5 days). After 5 days, the seed tissue thickens and becomes difficult for viewing.
- 4 Pipette  **10 µl** of diluted Hoyer's solution onto a glass slide and dissect developing ovules from the fruit directly onto the glass slide using sharp forceps

- 5 After dissection, pipette about  **30 µl** of the diluted Hoyer's over the developing ovules and place a coverslip on top. Then, set the slide flat, upright in a 4°C fridge.
- 6 Depending on how developed the ovules are, clear for at least 1 to 12 hours before viewing with a Differential Interference Contrast Microscope



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited