

# DNA EXTRACTION USING PHENOL-CHLOROFORM

Santiago Montero

## Abstract

**Citation:** Santiago Montero DNA EXTRACTION USING PHENOL-CHLOROFORM. **protocols.io**

dx.doi.org/10.17504/protocols.io.q2cdyaw

**Published:** 15 Jun 2018

## Guidelines

### REAGENT PREPARATION

★ **Laird's buffer (100ml)** (adjust pH to 8,5) [store at room temperature]

Tris	1,21 g
EDTA-NA <sub>2</sub>	0,19 g
NaCl	1,17 g
SDS	0,2% (= 1 ml if 20% SDS stock is used)
ddH <sub>2</sub> O	99 ml (depends of SDS's volume)

★ **Proteinase K (20mg/ml):** Dilute 10 mg of proteinase K in 0,5 ml of ddH<sub>2</sub>O. [store at -20°C]

★ **Phenol:** stock concentration

★ **P/C/I, Phenol-Chloroform-Isoamyl (25:24:1)** [store at 4°C]

To prepare a 250 ml solution:

Phenol	125 ml
Chloroform	120 ml
Isoamyl	5 ml

★ **C/I, Chloroform-Isoamyl (24:1)** [store at 4°C]

To prepare a 250 ml solution:

Chloroform	240 ml
Isoamyl	10 ml

★ **NaCl 5M:** 29,24 g of NaCl is dissolved in H<sub>2</sub>O up to 100 ml. Use it within 1 month.

Alternatively, **NaAc 3M:** 40,8 g of NaAc is dissolved in H<sub>2</sub>O up to 100 ml. [store at RT]

★ **Ethanol 100%** or alternatively: propanol. [store at room temperature]

★ **Ethanol 70%:** To prepare 50 ml: 36,5 ml Ethanol 100% + 13,5 ml ddH<sub>2</sub>O. [store at -20°C]

## Materials

✓ EDTA by Contributed by users

✓ Ethanol 100% by Contributed by users

Phenol by Sigma Aldrich

Proteinase K E00491 by Thermo Fisher Scientific

NaCl 53014 by Sigma Aldrich

Tris-HCl (Tris-Hydrochloride), 100gm H5121 by Promega

SDS SB0485.SIZE.500g by Bio Basic Inc.

✓ double distilled water (ddH<sub>2</sub>O) by Contributed by users

Phenol-chloroform-isoamyl alcohol 25:24:1 (PCI) 15593049 by Invitrogen - Thermo Fisher

✓ Ethanol 70% by Contributed by users

## Protocol

### DAY 1

#### Step 1.

Add 500 µl Laird's buffer in 1,5 ml eppendorf tube, one for each sample.

### DAY 1

#### Step 2.

Cut 10-30 µg of muscle tissue and put it in the tube with Laird's buffer.

#### DAY 1

##### **Step 3.**

Add 20 µl Proteinase K (20 mg/ml) to each tube

#### DAY 1

##### **Step 4.**

Incubate overnight in movement at 56°C (or at least 4 hours). \* If your samples are not completely solved, add more Proteinase K and incubate for longer time. DAY 2:

#### DAY 2

##### **Step 5.**

Add 500 µl Phenol to each tube and shake heavily during 10 min.

#### DAY 2

##### **Step 6.**

Centrifuge for 10 min (4°C) at 13000 rpm.

#### DAY 2

##### **Step 7.**

Pipette and keep all of the supernatant (top layer) into a new 1,5 eppendorf tube.

#### DAY 2

##### **Step 8.**

Add 500 µl Phenol-Chloroform-Isoamyl and shake heavily during 10 min.

#### DAY 2

##### **Step 9.**

Centrifuge for 10 min (4°C) at 13000 rpm.

#### DAY 2

##### **Step 10.**

Pipette and keep all of the supernatant (top layer) into a new 1,5 eppendorf tube.

#### DAY 2

##### **Step 11.**

Add 500 µl Chloroform-Isoamyl and shake heavily during 10 min.

#### DAY 2

##### **Step 12.**

Centrifuge for 10 min (4°C) at 13000 rpm.

## DAY 2

### Step 13.

Pipette and keep all of the supernatant (top layer) into a new 1,5 eppendorf tube. Be careful not to get anything of the down layer.

## DAY 2

### Step 14.

Add 0,1 volumes of 3M NaAC to each tube and mix it softly (do not use vortex). E.g.: 40 µl to 400 µl of supernatant.

## DAY 2

### Step 15.

Add 2 volumes of ice cold 95-100% ethanol (previously stored at -20°C) to each tube E.g.: 800 µl to 400 µl of supernatant. Mix it softly (turn the tubes upside down).

## DAY 2

### Step 16.

Leave at -20°C overnight (or at least 5-6 hours).DAY3:

## DAY 3

### Step 17.

Centrifuge for 30 min (4°C) at 13000 rpm.

## DAY 3

### Step 18.

Pour off ethanol.

## DAY 3

### Step 19.

Add 1000 µl ice cold ethanol 70%.

## DAY 3

### Step 20.

Centrifuge for 15 min (4°C) at 13000 rpm.

## DAY 3

### Step 21.

Pour off ethanol.

## DAY 3

**Step 22.**

Centrifuge for 5 min (4°C) at 13000 rpm.

---

DAY 3

**Step 23.**

Pour off residual ethanol.

---

DAY 3

**Step 24.**

Dry pellet completely by leaving the tube open at room temperature.

---

DAY 3

**Step 25.**

Dissolve pellet in 100 µl ddH<sub>2</sub>O.

---