

# Bacterial transformation

Stephen Floor

## Abstract

Standard protocol to transform bacteria with a plasmid using chemically competent E. coli and antibiotic resistance.

**Citation:** Stephen Floor Bacterial transformation. **protocols.io**

[dx.doi.org/10.17504/protocols.io.mnvc5e6](https://doi.org/10.17504/protocols.io.mnvc5e6)

**Published:** 16 Jan 2018

## Materials

- ✓ MACH1 or DH5a or TOP10 or NEB stable or etc for DNA purification/miniprep by Contributed by users
- ✓ BL21(DE3) or BL21-Star(DE3) or Rosetta2(DE3) or etc for protein purification by Contributed by users

## Protocol

### Step 1.

All steps in a microcentrifuge tube.

### Step 2.

Thaw cells on ice.

### Step 3.

Incubate 10-500 ng DNA with cells (10 to 50  $\mu$ L) on ice for 25 minutes in a microcentrifuge tube.

### Step 4.

Heat shock at 42 degrees for one minute in water bath.

### Step 5.

Recover for two minutes on ice.

### Step 6.

Add 180  $\mu$ L LB or SOC media, mix.

### Step 7.

Grow at 37 degrees for one hour with shaking. Warm the plates during this step.

### Step 8.

Plate 50 to 200  $\mu$ L of the transformation on plate with appropriate antibiotic.

### Step 9.

Leave at 37 degrees overnight or room temperature for two days.