



Jan 31,
2020

Post-IMS Autofluorescence Microscopy

Elizabeth Neumann¹, Jamie Allen¹, Heath Patterson¹, Danielle Gutierrez¹, Jeff Spraggins¹

¹Vanderbilt University

1 Works for me [dx.doi.org/10.17504/protocols.io.879hxr6](https://doi.org/10.17504/protocols.io.879hxr6)

VU Biomolecular Multimodal Imaging Center

Human BioMolecular Atlas Program (HuBMAP) Method Development Community



Jamie Allen ⚡

ABSTRACT

Scope:

Obtain autofluorescence microscopy images of tissues after IMS analysis

Expected Outcome:

A brightfield and green autofluorescence microscopy image of the tissue section that enables registration and correlation of different imaging modalities on a pixel by pixel basis.

- 1 Place microscope slide within adapter and insert into the Zeiss AxioScan Slide Scanner.
- 2 Perform Coarse focusing of the tissue using:
GFP filter set (ex. 450–490 nm; em. 500–550, green)
High lamp power (~90%) and moderate exposure times (~ 150 ms).
- 3 Perform Fine focusing of the tissue using:
GFP filter set (ex. 450–490 nm; em. 500–550, green)
Brightfield lamp at similar power and exposure time used in step 2 to build a focus map.
- 4 Define the imaging region that includes the tissue.
- 5 Acquire brightfield and autofluorescence image.
- 6 Export autofluorescence image as an OME-TIF File with the following options:
“BigTIFF”, “Use Tiles”, “Compress”, and “Convert to 8 Bit”.

Alternatively, other image file types, such as “BigTIFF” and “pngs” can also be useful, depending on the application.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited