

Simple PCR-SSP- rs6656401-In4 - CR1

Gabriela Kretzschmar, Luana Caroline Oliveira, Angelica Beate Winter Boldt

Abstract

This is a simple PCR-SSP for amplification of the rs6656401, located in intron 4 (CR1 gene).

Citation: Gabriela Kretzschmar, Luana Caroline Oliveira, Angelica Beate Winter Boldt Simple PCR-SSP- rs6656401-In4 - CR1. **protocols.io**



<https://www.protocols.io/view/simple-pcr-ssp-rs6656401-in4-cr1-p44dqyw>

Published: 18 Jun 2018

Before start

- 1- Wear clean gloves;
- 2- Clean pipettes and stand with hypochlorite and 70% alcohol;
- 3- Defreeze DNA samples and reagents (exception: Taq polymerase);
- 4- Gently mix the DNA samples and pass the reagents briefly on the vortex.
- 5- Centrifuge the DNA samples and reagents with a spin;
- 6- Keep reagents on ice;
- 7- Leave microtubes ready for preparation of the mix, keep the same ones identified and on ice.
- 8- Leave application plates properly identified.

Materials

-  1x Coral Buffer by [Qiagen](#)
- ✓ dNTP mix by Contributed by users
- ✓ custom made primers by Contributed by users
- ✓ 20 ng of genomic DNA by Contributed by users
- ✓ Ultrapure Water by Contributed by users
-  Taq DNA Polymerase by [Invitrogen - Thermo Fisher](#)

Protocol

Step 1.

The amounts described in this protocol are for one sample. Multiply them for the number of samples to be investigated.

Remember:

- Each sample shall be tested in two mix reactions 1 and 2 (2 tubes).
- Always make plus 10% mix.

For this PCR we use 2 microtubes because are 2 mixes.

1. Label properly the microtubes for the mixes, identifying them: mix 1 (with A allele of SNP of rs6656401) and mix 2 (with G allele SNP of rs6656401).
2. Keep these microtubes in ice.

Step 2.

Add Ultrapure Water to the microtube of each mix.

 AMOUNT

4.696 µl : for sample

 REAGENTS

✓ Ultrapure Water by Contributed by users


Step 3.

Add Coral Buffer to each mix

 AMOUNT

0.8 µl : for sample

 REAGENTS

 1x Coral Buffer by [Qiagen](#)

Step 4.

Add 0.2 mM dNTP

 AMOUNT

0.8 µl : for sample

 REAGENTS

✓ dNTP MIX (mM) 0.2 by Contributed by users

Step 5.

Add control primer (forward and reverse) - HGH:

0.08 µl for forward

0.08 µl for reverse

5'-3' Sequence:

HGH f

TGCCTTCCCAACCATTCCCTTA

HGH r

CCACTCACGGATTCTGTTGTGTTTC

The size of the fragment is 431bp

AMOUNT

0.08 µl : for sample/ 0.1 µM

REAGENTS

✓ HGH (Control Primer) - Forward and Reverse by Contributed by users

Step 6.

Add the specific primers:

- MIX 1 - In4_rs6656401_A_f + in4_G_r
- MIX 2 - In4_rs6656401_G_f + in4_G_r

0.24µl to each primer (F and R)

5'-3' Sequence - rs6656401 - primers forwards:

CTCTGTCTCCATCTTCTCA

CTCTGTCTCCATCTTCTCG

5'-3' Sequence - intron 4 - primer reverse:

CATAGTTGTAGTTGGGGATTG

The size of the fragment is 257bp

 **AMOUNT**

0.24 μl : for sample/ 0.3 μM

 **REAGENTS**

✓ Specific Primers Forward and Reverse by Contributed by users

Step 7.

Add an aliquot of DNA sample to each tube

Kepp the plate cold.

 **AMOUNT**

1 μl : for sample

 **REAGENTS**

✓ DNA sample by Contributed by users

Step 8.

Add the Taq DNA polymerase to the mix.

Briefly vortex and centrifuge for some seconds.

 **AMOUNT**

0.064 μl : 0.04 Units/ μl

 **REAGENTS**

 Taq DNA polymerase by [Invitrogen - Thermo Fisher](#)

Step 9.

Add an aliquot of the mix to each tube.

Remember, each sample shall be tested in two mix reactions 1 and 2 (2 tubes).

The final volume for each sample is 8 μl .

 **AMOUNT**

7 μl : for sample

Step 10.

Close well the tubes and centrifuge for some seconds

Put in the thermocycler

Step 11.

Use:

Thermal cycling begins with 94°C for 5 min and 20 s; followed by 33 cycles, where each cycle began with 94 °C for 20 s and ended with 72 °C for 40 s. The annealing temperatures were 58°C for the initial 11 cycles, 55°C for the following 11 cycles and 52°C for the last 11 cycles.

Warnings

Do not leave the Taq polymerase too long outside the freezer